



Product Specification

Huntington Disease GLHDDig2X GeneProber™ Digoxigenin labeled probe Huntington Disease Genotyping by Southern blot Analysis

For research use only. Not for use in diagnostic procedures for clinical purposes

Important Notice

For research use only. Not for use in diagnostic procedures for clinical purposes.

Product to be used by experienced researchers properly trained in performing molecular biology techniques following established safety procedures; and qualified and certified for research and appropriate interpretation of results and data.

Material supplied

Huntington Disease GeneProber™ Digoxigenin Labeled Probe

Huntington disease CAG triple repeat spanning region digoxigenin labeled probe for Southern blot detection of Pst I digested DNA

| Catalog No. | Description | Size |
|-------------|--|--------|
| 40-2025-41 | Huntington Disease GLHDDig2X GeneProber™ Digoxigenin labeled | 110 µl |

One tube containing 110 µl of *GeneProber™* GLHDDig2X probe at a concentration of ~40ng/ul. This probe is digoxigenin labeled for non-radioactive detection. The quantity supplied is sufficient for at least 5 20x20 cm blots using 20µl for each blot as probe.

Storage Instructions: Upon receipt immediately store at -20°C

Huntington Disease Genotyping

Background

Huntington disease (HD) is an autosomal dominant, progressive neurodegenerative disorder with a prevalence rate of about 5-10 affected persons per 100,000 in most western populations. The disorder presents with motor impairment, cognitive deterioration, and psychiatric symptoms.

HD is caused by a CAG trinucleotide expansion within the first exon of the IT15 gene on chromosome 4p16. The expanded CAG repeats are translated into a polyglutamine tract in the Huntington protein, which is believed to cause a dominant gain of function, leading to neuronal dysfunction and neurodegeneration.

The number of CAG repeats correlates inversely with the age of onset of symptoms. The American College of Medical Genetics/American Society of Human Genetics/ Huntington Disease Genetics Testing Working Group divided genotype/phenotype correlation in the following four categories for CAG repeat lengths:

- Normal allele, ≤ 26 CAG repeats, generating a normal phenotype;
- Intermediate allele, 27-35 CAG repeats, mutable normal allele generating a normal phenotype;
- HD allele with reduced penetrance, 36-39 CAG repeats, generating a normal or HD phenotype;
- HD allele, ≥ 40 CAG repeats, generating a HD phenotype.

The CAG trinucleotide expansion is unstable and can lengthen during transmission from parents to offspring. Thus, the stage of onset can decrease from one generation to the next, a phenomenon known as anticipation. HD anticipation is more intense in paternal transmission.

HD L34202 digested with: BamHI, PstI




Table 1. Trinucleotide Repeats in Human Genetic Disease

| Disease | Repeat ^a | Normal Length ^b | Intermediate Length (Premutation) ^{a,b} | Full Disease Length ^b |
|--|---------------------|----------------------------|--|-----------------------------------|
| Fragile XA (FRAXA) | (CGG) _n | 6-52 | 59-230 | 230-2,000 |
| Fragile XE (FRAXE) | (CCG) _n | 4-39 | ? (31-61) | 200-900 |
| Fragile XF (FRAXF) | (CGG) _n | 7-40 | ? | 306-1,008 |
| FRA16A | (CCG) _n | 16-49 | ? | 1,000-1,900 |
| Jacobsen Syndrome (FRA11B) | (CGC) _n | 11 | 80 | 100-1,000 |
| Kennedy Syndrome (SMBA) | (CAG) _n | 14-32 | ? | 40-55 |
| Myotonic Dystrophy (DM) | (CTG) _n | 5-37 | 50-80 | 80-1,000; congenital, 2,000-3,000 |
| Huntington disease (HD) | (CAG) _n | 10-34 | 36-39 | 40-121 |
| Spinocerebellar ataxia 1 (SCA1) | (CAG) _n | 6-39 | None Reported | 40-81 |
| Spinocerebellar ataxia 2 (SCA2) | (CAG) _n | 14-31 | None Reported | 34-59 |
| Spinocerebellar ataxia 3 (SCA3)/Machado Joseph disease (MJD) | (CAG) _n | 13-44 | None Reported | 60-84 |
| Spinocerebellar ataxia 6 (SCA6) | (CAG) _n | 4-18 | None Reported | 21-28 |
| Spinocerebellar ataxia 7 (SCA7) | (CAG) _n | 7-17 | 28-35 | 38-130 |
| Haw River syndrome (HRS; also DRPLA)) | (CAG) _n | 7-25 | ? | 49-75 |
| Friedreich ataxia (FRDA) | (GAA) _n | 6-29 | ? (>34-40) | 200-900 |

^a Typically, repeats tracts contain sequence interruptions. See Pearson and Sinden (1998a) for a discussion of the sequence interruptions.

^b No. of triplet repeats.

^c A question mark (?) indicates potential mutagenic intermediate length, and an ellipsis (...) indicates none. Not all diseases are associated with a premutation length repeats tract or premutation disease condition.


Molecular Analysis

The detection of expansion of a region of DNA sequence can be detected by PCR and Southern blotting procedures. These methods can be used for all disorders involving increase in size of a region of DNA. DNA analysis for direct detection of CAG expansion in Huntington Disease is based on enzymatic amplification of a fragment containing the CAG repeat sequence in exon I of the HD gene. This test detects the CAG expansion by the size of the amplified product; an increase in size is correlated with the corresponding number of CAG repeats and a calculated risk factor. Normal individuals have repeat numbers of up to 30, while individuals with a high probability of developing HD carry more than 37 repeats. Individuals with 30-37 repeats have a high probability of passing on repeats in the pathological size range.

Polymerase Chain Reaction (PCR) based methods are fundamentally similar. The two primers are constructed such that they span the region of the CAG trinucleotide repeat region. PCR is the most common method used to estimate the number of CAG repeats. Since the CAG repeats in the HD gene are immediately 5' of a CCG repeat which is also polymorphic in length, the PCR product of this primer pair excludes the known adjacent polymorphic CCG repeat that can contribute to an inaccurate determination of HD gene CAG repeat sizes in individuals who may have an HD gene CAG repeat allele close to the normal/affected boundary. Refer to HD GScan™ fluorescent PCR genotyping system [catalog No.: 40-2025-15FM]. Full mutations and PCR generated results that signify large expansions should be confirmed by Southern blot analysis.

Southern blot analysis for HD genotyping is performed by digestion of genomic DNA by Pst I followed by hybridization to a probe cognate to the Pst I fragment region that contains the amplified CAG repeats.

This manual details the Southern blot HD genotyping procedure using the GLHDDig2X GeneProber™.



Procedure: Chemiluminescent Southern Protocol

Material Supplied

One tube containing 110 µl of *GeneProber™* GLHDDig2X probe at a concentration of ~40ng/µl. This probe is digoxigenin labeled for non-radioactive detection. The quantity supplied is sufficient for at least five 20x20 cm blots using 20 µl for each blot as probe.

A. Genomic DNA Digestion

Important Note

-Digest genomic DNA with Pst I when using GLHDDig2X GeneProber™ as labeled probe.

| Restriction Digestion | |
|-----------------------------|-----------------|
| Component | Volume Quantity |
| Genomic DNA | 5 to 10µg |
| 10x Pst I Buffer | 10 µl |
| Pst I (~40 u/µl) | 4 µl |
| H₂O to | 100 µl |
| Overnight digestion at 37°C | |

Ethanol precipitate the digests; dissolve the pellets in 10 µl of 1x Loading buffer.

B. Electrophoresis and Transfer

1. Load samples to a 0.8% agarose gel. Electrophorese over night at 45mA for 20-24 hours. (1.6 kb fragment on the bottom of the gel).
2. Depurinate with 0.25N HCl (add 10 ml HCl to 500 ml H₂O) for 10 minutes.
3. Denature the DNA with 0.4N NaOH/0.6M NaCl for 30 min. at room temperature (RT).
4. Neutralize with 1.5M NaCl/0.5M Tris (pH 7.5) for 30 min. at RT.
5. Transfer overnight by Southern blot procedure to positively charged nylon membrane using 10xSSC.
6. Wash the membrane with 2x SSC and then bake at 80°C for 2 hours.

C. Hybridization

Gene Link recommends using Roche Biochemicals Digoxigenin based washing and detection system reagents.

1. Perform prehybridization at 55°C for 3 hours in 10 ml of Easy Hyb buffer (Roche Biochemicals).
2. Boil 20µl *GeneProber™* GLHDDig2X probe in 500µl of Easy Hyb for 10 minutes. Chill directly on ice.
3. Add the above probe to 10ml of Easy Hyb.
4. Discard the prehybridization buffer and replace it with the hybridization buffer containing the boiled probe. Hybridize overnight at 55°C.
5. Washing. Wash the membrane in 2xSSC/0.1% SDS at RT twice (5 min/wash), 0.5xSSC, 0.1%SDS twice at 65°C to 70°C (15 min/wash).
6. Warm the blocking reagent at this point. Prepare fresh 100 ml of Buffer MB by adding 10 ml of 10% blocking reagent [Gene Link Catalog no: 40-5026-10; Blocking solution for hybridization (10%)] and 10 ml of Maleic acid buffer 10X (Buffer M 10X) [Gene Link Catalog no: 40-5025-20; Maleic acid buffer 10X (Buffer M 10X)] to 80 mL of sterile water. Use 80 ml for blocking, the rest of 20 ml for making Anti-DIG-AP conjugate.

D. Anti-Dig Alkaline Phosphatase Binding

1. Equilibrate the membrane in 100ml of 1x washing buffer M for 1 minute.
2. Incubate the membrane in 80ml of Buffer MB (prepared in step 6 above) blocking solution at RT for 30 min.
3. Prepare 1:10000 Anti-DIG-AP conjugate at this point. Example, add 2 μ l to 20 ml Buffer MB (prepared in step 6 above).
4. Incubate the membrane in 20ml of Anti-DIG-AP conjugate solution at RT for 30 min.
5. Wash the membrane twice, 15 min/wash in 200ml of 1x washing buffer M at RT.
6. Equilibrate the membrane in 50ml of 1x Detection buffer for 2 min. Repeat 2 times. Total of three washes.

E. Detection

Detection with CDP star (Tropix) as substrate will yield reliable result by exposing to Kodak X-OMAT or XAR X-ray film for 1 hour to overnight at room temperature.

Transfer blot to a plastic sheet, (sheet protector cut from two sides to open up) and drain off excess buffer. Wipe off edges with paper towel. Blot should not be allowed to dry.

Spray CDP-star ready-to-use substrate evenly to cover the blot. DO NOT OVER SPRAY. Cover the blot with plastic sheet and wipe entire surface of the covered blot to expel any excess substrate and air bubbles. Expose the film at room temperature for 1 hr. or for shorter or longer time as required.


Luminescence continues for at least 24 hours and signal intensity remains almost constant during the first few hours. Multiple exposures can be taken to achieve the desired signal strength.

F. Stripping

Wash the membrane in water to remove the substrate. Then wash the membrane in 0.2N NaOH/0.1% SDS at 37°C for 30 minutes. Rinse the membrane in 2XSSC. Air dry.

G. Results & Interpretation

1. Normal hybridization pattern is ~1419 bp fragment with genomic DNA digested with Pst I using GLHDDig2X GeneProber™ as labeled probe.
2. Larger fragment size is attributable to expanded CAG repeat region. See HD probe region figure below.



Required reagents with recommended suppliers

Roche Applied Science

<http://www.roche-applied-science.com>

| Product Description | Catalog Number |
|---|----------------|
| Nylon Membranes, positively charged ; 20 sheets 10 x 15 cm | 11209272001 |
| DNA Molecular Weight Marker III, DIG-labeled ; 500 µl 10 µg/ml 5 µg | 11218603910 |
| DIG Easy Hyb ; 500 ml | 11603558001 |
| DIG Wash and Block Buffer Set ; 1 set 30 blots | 11585762001 |
| Anti-Digoxigenin-AP, Fab fragments from sheep; 200 µl 150 U | 11093274910 |
| CDP Star Ready to use; 2X 50 mL | 12041677001 |

Gene Link

<http://www.genelink.com/geneprodsite/category.asp?c=44>

Non-radioactive Southern Blot Reagents

| Product Description | Catalog No. | Unit Size |
|--|-------------|-----------|
| <u>Agarose LE Molecular Biology Grade 100 gms</u> | 40-3010-10 | 100 gms |
| <u>TAE Buffer 50 X Concentrate 1000 ml</u> | 40-3007-10 | 1 L |
| <u>TBE Buffer 5 X Concentrate; 1L</u> | 40-3008-10 | 1 L |
| <u>Loading buffer 10X BPB/XC non-denaturing; 1mL</u> | 40-3003-10 | 1 mL |
| <u>Loading buffer 10X BPB/XC non-denaturing ; 15 mL</u> | 40-3003-15 | 15 mL |
| <u>Lumisol II, Hybridization Solution; 200 mL</u> | 40-5023-20 | 200 mL |
| Depurination Solution (2X) for Southern Blotting; 1 L | 40-5034-10 | 1 L |
| Denaturation Solution (2X) for Southern Blotting; 1L | 40-5035-10 | 1 L |
| Neutralization Solution (2X) for Southern Blotting; 1L | 40-5036-10 | 1 L |
| Hybwash A; Hybridization Wash Solution Concentrate (20X SSC); 250 mL | 40-5020-25 | 250 mL |
| Hybwash B, Hybridization Wash Solution (10% SDS) ; 100 mL | 40-5021-10 | 100 mL |
| Maleic acid buffer 10X (Buffer M 10X); 200 mL | 40-5025-20 | 200 mL |
| 10% Blocking Reagent; 200 mL | 40-5026-20 | 200 mL |
| Detection Buffer 10X; Alkaline Phosphatase detection buffer; 100 mL | 40-5031-10 | 100 mL |
| CDP-Star® Substrate; Ready-to-Use 0.25 mM in spray bottle; 20 mL | 40-5010-20 | 20 mL |

Reagent Preparation

Most reagents with composition listed below are available in a molecular biology laboratory or these can be prepared in house. Gene Link catalog numbers are also listed if you like to purchase these common reagents.

| Depurination Solution (0.25M HCl) | | |
|--|--------------------|---------------|
| Product Description | Catalog No. | Volume |
| Depurination Solution (2X) for Southern Blotting | 40-5034-10 | 150 mL |
| Sterile water | | 150 mL |
| Total Volume | | 300 mL |

| Denaturation Solution (0.5M NaOH, 1.5M NaCl) | | |
|---|--------------------|---------------|
| Product Description | Catalog No. | Volume |
| Denaturation Solution (2X) for Southern Blotting | 40-5035-10 | 150 mL |
| Sterile water | | 150 mL |
| Total Volume | | 300 mL |

| Neutralization Solution (0.5M Tris-HCl pH 7.5, 1.5M NaCl) | | |
|--|--------------------|---------------|
| Product Description | Catalog No. | Volume |
| Neutralization Solution (2X) for Southern Blotting | 40-5036-10 | 150 mL |
| Sterile water | | 150 mL |
| Total Volume | | 300 mL |

| Hybwash I (2xSSC, 0.1% SDS) | | |
|--|--------------------|---------------|
| Product Description | Catalog No. | Volume |
| Hybwash A; Hybridization Wash Solution Concentrate (20X SSC) | 40-5020-25 | 35 mL |
| Sterile water | | 311 mL* |
| Hybwash B; Hybridization Wash Solution Concentrate (10% SDS) | 40-5021-10 | 4 mL* |
| Total Volume | | 350 mL |

* Volumes adjusted to whole numbers

| Hybwash II (0.5xSSC, 0.1%SDS) | | |
|--|--------------------|---------------|
| Product Description | Catalog No. | Volume |
| Hybwash A; Hybridization Wash Solution Concentrate (20X SSC) | 40-5020-25 | 9 mL* |
| Sterile water | | 337 mL |
| Hybwash B; Hybridization Wash Solution Concentrate (10% SDS) | 40-5021-10 | 4 mL* |
| Total Volume | | 351 mL |

* Volumes adjusted to whole numbers

1X Maleic Acid Buffer (Buffer M 1X)
(100 mM Maleic acid, 150 mM NaCl pH7.5)

| Product Description | Catalog No. | Volume |
|---------------------------------------|-------------|---------------|
| Maleic acid buffer 10X (Buffer M 10X) | 40-5025-20 | 10 mL |
| Sterile water | | 90 mL |
| Total Volume | | 100 mL |

Buffer MB
(1 x Maleic acid buffer (Buffer M) with Blocking Reagent)

Always prepare fresh!

| Product Description | Catalog No. | Volume |
|---------------------------------------|-------------|---------------|
| Maleic acid buffer 10X (Buffer M 10X) | 40-5025-20 | 10 mL |
| Sterile water | | 80 mL |
| 10% Blocking Reagent* | 40-5026-10 | 10 mL |
| Total Volume | | 100 mL |

The prepared reagent will be turbid yellow in color

* The 10% Blocking Reagent is turbid yellow in color and will form precipitates on storage.

Warm to 50°C and shake well before aliquoting. **DO NOT SHAKE VIGOROUSLY**

1X Detection Buffer, Alkaline phosphatase detection buffer
(100mM Tris-HCl pH 9.5, 100mM NaCl)

| Product Description | Catalog No. | Volume |
|---|-------------|---------------|
| Detection Buffer 10X; Alkaline phosphatase detection buffer | 40-5031-10 | 10 mL |
| Sterile water | | 90 mL |
| Total Volume | | 100 mL |

Appendix: Protocols

Genomic DNA Purification

Genomic DNA is usually extracted from blood. A simple procedure is given below that purifies ~10 µg DNA from 300 µl blood using a 30 minute procedure.

Omni-Pure™ Genomic DNA Purification System Catalog Number: 40-4010-01
Rapid DNA Purification Protocol for 300 µl Whole Blood

A. Initial Preparation

1. Label two sets of 1.5 ml tubes per sample.
2. Add 900 µl GD-1 solution (RBC Lysis Solution) to one tube for each sample.
3. Add 300 µl Isopropanol (2-propanol) to one tube for each sample. Cap the tubes.

B. Cell Lysis

1. To the tube containing 900 µl GD-1 solution (RBC Lysis Solution) using a filter tip pipet transfer 300 µl whole blood. Cap and gently mix by inversion. Incubate for 1-3 minutes at room temperature. Mix by inversion a few times during this incubation period. Incubate longer for fresh blood cells as they are intact and not lysed already.
2. Centrifuge at 3 K rpm for 20 seconds to pellet the white blood cells. A reddish white pellet should be clearly visible. Decant and discard supernatant leaving behind the last few droplets. Do not totally remove the supernatant.
3. Completely resuspend the white blood cell pellet by vigorously vortexing the tube. Ensure that the pellet is completely resuspended.
4. To the resuspended cells add 300 µl GD-2 solution (Cell Lysis Solution). Mix by gentle vortexing. You will notice release of DNA by the thickening of the liquid in the sample. Samples may be stored at this stage for processing later. It has been shown that the samples are stable in Cell Lysis Solution for at least 2 years at room temperature.

C. Protein Precipitation

1. Add 100 µl GD-3 solution (Protein Precipitation Solution) to the sample in cell lysis solution.
2. Vortex vigorously for 20 seconds. Small particles of brown color will appear and be visible at this stage.
3. Centrifuge at 5 K rpm for 1 minute to pellet the precipitated proteins. A clearly visible brown pellet containing proteins should be collected at the bottom of the tube.

D. DNA Precipitation

1. Decant the supernatant containing the DNA to a new appropriately labeled tube (see initial preparation above) containing 300 µl 100% Isopropanol (2-propanol).
2. Mix the sample by inversion until a visible white floating DNA strand-particle is identified. Mixing by inversion 30-40 is usually sufficient.
3. Centrifuge at 6 K rpm for 1 minute to collect the DNA as a pellet. A white DNA pellet should be clearly visible.
4. Decant supernatant and place tube inverted on a clean Kimwipe™ tissue paper to drain the remaining supernatant.
5. To remove residual salts, add 300 µl of 70% ethanol. Vortex gently.
6. Centrifuge at 6 K rpm for 1 minute to collect the DNA as a pellet. Gently take out the tubes so that the pellet is not dislodged. While holding the tube, rotate tube so that you can watch the pellet. Now carefully decant the ethanol, keeping an eye on the pellet so that it does not flow away.
7. Place tube inverted on a clean Kimwipe™ tissue paper to drain the remaining ethanol.
8. Air dry the DNA pellet. Do not use vacuum.

E. DNA Reconstitution & Use

1. Add 100 µl of GD-4 solution (DNA Reconstitution Solution). Vortex gently. Incubate at 60°C for 5 minutes to facilitate dissolution or keep overnight at room temperature.
2. Store DNA at 4 °C. For long-term storage, place sample at -20 °C or -80 °C.
3. Average yield of 10 µg is expected from 300 µl blood DNA. The range is between 5 µg to 15 µg.
4. The 100 µl of purified DNA obtained will have an average concentration of ~100 ng/µl.
5. For PCR amplification use 1-2 µl.
6. Use 100 µl for restriction digestion followed by Southern blot analysis.
7. It is convenient to perform multiple 300 µl blood DNA purification instead of scaling up the procedure.

Gel Electrophoresis of DNA

Gel electrophoresis of PCR products is the standard method for analyzing reaction quality and yield. PCR products can range up to 10 kb in length, but the majority of amplifications are at 1 kb and below. Agarose electrophoresis is the classical method to analyze amplification products from 150 bp to greater than 10 kb. Polyacrylamide gel electrophoresis should be used for resolution of short fragments in the range of 100 bp to 500 bp when discrimination of as small as a 10 bp difference is required.

PAGE gels for PCR products can be formulated with the amount of cross-linker chosen to give pore sizes optimal for the size of DNA fragment desired. Gels are most often stained in ethidium bromide, even though the fluorescence of this stain is quenched by polyacrylamide, which decreases sensitivity 2-5 fold. This decrease in sensitivity generally does not present a problem, because most PCR reactions yield product levels in the microgram range, and ethidium will detect as little as 1/10 of this amount. Polyacrylamide gels can be stained by silver staining for more sensitive detection.

Agarose Gel Electrophoresis of DNA

Agarose gels are typically run at 20 to 150V. The upper voltage limit is the amount of heat produced. At room temperature about 5 Watts is correct for a minigel (Volts x Amps = Watts). At low voltages migration is linearly proportional to voltage, but long DNA molecules migrate relatively faster in stronger fields. Migration is inversely proportional to the log of the fragment length; a log function also governs migration rate and gel concentration (0.5 to 2% for most purposes). Furthermore, supercoiled / circular DNA molecules migrate at different rates from linear molecules; single-stranded DNA and RNA migrate at similar rates, but usually faster than double-stranded DNA of the same length. Salt in the samples increases conductivity and, hence, migration rate.

The buffers used for most neutral agarose gels (the gel itself and the solution in which it lies) is 1 x TAE or 1 x TBE. Agarose powder is added to the buffer at room temperature, heated in a microwave and boiled slowly until the powder has dissolved. Cast the gel on a horizontal surface once the agarose has been cooled to ca. 60° C (just cool enough to hold) and add 0.1 µg of ethidium bromide solution for each ml of gel volume. At times, during removal of the comb, it is possible to tear the bottom of the sample wells gels, which results in sample leakage upon loading. This can be avoided by removing the comb after the gel has been placed in the running buffer.

- Use TAE buffer for most molecular biology agarose gel electrophoresis.

Recipe


| 1 X TAE Buffer |
|---|
| Agarose Gel Electrophoresis Buffer |
| 40 mM Tris-Acetate pH 7.8 |
| 1 mM EDTA |

| 1 X TBE |
|--|
| Agarose and Polyacrylamide Gel Electrophoresis Buffer |
| 0.089 M Tris |
| 0.089 M Boric Acid |
| 0.002 M EDTA |

Spectrophotometric Determination of DNA Concentration & Estimation by Agarose Gel Electrophoresis

Measuring the optical density (OD) or absorbance at 260 nm (A_{260}) in a UV spectrophotometer is a relatively accurate method for calculating the concentration of DNA in an aqueous solution if a standard curve is meticulously prepared. An A_{260} of 1, using a 1 cm path length, corresponds to a DNA concentration of 50 $\mu\text{g/ml}$ for double stranded DNA, 40 $\mu\text{g/ml}$ for RNA and 33 $\mu\text{g/ml}$ for oligonucleotides. However, this method is not suitable for determining concentrations of dilute solutions of DNA, as the sensitivity of this method is not very high. For reliable readings, the concentration of double stranded DNA must be greater than 1 $\mu\text{g/ml}$. A simple, inexpensive method for the estimation of nanogram quantities of DNA is described in the following section. We recommend the use of agarose gel electrophoresis for routine approximate determination of DNA concentration.

The amount of DNA in a sample may be estimated by running the sample alongside standards containing known amounts of the same-sized DNA fragment. In the presence of ethidium bromide staining, the amount of sample DNA can be visually estimated by comparing the band intensity with that of the known standards.



An unknown amount of a 4 kb DNA fragment (unknown) was run alongside known quantities (indicated in nanograms) of the same DNA fragment. As estimated by visual comparison with the known standards, the unknown sample contained 240-320 ng of DNA.



Ethidium bromide is a carcinogen. Follow Health and Safety Procedures established by your institution. Follow proper Hazardous Material Disposal procedures established by your institution.

- Use 0.1 μg of ethidium bromide solution for each ml of gel volume.

Polymerase Chain Reaction

PCR Components and Analysis

PCR buffer conditions vary and it is imperative to optimize buffer conditions for each amplification reaction. At Gene Link most amplification reactions have been optimized to work with the following standard buffer condition, unless otherwise indicated.

dNTP Concentration

Standard dNTP concentration of 0.2 mM of each base is used. See section on PCR additives when dNTP concentration is changed.

MgCl₂ Concentration

The concentration of Mg²⁺ will vary from 1-5 mM, depending upon primers and substrate. Since Mg²⁺ ions form complexes with dNTPs, primers and DNA templates, the optimal concentration of MgCl₂ has to be selected for each experiment. Low Mg²⁺ ions result in a low yield of PCR product, and high concentrations increase the yield of non-specific products and promote mis-incorporation. Lower Mg²⁺ concentrations are desirable when fidelity of DNA synthesis is critical. The recommended range of MgCl₂ concentration is 1-4 mM, under the standard reaction conditions specified. At Gene Link, using the standard PCR buffer with KCl, a final dNTP concentration of 0.2 mM, a MgCl₂ concentration of 1.5 is used in most cases. If the DNA samples contain EDTA or other chelators, the MgCl₂ concentration in the reaction mixture should be raised proportionally. Given below is an MgCl₂ concentration calculation and addition table using a stock solution of 25 mM MgCl₂.

| MgCl ₂ Concentration & Addition Table | | | | | | | | |
|--|-----|------|-----|------|-----|-----|-----|-----|
| Final concentration of MgCl ₂ in 50 µl reaction mix, (mM) | 1.0 | 1.25 | 1.5 | 1.75 | 2.0 | 2.5 | 3.0 | 4.0 |
| Volume of 25 mM MgCl ₂ , µl | 2 | 2.5 | 3 | 3.5 | 4 | 5 | 6 | 8 |

Primer Concentration

The final concentration of primers in a PCR reaction is usually 0.5 to 1 µM (micromolar). This is equivalent to 0.5 to 1 pmol/µl. For a 100 µl reaction you would add 50 to 100 pmols. At Gene Link we use 0.5 pmol/µl in the final PCR.

Genemer™ Reconstitution

Stock Primer Mix: Dissolve the supplied 10 nmols of lyophilized Genemer™ in 100 µl sterile TE. The 10 nmols of primers when dissolved in 100 µl will give a solution of 100 µM i.e. 100 pmols/µl.

Primer Mix: Prepare a 10 pmols/µl Primer Mix solution by a ten fold dilution of the stock primer mix.

Example: Add 180 µl sterile TE to a new tube, to this tube add 20 µl of primer stock solution. Label this tube as Primer Mix 10 pmols/µl.

Amplification Thermal Cycling

Hot Start: It is essential to have a 'Hot Start' profile for amplification of any

● Recipe

| Standard Gene Link PCR Buffer Composition | |
|---|----------------|
| 10 X PCR buffer | 1 X PCR buffer |
| 100 mM Tris-HCl pH 8.3 | 10 mM |
| 500 mM KCl | 50 mM |
| 15 mM MgCl ₂ | 1.5 mM |
| 0.01% Gelatin | 0.001% |

● Recipe

| 2.0 mM dNTP Stock Solution Preparation* | |
|---|--------|
| Component | Volume |
| 100 mM dGTP | 100 µl |
| 100 mM dATP | 100 µl |
| 100 mM dTTP | 100 µl |
| 100 mM dCTP | 100 µl |
| Water | 4.6 ml |
| Total Volume | 5 ml |

*Aliquot and freeze



Always use filter barrier pipette tips to prevent cross contamination

● Recipe

| TE Buffer pH 7.5 Composition |
|------------------------------|
| 1 X TE Buffer pH 7.5 |
| 10 mM Tris-HCl pH 7.5 |
| 1 mM EDTA |



• Program your thermal cycler instrument with an amplification profile

fragment from a complex template like human genomic DNA. Taq polymerase has low activity at room temperature and it is essential to minimize any mis-priming in the first cycle of amplification. A typical hot start profile is given below. Various enzyme preparations are available which are activated by heat in the first cycle. A simple hot start protocol is given below that can be used with regular Taq polymerase. See the section on PCR additives for amplification of products from high GC content templates.

| Hot Start | | |
|---|---------------------|--------|
| Step | Time & Temperature | Cycles |
| Initial Denaturation | 95 °C for 5 minutes | 1 |
| Annealing | 60 °C Hold Infinity | Hold |
| Comments: Add Taq premix while on hold. | | |

Amplification File

The initial denaturation step at 94 °C for 30 seconds is sufficient for all templates. The number of cycles is usually set to 30 and is sufficient to amplify 1-10 µg of product depending on the initial concentration of template. A higher number of cycles from 35-45 cycles may be used, but internal priming on the product and over amplification of unwanted bands often result from over-cycling. Generally, it is better to focus on optimizing reaction conditions than to go beyond 35 cycles.

| Typical Amplification File | | | |
|---|-------------|-----------|--------|
| Step | Temperature | Time | Cycles |
| Denaturation | 94 °C | 30 sec. | 30 |
| Annealing | * | 30 sec. | |
| Elongation | 72 °C | 30 sec. | |
| Fill in Extension | 72 °C | 7 minutes | 1 |
| Hold | 4 °C | Infinity | Hold |
| *Based on the Tm of the primers. Usually varies from 50 °C to 65 °C | | | |

prior to beginning the amplification protocol. Consult your appropriate instrument manufacturer's manual.

● Recipe

| Typical PCR Premix (/50µl) | |
|--|--------------|
| Component | Volume |
| 10 x PCR Buffer | 5 µl |
| 2.0 mM dNTP mix (each) | 5 µl |
| Primer Mix (10 pmol/µl each) or 2.5µl of 10 pmol/µl of individual primer (final 25 pmol of each primer/50µl) | 2.5 µl |
| H ₂ O | 37.5 µl |
| Total Volume | 50 µl |

● Recipe

| PCR reaction (/50µl) | |
|------------------------|--------|
| Component | Volume |
| PCR premix | 45 µl |
| 100ng/µl diluted DNA | 1 µl |
| Hot start and then add | |
| Taq premix | 5 µl |

● Recipe

| Taq Premix EM (/50µl) | |
|--|--------|
| Component | Volume |
| PCR Premix | 6 µl |
| Taq polymerase (5 u/µl) | 0.25µl |
| Add 5 µl/50 µl rxn after initial denaturation. | |
| Use 2.5 units of Taq for 100 µl reactions. Taq is usually supplied at a concentration of 5 units/µl | |

| PCR Premix Preparation (PP) | | |
|---|----------------|-----------------------------|
| Component | 1 X 50 µl Rxn. | 10 X 50 µl Rxns. |
| Sterile Water | 32 µl | 320 µl |
| 10 X PCR Buffer | 4.5 µl | 45 µl |
| 2.0 mM dNTP | 5 µl | 50 µl |
| 10 pmol/µl Primer Mix | 2.5 µl | 25 µl |
| Taq Enzyme Mix (EM) | 5 µl | 50 µl |
| Template DNA (~500 ng) | 1-2 µl | Add 1-2 µl DNA to each tube |
| Total Volume | 50 µl | |
| Keep on ice during set up. After adding template start PCR File | | |

i • The PCR premix preparation protocol is written considering that more than one amplification reaction will be performed at the same time. If only one reaction is planned then there is no need to prepare the Taq Enzyme Mix (EM).


Recipe

| Gene Link PCR Buffer |
|--------------------------|
| 1 X PCR Buffer |
| 10 mM Tris-HCl pH 8.3 |
| 50 mM KCl |
| 1.5 mM MgCl ₂ |
| 0.001% Gelatin |

Yield and Kinetics

The target will be amplified by up to 10⁶ fold in a successful reaction, but the amplification will usually plateau at 1-10 µg. Thus, 1 pg of target sequence in the reaction is a good place to begin.

PCR reactions produce product in a nonlinear pattern. Amplification follows a typical exponential curve until some saturation point is reached. Generally products will not be further amplified once 1-5 µg has been generated. Saturation by one product of a reaction does not always prevent further amplification of other generally unwanted products. Over-cycling may decrease the quality of an otherwise good reaction. When first optimizing a reaction, it is advisable to take samples every 5 or 10 cycles to determine the number of cycles actually needed.



Purification of PCR Product

Various purification methods are available for the purification of PCR products. The selection of a particular method over another is based on the downstream application and the initial robustness of the amplification. Usually no further purification is required for most cloning experiments if a single fragment is amplified, whereas for sequencing applications the amplified product should be purified from the primers and any other minor amplification products.

The preferred method of purification of an amplified fragment is the excision of the fragment band after agarose gel electrophoresis. This method yields the purification of a single fragment; as such care should be taken to excise a gel piece containing a single electrophoretically resolved fragment. The Omni-Clean™ Purification System available from Gene Link can be used for this purpose. Catalog No. 40-4110-10 for bead based system; 40-4120-10 for spin column based system and 40-4130-10 for DNA concentration. Please refer to product insert for detailed protocol or visit www.genelink.com.

A. Purification of DNA from gel slices using glass beads. Provides purified single fragment.

[Omni-Clean™ Gel DNA Beads Purification System; Catalog No. 40-4110-10]

Protocol

1. By weight, determine the volume of the excised DNA fragment.
2. Add 3 volumes of NaI solution and heat to 55 °C. Visually determine the dissolution of gel pieces.
3. Add 1 µl of glass bead suspension per µg of DNA and vortex.
4. Centrifuge at 2K rpm for 20 seconds to pellet glass bead/DNA complex. Discard supernatant.
5. Re-suspend pellet in 400 µl Omni-Clean™ wash buffer. Centrifuge at 2K rpm for 20 seconds and discard wash buffer.
6. Pipet out any remaining buffer in the tube.
7. Add 25 µl water or TE; re-suspend pellet and centrifuge at 2K rpm for 20 seconds.
8. The supernatant contains the purified DNA. Using a pipet, collect the supernatant and transfer to a new appropriately labeled tube.

B. Purification of DNA from gel slices using spin column. Provides purified single fragment.

[Omni-Clean™ Gel DNA Spin Column Purification System; Catalog No. 40-4120-50]

Protocol

1. By weight, determine the volume of the excised DNA fragment.
2. Add 3 volumes of NaI solution and heat to 55 °C. Visually determine the dissolution of gel pieces.
3. Add the above solution to the spin column assembled on a collection tube.
4. Let the solution flow by gravity or centrifuge at 2K rpm for 20 seconds. Discard flow through collected in the collection tube.
5. Add 400 µl Omni-Clean™ wash buffer to the spin column. Centrifuge at 2K rpm for 2 minutes and discard wash buffer collected in the collection tube.
6. Replace the collection tube with a new appropriately labeled 1.5ml tube.
7. Add 25 µl water or TE to the spin column. Let sit for 3 minutes.
8. Centrifuge at 2K rpm for 2 minutes.
9. The collection tube contains the purified DNA.

C. Purification of DNA from solution using glass beads. Provides removal of salts, primers and dNTP.

[Omni-Clean™ DNA Beads Concentration System; Catalog No. 40-4130-10]

Protocol

1. Determine volume of DNA solution and add 3 volumes of NaI solution.
2. Add 1 µl of glass bead suspension per µg of DNA.
3. Centrifuge at 2K rpm for 20 seconds to pellet glass bead/DNA complex. Discard supernatant.
4. Re-suspend pellet in 400 µl Omni-Clean™ wash buffer.
5. Centrifuge at 2K rpm for 20 seconds and discard wash buffer.
6. Pipet out any remaining buffer in the tube.
7. Add 25 µl water or TE; re-suspend pellet and centrifuge at 2K rpm for 20 seconds.
8. The supernatant contains the purified DNA. Using a pipet, collect the supernatant and transfer to a new appropriately labeled tube.

D. Purification of DNA from solution using spin column. Provides removal of salts, primers and dNTP.

[Omni-Clean™ DNA Spin Column Concentration System; Catalog No. 40-4140-10]

Protocol

1. Determine volume of DNA solution and add 3 volumes of NaI solution.
2. Add the above solution to the spin column assembled on a collection tube.
3. Let the solution flow by gravity or centrifuge at 2K rpm for 20 seconds. Discard flow through collected in the collection tube.

4. Add 400 µl Omni-Clean™ wash buffer to the spin column. Centrifuge at 2K rpm for 2 minutes and discard wash buffer collected in the collection tube.
5. Replace the collection tube with a new appropriately labeled 1.5ml tube.
6. Add 25 µl water or TE to the spin column. Let sit for 3 minutes.
7. Centrifuge at 2K rpm for 2 minutes.
8. The collection tube contains the purified DNA.

PEG Precipitation

Primers and salts are efficiently removed by a simple PEG precipitation. This method is recommended for downstream DNA sequencing application. This method is generally used for plasmid DNA.

Protocol

1. To 50 µl of amplified PCR reaction add 6.0 µl of 5 M NaCl and 40 µl of 13% (w/v) PEG 8000. Incubate the mixture on ice for 20-30 minutes.
2. Collect the DNA precipitate by centrifugation at maximum speed for 15 minutes at 4 °C in a microfuge. Carefully remove the supernatant by gentle aspiration.
The pellet of DNA is translucent and generally invisible at this stage.
3. Rinse the pellet with 500 µl of 70% ethanol.
The precipitate changes to a milky-white color and becomes visible.
4. Carefully pour off the 70% ethanol. Rinse the DNA pellet once more with 70% ethanol. Store the tube in an inverted position at room temperature until the last visible traces of ethanol have evaporated.
5. Dissolve the DNA in 20 µl of H₂O.
6. Run an aliquot on an agarose gel to confirm the presence of the correct amplified product. The purified DNA is sequence grade and can be used directly for sequencing.

Gel Filtration

Primers and salts are efficiently removed by gel filtration using Sephadex G-50. This method is recommended for downstream DNA sequencing application.

Protocol

1. Hydrate Sephadex G-50 ahead of time in sterile water or TE (10mM Tris pH 8, 1 mM EDTA). Take out from fridge if already stored hydrated. Bring to room temperature.
2. Assemble a spin column on a collection tube.
3. Add 700 µl of hydrated Sephadex G-50 to each spin column, initiate flow using rubber bulb or any other method.
4. Allow flowing by gravity till there is no more fluid left above the Sephadex G-50 bed. Discard flow through from the collection tube.
5. Spin the spin column placed inside the collection tube for 2 minutes at 3 K rpm.
6. Change collection tube to new 1.5 ml tube appropriately labeled with sample name.
7. Apply up to 50 µl sample gently to the G-50 bed of the column.
8. Spin for 2 minutes at 3 K rpm.
9. Purified sample is collected in the collection tube. The eluent collected in the 1.5 ml tube is free of salts and primers shorter than 35-40mer.

References

1. Kovarova, M; and Draber, P; (2000) New Specificity and yield enhancer for polymerase chain reactions (2000) Nucl. Acids. Res. 28: e70.
2. Henke, W., Herdel, K., Jung, K., Schnorr, D. and Stefan A. Loening, S. (1997) Betaine improves the PCR amplification of GC-rich DNA sequences. Nucl. Acids Res. 25: 3957-3958.
3. Daniel S. Mytelka, D.S., and Chamberlin, M.J.,(1996) Analysis and suppression of DNA polymerase pauses associated with a trinucleotide consensus. Nuc. Acids Res., 24:2774–278.
4. Keith, J. M., Cochran, D.A.E., Lala, G.H., Adams, P., Bryant, D.and Mitchelson, K.R. (2004) Unlocking hidden genomic sequence. Nucl. Acids Res. 32: e35.
5. Owczarzy, R., Dunietz, I., Behlke, M.A., Klotz, I.M. and Joseph A. Walder (2003) Thermodynamic treatment of oligonucleotide duplex–simplex equilibria. PNAS, 100:14840-14845.

GeneProber™ Product Ordering Information

The GeneProber™ product line is based on the chemiluminescent Southern blot detection method. Gene Link's non-radioactive detection systems for genotyping of triple repeat disorders are rapid, reliable and as sensitive as the ³²P labeled southern blots. No more decayed probes and radioactive exposure. Kits are available for reliable genotyping of the fragile X, myotonic dystrophy and other triple repeat mutation group disorders.

Unlabeled GeneProber™ probes are also available for radio labeling and radioactive based detection. Gene Link strongly recommends the use of non-radioactive gene detection systems. Consider switching to Gene Link's product line of non-radioactive detection systems

| Product | Unit Size | Catalog No. |
|---|-----------|-------------|
| Fragile X GeneProber™ GLFX1 Probe unlabeled | 500 ng | 40-2004-40 |
| Fragile X GeneProber™ GLFXDig1 Probe Digoxigenin labeled | 110 µL | 40-2004-41 |
| Huntington's Disease GeneProber™ GLHD14 Probe unlabeled | 500 ng | 40-2025-40 |
| Huntington's Disease GeneProber™ GLHDDig2X Probe Digoxigenin labeled | 110 µL | 40-2025-41 |
| Myotonic Dystrophy GeneProber™ GLDM1 Probe unlabeled | 500 ng | 40-2026-40 |
| Myotonic Dystrophy GeneProber™ GLDMDig2 Probe Digoxigenin labeled | 110 µL | 40-2026-41 |
| Friedreich's Ataxia GeneProber™ GLFRDA21 Probe unlabeled | 500 ng | 40-2027-40 |
| Friedreich's Ataxia GeneProber™ GLFRDADig21 Probe Digoxigenin labeled | 110 µL | 40-2027-41 |

GScan™ Products Product Ordering Information

Gene Link's GScan™ gene detection products are safe, convenient and sensitive, and afford automated compilation of data. The kits contain optimized PCR amplification reagents and a wide array of fluorescent-labeled primers for genotyping after PCR using fluorescent genetic analyzer instrument(s). Included in these kits are ready-to-run control samples of various repeats of the triple repeat disorder kit. These control samples are for calibration with the molecular weight markers for accurate size determination of the amplified fragments.

The GScan™ kits are simple and robust for routine triple-repeat detection of greater than 100 repeats of all triple repeat disorders listed, except Fragile X. The CGG repeat in Fragile X can be detected up to ~50 repeats.

| Product | Unit Size | Catalog No. |
|--|-----------|---------------|
| Fragile X GScan™ Kit for fluorescent detection; 100 reactions kit | 1 kit | 40-2004-15XX |
| Fragile X GScan™ Kit for fluorescent detection; 20 reactions kit | 1 kit | 40-2004-15FMS |
| Huntington's Disease GScan™ Kit for fluorescent detection; 100 reactions kit | 1 kit | 40-2025-15XX |
| Huntington's Disease GScan™ Kit for fluorescent detection; 20 reactions kit | 1 kit | 40-2025-15FMS |
| Myotonic Dystrophy GScan™ Kit for fluorescent detection; 100 reactions kit | 1 kit | 40-2026-15XX |
| Myotonic Dystrophy GScan™ Kit for fluorescent detection; 20 reactions kit | 1 kit | 40-2026-15FMS |
| Friedreich's Ataxia GScan™ Kit for fluorescent detection; 100 reactions kit | 1 kit | 40-2027-15XX |
| Friedreich's Ataxia GScan™ Kit for fluorescent detection; 20 reactions kit | 1 kit | 40-2027-15FMS |

All Gene Link products are for research use only

Current pricing are posted at <http://www.genelink.com/>

Huntington's Disease Product Ordering Information

| Product | Unit Size | Catalog No. |
|--|----------------------|--------------|
| Huntington's Disease Genemer™ Primer pair Primers for amplification of CAG triple repeat spanning region. The quantity supplied is sufficient for 400 regular 50 µl PCR reactions. | 10 nmols | 40-2025-10 |
| Huntington's Disease PCRProber™ AP labeled probe Alkaline phosphatase labeled probe | 12 µL | 40-2025-31 |
| Huntington's Disease PCRProber™ Kit for chemiluminescent detection Kit for performing PCR amplification and chemiluminescent based detection. | 5 blots [50 rxns] | 40-2025-32 |
| Huntington's Disease Genemer™ Kit for Radioactive Detection Kit for amplification and radioactive detection of Huntington's Disease CAG triple repeat region amplified PCR products using ³⁵ S or ³² P. 100 Reactions. | 1 Kit [100 rxns] | 40-2025-20 |
| Huntington's Disease GScan Kit for fluorescent detection Kit for performing fluorescent PCR amplification based detection. Various dye kits. XX=FM for 6-Fam; HX for Hex; TT for Tet; C3 for Cy3 and C5 for Cy5. | 1 Kit [100 rxns] | 40-2025-15XX |

Genemer™ GScan Control DNA Cloned fragment of the mutation region of a particular gene. These control DNAs are ideal genotyping templates for optimizing and performing control amplification with unknown DNA. The size of the triple repeats has been determined by sequencing and gel electrophoresis. The stability of size repeats upon cloning and amplification has NOT been determined. Thus, the size should be considered approximate and there is no claim for each fragment to contain the exact number of triple repeats. These control DNAs are sold with the express condition that these NOT be used for exact triple repeat size determination of DNA of unknown genotype. The control DNA should be used for determining the performance of specific Genemer™ and PCRProber™ Gene Link products.

| | | |
|--|-------|--------------|
| Huntington's Disease 7 ~CAG repeat GScan Genemer Control DNA; HEX labeled | 25 µL | 40-2025-05HX |
| Huntington's Disease 18 ~CAG repeat GScan Genemer Control DNA; HEX labeled | 25 µL | 40-2025-01HX |
| Huntington's Disease 31 ~CAG repeat GScan Genemer Control DNA; HEX labeled | 25 µL | 40-2025-07HX |
| Huntington's Disease 34 ~CAG repeat GScan Genemer Control DNA; HEX labeled | 25 µL | 40-2025-02HX |
| Huntington's Disease 37 ~CAG repeat GScan Genemer Control DNA; HEX labeled | 25 µL | 40-2025-08HX |
| Huntington's Disease 44 ~CAG repeat GScan Genemer Control DNA | 25 µL | 40-2025-03HX |
| Huntington's Disease 49 ~CAG repeat GScan Genemer Control DNA | 25 µL | 40-2025-09HX |
| Huntington's Disease 89 ~CAG repeat GScan Genemer Control DNA | 25 µL | 40-2025-04HX |
| Huntington's Disease 116 ~CAG repeat GScan Genemer Control DNA | 25 µL | 40-2025-06HX |
| Huntington's Disease 134 ~CAG repeat GScan Genemer Control DNA | 25 µL | 40-2025-61HX |
| Huntington's Disease 182 ~CAG repeat GScan Genemer Control DNA | 25 µL | 40-2025-62HX |

All Gene Link products are for research use only

Current pricing are posted at <http://www.genelink.com/>

Related Products Ordering Information

Omni-Pure™ DNA & RNA Purification Systems

| Product | Catalog No. | Unit Size*(Purifications) |
|--|-------------|---------------------------|
| Omni-Pure™ Blood DNA Purification System | 40-4010-01 | 100 |
| Omni-Pure™ Blood DNA Purification System | 40-4010-05 | 500 |
| Omni-Pure™ Blood DNA Purification System | 40-4010-10 | 1000 |
| Omni-Pure™ Tissue DNA Purification System | 40-4050-01 | 100 |
| Omni-Pure™ Tissue DNA Purification System | 40-4050-05 | 500 |
| Omni-Pure™ Tissue DNA Purification System | 40-4050-10 | 1000 |
| Omni-Pure™ Plant DNA Purification System | 40-4060-01 | 100 |
| Omni-Pure™ Plant DNA Purification System | 40-4060-05 | 500 |
| Omni-Pure™ Plant DNA Purification System | 40-4060-10 | 1000 |
| Omni-Pure™ Viral DNA Purification System | 40-3720-01 | 100 |
| Omni-Pure™ Viral DNA Purification System | 40-3720-05 | 500 |
| Omni-Pure™ Microbial DNA Purification System | 40-3700-01 | 100 |
| Omni-Pure™ Microbial DNA Purification System | 40-3700-05 | 500 |
| Omni-Pure™ Viral RNA Purification System | 40-3650-01 | 100 |
| Omni-Pure™ Viral RNA Purification System | 40-3650-05 | 500 |

*Sample volume for each purification system varies. Each purification yields sufficient quantity for desired applications.

Omni-Clean™ Gel DNA Purification and Concentration Systems

| Product | Catalog No. | Unit Size*(Purifications) |
|---|-------------|---------------------------|
| Omni-Clean™ Gel DNA Beads Purification System | 40-4110-10 | 100 |
| Omni-Clean™ Gel DNA Beads Purification System | 40-4110-50 | 500 |
| Omni-Clean™ Gel DNA Spin Column Purification System | 40-4120-10 | 100 |
| Omni-Clean™ Gel DNA Spin Column Purification System | 40-4120-50 | 500 |
| Omni-Clean™ DNA Beads Concentration System | 40-4130-10 | 100 |
| Omni-Clean™ DNA Beads Concentration System | 40-4130-50 | 500 |
| Omni-Clean™ DNA Spin Column Concentration System | 40-4140-10 | 100 |
| Omni-Clean™ DNA Spin Column Concentration System | 40-4140-50 | 500 |

*Sample volume for each purification system varies. Each purification yields sufficient quantity for desired applications.

Omni-Pure™ Plasmid DNA Purification Systems

| Product | Catalog No. | Unit Size*(Purifications) |
|--|-------------|---------------------------|
| Omni-Pure™ Plasmid DNA Purification System | 40-4020-01 | 100 |
| Omni-Pure™ Plasmid DNA Purification System | 40-4020-05 | 500 |

*Sample volume for each purification system varies. Each purification yields sufficient quantity for desired applications.

All Gene Link products are for research use only

Current pricing are posted at <http://www.genelink.com/>

Southern Blot Buffers & Reagents

| Product | Catalog No. | Unit Size |
|--|-------------|-------------|
| Agarose Tablets, 0.5 gm each; 100 tablets | 40-3011-10 | 100 tablets |
| Agarose LE Molecular Biology Grade; 100 g | 40-3010-10 | 100 g |
| Agarose LE Molecular Biology Grade; 500 g | 40-3010-50 | 500 g |
| Hybwash A, Hybridization Wash Solution (20X SSC); 200 mL | 40-5020-20 | 200 mL |
| Hybwash B, Hybridization Wash Solution (10% SDS); 100 mL | 40-5021-10 | 100 mL |
| TAE Buffer; 50 X Concentrate; 100 mL | 40-3007-01 | 100 mL |
| TAE Buffer; 50 X Concentrate; 1 L | 40-3007-10 | 1 L |
| TBE Buffer; 5 X Concentrate; 1 L | 40-3008-10 | 1 L |
| Buffer M 10X (Maleic Acid buffer); 100 mL | 40-5025-10 | 100 mL |
| 10% Blocking solution; 100 mL | 40-5026-10 | 100 mL |
| Loading Buffer 2X BPB/XC Denaturing for Sequencing; 1 mL | 40-5027-10 | 1 mL |
| 10x AP Detection buffer (alkaline phosphatase detection buffer); 100 mL | 40-5031-10 | 100 mL |
| Lumisol™ I Hybridization Solution; contains formamide; 200 mL | 40-5022-20 | 200 mL |
| Lumisol™ II Hybridization Solution; for non-toxic hybridizations; 200 mL | 40-5023-20 | 200 mL |
| Lumisol™ III Hybridization Solution; for oligo probes; 200 mL | 40-5024-20 | 200 mL |
| CDP-Star® Substrate; Ready-to-Use 0.25 mM in spray bottle; 10 mL | 40-5010-10 | 10 mL |

Loading Buffers

| Product | Catalog No. | Size |
|---|-------------|-------|
| Gel Loading Buffer 5X BPB/XC non-denaturing; 1 mL | 40-3002-10 | 1 mL |
| Gel Loading Buffer 5X BPB/XC non-denaturing; 15 mL | 40-3002-15 | 15 mL |
| Gel Loading Buffer 10X BPB/XC non-denaturing; 1 mL | 40-3003-10 | 1 mL |
| Gel Loading Buffer 10X BPB/XC non-denaturing; 15 mL | 40-3003-15 | 15 mL |
| Gel Loading Buffer 5X Orange G/XC non-denaturing; 1 mL | 40-3004-10 | 1 mL |
| Gel Loading Buffer 5X Orange G/XC non-denaturing; 15 mL | 40-3004-15 | 15 mL |
| Gel Loading Buffer 2X BPB/XC Denaturing for Sequencing; 1 mL | 40-5027-10 | 1 mL |
| Gel Loading Buffer 2X BPB/XC Denaturing for Sequencing; 15 mL | 40-5027-15 | 15 mL |
| DNA SDS Gel Loading Buffer 5X BPB/XC DNA binding protein denaturing buffer ; 1 mL | 40-5028-10 | 1 mL |
| DNA SDS Gel Loading Buffer 5X BPB/XC DNA binding protein denaturing buffer; 15 mL | 40-5028-15 | 15 mL |
| RNA Gel Loading Buffer 2X BPB/XC with ethidium bromide; 1 mL | 40-5029-10 | 1 mL |
| RNA Gel Loading Buffer 2X BPB/XC with ethidium bromide; 15 mL | 40-5029-15 | 15 mL |
| RNA Gel Loading Buffer 2X BPB/XC without ethidium bromide ; 1 mL | 40-5030-10 | 1 mL |
| RNA Gel Loading Buffer 2X BPB/XC without ethidium bromide; 15 mL | 40-5030-15 | 15 mL |

Omni-Marker™

| Product | Catalog No. | Size* |
|--|-------------|--------|
| Omni-Marker™ Universal unlabeled; 1 mL | 40-3005-10 | 1 mL |
| Omni-Marker™ Low unlabeled; 1 mL | 40-3006-10 | 1 mL |
| Omni-Marker™ GScan™-2 Tamra labeled 50 bp - 600 bp; 500 µL | 40-3062-05 | 500 µL |

All Gene Link products are for research use only

Current pricing are posted at <http://www.genelink.com/>

