



Product Specifications

DNA Molecular Weight Markers, Electrophoresis Reagents
Polymerase Chain Reaction Reagents
Custom Primers and Probes
Hybridization and Detection Reagents

DNA & RNA Gel Loading Buffers

	Catalog No.	Product	Size
<input type="checkbox"/>	40-3002-10	Gel Loading Buffer 5X BPB/XC non-denaturing	1 ml
<input type="checkbox"/>	40-3002-15	Gel Loading Buffer 5X BPB/XC non-denaturing	15 ml
<input type="checkbox"/>	40-3003-10	Gel Loading Buffer 10X BPB/XC non-denaturing	1 ml
<input type="checkbox"/>	40-3003-15	Gel Loading Buffer 10X BPB/XC non-denaturing	15 ml

Storage

Store at room temperature or at +4°C up to 12 months.
For longer periods, store at -20°C.

Product Description & Application

Gel Loading Buffer 5X BPB/XC and 10X BPB/XC are concentrates that contain bromophenol blue and xylene cyanol as tracking dyes. The migration of the tracking dyes differs in different percentage gels. See table below for approximate co-migration of dyes with DNA base pair sizes. In agarose gels over a range of 0.7% to 1.4% (w/v), xylene cyanol co-migrates with ~4 kb linear double strand DNA and bromophenol blue co-migrates with ~300 bp double strand DNA. In a 12% polyacrylamide gel xylene cyanol co-migrates with approximately 70 bp DNA and bromophenol blue co-migrates with ~20 bp DNA.

Suitable for use in nucleic acid gel electrophoresis
DNase, RNase: None detected

Application

Gel Loading Buffer 5X BPB/XC non-denaturing

Add 1 volume of the Gel Loading Solution to 4-5 volumes of sample, mix well and load to gel.

Gel Loading Buffer 10X BPB/XC non-denaturing

Add 1 volume of the Gel Loading Solution to 9-10 volumes of sample, mix well and load to gel.

10X Composition

0.25%(w/v) Bromophenol blue
0.25%(w/v) Xylene cyanole FF
15%(w/v) Ficoll in water

Reference

1. Sambrook, J., *et al.*, Molecular Cloning, A Laboratory Manual, Cold Spring Harbor Laboratory, p. 6.13-6.62 (1989)

Recommended Gel Percentages for Separation of Linear DNA*			
Agarose gel, %	Range of separation, bp	Polyacrylamide gel, %	Range of separation, bp
0.5	1,000-30,000	3.5	100-1,000
0.7	800-12,000	5.0	80-500
1.0	500-10,000	8.0	60-400
1.2	400-7,000	12.0	40-200
1.4	200-4,000	20.0	5-100
2.0	50-2,000		

* Approximate electrophoretic resolution of DNA fragments sizes.

Migration Rates of the Marker Dyes in Polyacrylamide Gels*		
Polyacrylamide gel, %	Bromophenol blue*	Xylene cyanol FF*
Non-denaturing gels		
3.5	100 bp	460 bp
5.0	65 bp	260 bp
8.0	50 bp	160 bp
12.0	45 bp	70 bp
15.0	20 bp	60 bp
20.0	12 bp	45 bp
Denaturing gels		
5.0	35 bases	130 bases
6.0	29 bases	106 bases
8.0	26 bases	76 bases
10.0	12 bases	55 bases
20.0	8 bases	28 bases

* The approximate sizes of DNA fragments with which the indicated marker dye co-migrates.

Migration Rates of the Marker Dyes in Agarose Gels			
Agarose concentration, %	Xylene cyanol FF	Bromophenol blue	Orange G
0.7-1.7	~4000bp	~300bp	~50bp
2.5-3.0	~800bp	~100bp	~30bp

* The approximate sizes of DNA fragments with which the indicated marker dye co-migrates.

Reference

Sambrook, J., et al., Molecular Cloning. A Laboratory Manual, Cold Spring Harbor Laboratory, Cold Spring Harbor, N.Y., 12.89, 5.42, 2001

Ordering Information

Loading Buffers

Product	Catalog No.	Size	Price \$
Gel Loading Buffer 5X BPB/XC non-denaturing	40-3002-10	1 ml	10.00
Gel Loading Buffer 5X BPB/XC non-denaturing	40-3002-15	15 ml	80.00
Gel Loading Buffer 10X BPB/XC non-denaturing	40-3003-10	1 ml	16.00
Gel Loading Buffer 10X BPB/XC non-denaturing	40-3003-15	15 ml	95.00
Gel Loading Buffer 5X Orange G/XC non-denaturing	40-3004-10	1 ml	10.00
Gel Loading Buffer 5X Orange G/XC non-denaturing	40-3004-15	15 ml	80.00
Gel Loading Buffer 2X BPB/XC Denaturing for Sequencing	40-5027-10	1 ml	10.00
Gel Loading Buffer 2X BPB/XC Denaturing for Sequencing	40-5027-15	15 ml	80.00
DNA SDS Gel Loading Buffer 5X BPB/XC DNA binding protein denaturing buffer	40-5028-10	1 ml	10.00
DNA SDS Gel Loading Buffer 5X BPB/XC DNA binding protein denaturing buffer	40-5028-15	15 ml	60.00
RNA Gel Loading Buffer 2X BPB/XC with ethidium bromide	40-5029-10	1 ml	36.00
RNA Gel Loading Buffer 2X BPB/XC with ethidium bromide	40-5029-15	15 ml	82.00
RNA Gel Loading Buffer 2X BPB/XC without ethidium bromide	40-5030-10	1 ml	26.00
RNA Gel Loading Buffer 2X BPB/XC without ethidium bromide	40-5030-15	15 ml	72.00

Related Products Ordering Information

Buffers & Reagents

Product	Catalog No.	Size	Price \$
Agarose LE Molecular Biology Grade; 100 gms	40-3010-10	100 gms	120.00
Agarose LE Molecular Biology Grade; 500 gms	40-3010-50	500 gms	410.00
Agarose Tablets, 0.5 gm each	40-3011-10	100 tablets	100.00
TAE Buffer; 50 X Concentrate	40-3007-01	100 ml	32.00
TBE Buffer; 5 X Concentrate	40-3008-10	1000 ml	35.00
Hybwash A, Hybridization Wash Solution	40-5020-20	200 ml	65.00
Hybwash B, Hybridization Wash Solution	40-5021-10	100 ml	50.00
10x Washing buffer	40-5025-20	200 ml	125.00
10% Blocking solution	40-5026-10	100 ml	75.00
Seq. Loading buffer	40-5027-00	1 ml	10.00
10x AP Detection buffer	40-5031-10	100 ml	65.00
Lumisol™ I Hybridization Solution; contains formamide	40-5022-20	200 ml	75.00
Lumisol™ II Hybridization Solution; for non-toxic hybridizations	40-5023-20	200 ml	75.00
Lumisol™ III Hybridization Solution; for oligo probes	40-5024-20	200 ml	75.00

Omni-Marker™

Product	Catalog No.	Size*	Price \$
Omni-Marker™ Universal unlabeled	40-3005-10	1 ml	90.00
Omni-Marker™ Low unlabeled	40-3006-10	1 ml	90.00
Omni-Marker™ GScan-2 Tamra labeled 50 bp - 600 bp	40-3062-05	500 µl	325.00

Prices subject to change without notice.

All Gene Link products are for research use only