

### Custom Oligo Specifications

Gene Link custom oligonucleotides are supplied desalted and lyophilized. They are ready to use after appropriate reconstitution. Dry oligonucleotides are stable at room temperature for an extended period of time.

#### Storage & Reconstitution

The oligonucleotide should preferably be frozen upon receipt. TE buffer (10mM Tris, 1mM EDTA, pH 7.5) is recommended for dissolving the oligonucleotides; EDTA inhibits the activity of the nucleases. Further dilution can be made in distilled sterile water. After reconstitution store the stock solution at -80°C or -20°C.

#### Gel Photo Documentation

An actual gel picture of the synthesized custom oligonucleotide is supplied. Polyacrylamide gels of 12 to 15% are run, depending upon the length of the custom oligonucleotide. A major single band represents high purity of the crude oligonucleotide.

#### Purity & Usage

The crude, desalted oligonucleotide supplied is suitable for all amplification and sequencing protocols. Gel purification is advised for all oligos used for cloning applications and for oligos longer than 50mer.

#### Biophysical Data

Each oligo after desalting is quantified by recording  $A_{260}$ . Exact nmols and  $\mu\text{g}$  is determined by the extinction coefficient and molecular weight of the oligo.

### Oligo Scale of Synthesis and Typical Yield

Scale	Crude Desalted			RPC Purified**			Gel Purified		
	20mer oligo* Typical yield			30mer oligo* Typical yield			50mer oligo* Typical yield		
	A <sub>260</sub> Units	nmols	mg	A <sub>260</sub> Units	nmols	mg	A <sub>260</sub> Units	nmols	mg
50 nmol	8-10	30+	0.2-0.3	4-5	12+	0.1-0.16	NR* [1-2]	NR* [2-4]	NR* [0.03-0.06]
200 nmol	20-25	80+	0.6-0.8	8-12	24+	0.26-0.40	4-6	8+	0.13-0.2
1 $\mu\text{mol}$	100-120	400+	3-4	40-50	30+	1.3-1.6	20-25	40+	0.6-0.8
<b>Purity &amp; Yield</b>	Purity is more than 80% depending on oligo sequence and structure. Refer to coupling efficiency table for oligo length dependent purity and yield.  No further purification required for PCR and sequencing applications.  Gel purification recommended for oligos above 50mer and all applications involving cloning and mutagenesis.			Purity 85% to 95% depending on oligo sequence and structure.  Yield and purity will be lower for sequences with high GC content.  Not recommended for oligos longer than 35mer.  **RPC is reverse phase purification using a cartridge; a substitute for HPLC.			Purity 98% to ~100% depending on oligo sequence and structure.  Yield will gradually decrease as length of oligo increases. Palindromes, hairpins and high GC content oligos and oligos containing stretches of 3 or more G's induces strong secondary structure and base stacking thus decreasing purity and yield.  NR* Not Recommended		
*Yield of 30 $\mu\text{g}/A_{260}$ unit for oligos is calculated for an ~equimolar base composition. Long stretches of a single base or homopolymers will have variable yields. Example for homopolymeric 50mer: A(50) = ~20/ $A_{260}$ Unit; G(50) = ~28/ $A_{260}$ Unit; T(50) = ~35/ $A_{260}$ Unit and C(50) = ~39/ $A_{260}$ Unit.									

### Oligo Reconstitution and Use

Gene Link oligos are supplied lyophilized. These are stable at room temperature for an extended period of time. TE buffer (10mM Tris, 1mM EDTA, pH 7.5) is recommended for dissolving the oligonucleotides; EDTA inhibits the activity of the nucleases. Further dilution can be made in distilled sterile water. After reconstitution store the stock solution at -80°C or -20°C.

Standard PCR Set Up		
Reagent	Final Concentration	Quantity/ 50 $\mu\text{l}$ Reaction
Sterile deionized water	-	variable
10X *PCR buffer	1X	5 $\mu\text{l}$
2mM dNTP mix	0.2mM of each	5 $\mu\text{l}$
Primer I, 10 $\mu\text{M}$ (10pmol/ $\mu\text{l}$ )	0.5 $\mu\text{M}$	2.5 $\mu\text{l}$
Primer II, 10 $\mu\text{M}$ (10pmol/ $\mu\text{l}$ )	0.5 $\mu\text{M}$	2.5 $\mu\text{l}$
Taq DNA Polymerase, 5U/ $\mu\text{l}$	1.25u/50 $\mu\text{l}$	0.25 $\mu\text{l}$
Template DNA	10pg-1 $\mu\text{g}$	variable
*Final MgCl <sub>2</sub> concentration is 1.5mM		

#### Oligo Reconstitution

##### Stock solution of 500pmols/ $\mu\text{l}$ [500 $\mu\text{M}$ ]

Gene Link provides the exact amount of nmols of each oligo supplied on the tube and on the Oligo Report. Multiply the 'nmol' amount by 2 to arrive at the volume of TE to be added.

**Example:** 45.10 nmols  $\times$  2 = 90.2 $\mu\text{l}$

Dissolve the oligo in 90.2 $\mu\text{l}$  to get 500pmols/ $\mu\text{l}$  stock solution. Use as required.

Dilute 10 fold to prepare a 50pmols/ $\mu\text{l}$  [50 $\mu\text{M}$ ]. Use as required.

##### Stock solution of 100 pmols/ $\mu\text{l}$ [100 $\mu\text{M}$ ]

Gene Link provides the exact amount of nmols of each oligo supplied on the tube and on the Oligo Report. Multiply the 'nmol' amount by 10 to arrive at the volume of TE to be added.

**Example:** 45.10 nmols  $\times$  10 = 451 $\mu\text{l}$

Dissolve the oligo in 451 $\mu\text{l}$  to get 100pmols/ $\mu\text{l}$  stock solution. Use as required.

Dilute 10 fold to prepare a 10pmols/ $\mu\text{l}$  [10 $\mu\text{M}$ ]. Use as required.

#### Examples of Use

##### Polymerase Chain Reaction (PCR)

The final concentration of primers in a PCR reaction is 0.2–1.0 $\mu\text{M}$ . This is equivalent to 0.2–1 pmol/ $\mu\text{l}$ . At Gene Link, for a standard PCR we use 0.5 pmol/ $\mu\text{l}$ .

##### Sequencing

The final concentration of primer in automated sequencing is from 4 to 10 pmols (~0.05 - 0.1 $\mu\text{g}$ ). Use the oligo reconstitution protocol to prepare a 100 pmols/ $\mu\text{l}$  [100 $\mu\text{M}$ ] solution and then dilute 10 fold to get 10 pmol/ $\mu\text{l}$  solution. Use 1 $\mu\text{l}$  (10 pmols).

##### Quick Conversion Table

1 $\mu\text{M}$  ( $\mu\text{Molar}$ ) = 1 pmol/ $\mu\text{l}$  (picomoles/ $\mu\text{l}$ )  
1 mM (milliMolar) = 1 nmols/ $\mu\text{l}$  (nanomoles/ $\mu\text{l}$ )

**Example:** 20 $\mu\text{Molar}$  primer solution is 20 pmol/ $\mu\text{l}$