

### Custom Oligo Specifications

Gene Link custom oligonucleotides are supplied desalted and lyophilized. They are ready to use after appropriate reconstitution. Dry oligonucleotides are stable at room temperature for an extended period of time.

#### Storage & Reconstitution

The oligonucleotide should preferably be frozen upon receipt. TE buffer (10mM Tris, 1mM EDTA, pH 7.5) is recommended for dissolving the oligonucleotides. After reconstitution store the stock solution at -80°C or -20°C.

#### Gel Photo Documentation

An actual gel picture of the synthesized custom oligonucleotide is supplied. A major single band represents high purity of the crude oligonucleotide.

#### Purity & Usage

The crude, desalted oligonucleotide supplied is suitable for all amplification and sequencing protocols. Gel purification is advised for all oligos used for cloning applications and for oligos longer than 50mer.

#### Biophysical Data

Each oligo after desalting is quantified by recording  $A_{260}$ . Exact nmols and  $\mu\text{g}$  is determined by the extinction coefficient and molecular weight of the oligo.

### Oligo Scale of Synthesis and Typical Yield of Unmodified Oligos\*

Scale	Crude Desalted		RPC Purified***		Gel Purified	
	20mer oligo**		30mer oligo**		50mer oligo**	
	$A_{260}$ Units	nmols	$A_{260}$ Units	nmols	$A_{260}$ Units	nmols
50 nmol	8-10	30+	4-5	12+	NR* [1-2]	NR* [2-4]
200 nmol	20-25	80+	8-12	24+	4-6	8+
1 $\mu\text{mol}$	100-120	400+	40-50	30+	20-25	40+
Purity & Yield	Purity is more than 80% depending on oligo sequence and structure.		Purity 85% to 95% depending on oligo sequence and structure. Not recommended for oligos longer than 35mer.		Purity 98% to ~100% depending on oligo sequence and structure. Yield will gradually decrease as length of oligo increases.	

\*The yield of modified oligos varies based on modification.

\*\*Yield of 30  $\mu\text{g}/A_{260}$  unit for oligos is calculated for an ~equimolar base composition. Long stretches of a single base or homopolymers will have variable yields.

Example for homopolymeric 50mer: A(50) = ~20/ $A_{260}$  Unit; G(50) = ~28/ $A_{260}$  Unit; T(50) = ~35/ $A_{260}$  Unit and C(50) = ~39/ $A_{260}$  Unit.

\*\*\*RPC is reverse phase purification using a cartridge; a substitute for HPLC.

NR\* Not Recommended.

### Fluorescent Molecular Primers & Probes

The use of fluorescent dyes in molecular biology has rapidly transformed from just single dye labeled primers for fragment analysis to the use of multiple labeled dyes and quenchers as probes for real time quantitative PCR analysis. Fluorescence based detection offers a safe and sensitive method for quantitative detection. Gene Link offers synthesis of all different forms of molecular primers and probes. We provide technical service in the design of novel probes and synthesize numerous combinations of dyes, quenchers, RNA, phosphorothioate, 2' O methyl and chimeric probes.

#### Excitation and Emission

The excitation level of molecules varies at different wavelengths. Molecules exposed to a beam of light absorb more at a particular wavelength. This specific

wavelength is termed as the Excitation Maxima. The emission maximum is the wavelength at which the maximum amount of light is released. The molecule stays in the excited state for a finite time, usually <1-10 nanoseconds and returns to the relaxed state upon emission of energy. Excitation and Emission is a cyclic process and consequently can be repeated to an extent before it starts to fade, termed as photo-bleaching.

Different fluorescent dyes are used for molecular probes and primer design. The dyes are selected based on the excitation and emission wavelengths, bleaching, quenching and various other biophysical factors.

#### Quenching

Reduction in the expected fluorescence emission is termed as quenching. The phenomenon of quenching forms the basis of the mode of action of molecular probes; the designed and controlled fluorescence based on hybridization to the target sequence.

#### Quenchers

Dye	Absorbance max (nm)
Dabcyl	453
*BHQ-1	534
*BHQ-2	579
*BHQ-3	672

Placing a molecule that absorbs light in close proximity to the fluorophore can induce quenching. The quenching effect is exhibited by fluorescent as well as non-fluorescent molecules. A non-fluorescent quencher

is the basis of the design of Molecular Beacons.

#### Fluorescence Resonance Energy Transfer (FRET)

Resonance energy transfer, often known as fluorescence resonance energy transfer (FRET) or Förster energy transfer. It is the radiation-less transfer of excitation energy from a donor to an acceptor. An important consequence of this transfer is that there is no emission of light by the donor. The acceptor may or may not be fluorescent.

FRET varies based on the degree of spectral overlap of the donor and acceptor. This is called the "spectral overlap" or sometimes the "Förster overlap integral". This describes the amount of overlap where resonance can occur, i.e. where the donor and acceptor have the same frequencies.

#### TaqMan Probes

TaqMan (also known as Fluorogenic 5' nuclease assay) probes contain two dyes, a reporter dye (e.g. 6-FAM) at the 5' end and a 3' acceptor dye, usually TAMRA. Recent designs substitute the 3' TAMRA fluorescent acceptor quencher dye with non-fluorescent quencher, e.g. BHQ-1. The proximity of the quencher to the reporter in an intact probe allows the quencher to suppress, or "quench" the fluorescence signal of the reporter dye through FRET.

#### Molecular Beacons

Molecular beacons are hairpin shaped oligos with a fluorophore and a quencher at either ends. The loop serves as the specific target sequence. The stem is formed by the annealing of complementary arm

### Fluorophore Absorbance and Emission Data\*

Dye	Color	Absorbance Max (nm)	Emission Max (nm)	Extinction Coefficient
6-FAM (Fluorescein)	Green	494	525	74850
TET	Orange	521	536	85553
HEX	Pink	535	556	95698
Cy 5	Violet	646	667	250000
Cy 5.5	Blue	683	707	190000
Cy 3	Red	552	570	150000
Cy 3.5	Purple	588	604	150000
Cy 7	Near IR	743	767	200000
Tamra	Rose	565	580	87000
ROX	Purple	587	607	105000
JOE	Mustard	528	554	105000
Alexa Dye Series	Varies	Varies	Varies	Varies

\*Color and fluorescence data vary with pH. Consult appropriate dye manufacturer for details.

sequences on the ends of the probe sequence. The stem keeps the fluorophore and the quencher in close proximity to each other, causing the fluorescence of the fluorophore to be quenched by energy transfer. When the probe encounters a target molecule, it forms a hybrid that is longer and more stable than the stem leading to the restoration of fluorescence.

\*BHQ is a registered trademark of Biosearch Technologies. Complete disclaimer of license statement for Molecular Beacons products, PHRI Molecular Beacon and BHQ license agreement can be viewed at the following link: [www.genelink.com/newsite/products/MBPriceList.asp](http://www.genelink.com/newsite/products/MBPriceList.asp)

