



Premium Oligonucleotide Synthesis

QUALITY • CONSISTENCY • CONFIDENCE

Gene Link oligos are for demanding applications and consistent results. We believe that investigators who value time and have no room for an experiment to fail due to oligo quality should consider Gene Link. Our numerous quality control steps for each oligo assure confidence.



GOLD STANDARD

Actual Gel Photo

An actual gel photo of each oligo is affixed on the oligo report. An absolute testimony of quality.

Gene Link has raised the standard since inception over a decade ago.

We have the pictures to prove it!

Superior to “Mass-Produced Factory Oligos”

Gene Link is not an oligo factory. Each oligo is synthesized, processed and quality assured to Gene Link’s absolute standards. This includes coupling efficiency monitoring of each base during synthesis and electrophoretic analysis of each oligo on a polyacrylamide gel to visually assess quality.



Coupling Efficiency

We maintain a coupling efficiency threshold of greater than 99.5% for all oligos by using premium reagents of exacting specifications, membrane synthesis, state-of-the-art instruments and optimized software-driven protocols. This may not be evident when comparing short oligos, as PCR and sequencing reactions are very robust and can tolerate up to 50% failure/truncated sequence oligos. However, you are clearly taking a chance by using long oligos synthesized at anything below 99.5% coupling efficiency.

Trityl Monitoring

All Gene Link DNA synthesizers are equipped with trityl monitors for monitoring coupling efficiency of each added base. The instruments are programmed to halt when it falls below the threshold.

See example of routine trityl bars.



Actual trityl coupling efficiency of a 210 mer.

The Gene Link Advantage

- Stringent Quality Control Measures
- Specializing in Long Oligos up to 250 mer
- Trityl Monitoring of All Oligos
- Polyacrylamide Gel Photograph of Each Oligo
- All Modifications Available
- All Oligo Types Available
- Easy Online Ordering System
- Online Design and Analysis Tools
- Knowledgeable Technical Support
- Personalized, Friendly Customer Service



Long Oligos

Ask our competitors how often they synthesize 200 to 250 mer oligonucleotides. Gene Link specializes in long oligos.

You are invited to compare.

Gene Link has been leading the way by providing the most user friendly online experience in oligo ordering. From oligo design and analysis, to the convenient ordering system and the assurance of a secured transaction, Gene Link provides the most comprehensive web resource in the industry.

Features include convenient NCBI blasting and secondary structure analysis, simple import tools for large orders in spreadsheet or text file format, and three levels of review and editing.

Applications include RNAi Explorer™, a robust siRNA search and design tool, and a standalone Oligo Explorer™ application for online acquisition of sequences and oligo design.

Custom Oligo Ordering System

- Classic ordering system with extensive analysis features
- Timesaver multi oligo import from spreadsheet and text files
- Ability to handle mixed oligo types, purity and modifications in a single order
- Selection of oligo type (DNA, RNA, Phosphorothioate, Chimeric, etc.)
- Simple drop down menu selection for 5', internal or 3' modifications
- Analyze for oligo hairpin and loops
- Integrates with NCBI Blast for homology checks
- Flip 3' to 5' and reverse complement

Online Oligo Analysis

- Simulate annealing, loops and hairpin formations
- Calculate MW, EC, T_m , A_{260} , etc.

SINGLE OLIGO DESIGN — RESULTS —

Click for secondary structure

All thermodynamic values including T_m and secondary structures calculated and displayed solely indicate the relative stability of the secondary structures. They should only be used to compare the relative stability of the structures.

Name : Ex 2871
 Sequence : 5'-AACAGAACGGGTGTACAGCCGCTGACCCGTGA-3'
 Length : 36 nt

T_m (Chemical) : 88.3 °C
 T_m (Local) : 88.3 °C
 T_m (BLI) : 88.3 °C

GC : 50.8 %

ΔG : -76.9 kCal/mol
 ΔH : -289.8 kCal/mol
 ΔS : -905.3 kCal/mol

3'-tail GC % : 51.4 %
 3'-tail ΔG : -11.6 kCal/mol

Secondary structure results are truncated to show the top 4 most stable structures.
[Click here to see all possible associations](#)

Ex 2871 sequence alignment

5'-AACAGAACGGGTGTACAGCCGCTGACCCGTGA-3'
 3'-AATGACCCCGTGCACACATGATGAGGACGAC-5'

Block of 22 bp IF loop

5'-AACAGAACGGGTGTACAGCCGCTGACCCGTGA-3'
 3'-AATGACCCCGTGCACACATGATGAGGACGAC-5'

Save Session

Too busy to order all of your oligos in one session? Gene Link's answer to the multitasking researcher with endless interruptions is the "Save Session" feature. Enter as many oligos as you wish, click the "Save Session" button and resume at your will. Your oligos will be saved. What's more, you'll save money on shipping by consolidating your multiple orders into one.

Gene Link™

CUSTOM OLIGOS ORDER ENTRY

Sequence Number: Ex 2859 (Link to 22 (AACAGAAC))

Name Oligo: Name day dependent if ordered by 22 noon EDT. Extra charges apply.

Quantity: 1000 nmol

Purity: Gel Purified

5' modification: (None)

3' modification: (None)

Internal modification: (None)

Please Note:
 Do not include 5' and 3' positions and do not use 3' only. (A=AT; T=TA; C=GC; G=CG; G=GT; C=CC; A=AA; T=TT; C=CT; G=GG; T=TT; A=AA; C=CC; G=GG; T=TT)
 For designations use T. For all other standard nucleotides use N.
 For oligos longer than 33 bases, 3' end must include 4 bases and gel purification is strongly recommended due to an increasing number of tails as sequence above that length.

SEQUENCE: 5' (sequences in spaces between bases, all CAPITALS, no hard returns): 3'

AACAGAACGGGTGTACAGCC

ANALYZE (more similar)
 NCBI BLAST (see below)
 FLIP 3' TO 5'
 REVERSE COMPLEMENT

Seq: 22 | MW: 46 | MW: 6770.40 | Tr: 56.3 | EC: 326.8 | ΔG_{260} : 95.46 | ΔG_{260} : 5.76

Comments regarding this oligo

SAVE OLIGO NEXT OLIGO | REVIEW SAVED OLIGOS

(The screen will refresh when set up for next oligo) (NOTE: Oligos will be SAVED with the "Save Oligo" button a class before it can occur in the registry list)

Go Back | Home | Support | About Us



Multiple Oligo NCBI Blast

Click [NCBI BLAST](#) to ascertain homologies to other sequences.

Perform NCBI Blast of multiple sequences at once by using Gene Link's online MultiBlast application. Import all the sequences using a spreadsheet or a text file. All of your sequences will be blasted and results retrieved. Gene Link offers a very convenient approach to perform multiple blast searches.





OLIGO EXPLORER 1.2

Oligo Explorer™

A PC-based application for standalone DNA sequence retrieval and oligo design.

Oligo Explorer™ was developed to design PCR and sequencing primers. Oligo Explorer™ is an efficient easy-to-use tool to determine primer properties like T_m , GC%, primer loops and primer dimers.

Oligo Explorer™ also includes a powerful "Primer Wizard" tool that helps you to find suitable primer pairs. You can set your own parameters for the primer pair search engine or use the default parameters. "Primer Wizard" suggests primer pairs that amplify PCR products of the given length. Individual primer pairs are suggested that theoretically will not form stable primer dimers or primer loops.



Molecular Biology Convenience Applets

Gene Link has numerous online applets for quick calculations. The BioCalculator is a series of applets for simplifying the routine laboratory calculations. The following convenient calculators are available:

- Oligo Resuspension
- Oligo Dilution
- Oligo T_m
- Reagent Dilution
- Molarity Determination
- Ligation
- Base/Dye Ratio

Oligo Specifications Report

Gene Link's Custom Oligonucleotide Synthesis Report specifies each oligo name and sequence along with its pertinent physical properties such as MW, %GC, T_m , A_{260} units, etc. Our report is also unique in that we affix an actual polyacrylamide gel electrophoresis photograph onto each report, so that you also may visually attest to the quality of our product.

From your custom oligo to the presentation of our oligo synthesis report, not a step of quality is overlooked. *You are invited to compare.*

Custom Oligo Specifications

Gene Link custom oligonucleotides are supplied desalted and lyophilized. They are ready to use after appropriate reconstitution. Dry oligonucleotides are stable at room temperature for an extended period of time.

Storage & Reconstitution

The oligonucleotide should preferably be frozen upon receipt. TE buffer (10 mM Tris, 1 mM EDTA, pH 7.5) is recommended for dissolving the oligonucleotides. After reconstitution store the stock solution at -80°C or -20°C .

Purity & Usage

The crude, desalted oligonucleotide supplied is suitable for all amplification and sequencing protocols. Gel purification is advised for all oligos used for cloning applications and for oligos longer than 50 mer.

Biophysical Data

Each oligo after desalting is quantified by recording A_{260} . Exact nmols and μg are determined by the extinction coefficient and molecular weight of the oligo.

Gel Photo Documentation

An actual gel picture of the synthesized custom oligonucleotide is supplied. A major single band represents high purity of the crude oligonucleotide.

Customer Name: Alyson Rodgers
Customer Number: 3053287
Order Number: 136038
Date: June 23, 2004

Line	Oligo Name	Sequence (5'-3')	Size	MW	%GC	T_m	nmols	μg	A_{260} Units
1	Primer 1	CAFDCTCAGACACAGCTAGCTCAAGACCTTGGCGGT-CAMTTAGGATACCTAGG	51	16,715	74.8	48.7	785.3	35.61	
2	Primer 2	GGTCTCTTAGTCAAGGAGCTTCCGCGAGTCCCGGTGGGSAFACCTGATCAGCTACTGCTGAGTCA	54	18,710	77.9	47.6	941.8	32.11	
3	Primer 3	CAFDCTCAGACACAGCTAGCTCAAGACCTTGGCGGT-CAMTTAGGATACCTAGG	46	14,800	74.8	45.5	672.7	25.10	
4	Primer 4	CTCAGGACGAAATCGGGAGCGGGCACTTCGTCGGGGCGTTC	40	13,850	77.0	49.0	603.6	21.00	
5	Primer 5	GGGATTCGGTGCACAGGCTTGGTGA	20	7,888	83.8	48.2	275.7	12.70	
6	Primer 6	GGTGTGTGTGGGGTCCCA	10	5,957	58.6	52.0	288.8	0.40	
7	Primer 7	AGAGGAAAGGATAGGAGAC	30	8,287	53.4	40.2	281.8	10.37	
8	Primer 8	CGACCTGCTGTGACAGGCTTCTTGGTGGATGTCATCTCGCGCGGTTCAGTTCGCTGCTGAGGACTGSAK	36	12,869	77.7	5.0	222.7	7.35	
9	Primer 9	TGGTCGAGTCTAGCGCTTCGTCGAGGAAATTTAGGAAAGAGAAAGGCTTCTGATGATTTAGCA	56	20,510	79.2	6.5	252.6	0.14	
10	Primer 10	ATTTCTCAGTACTGTGTTTCAGGAGGAGGATCTTNGATGFGTGGGGTGTACAGACTGAAAGTCAAAAGAGATTTGATAGGATTTTCAAGGAGGATCAAAAGGTTCTCTTTGATTCGAGTTTCTGCTTACTGCTCTTCTGCTGCTTNGGT	101	48,730	79.2	6.2	308.4	10.07	
11	Primer 11	CTCAGGACGAAATCGGGAGCGGGCACTTCGTCGGGGCGTTC	40	13,850	77.0	38.6	388.8	13.00	

Substitution: These units of the label per oligo plate must identify. One can be placed at each end and the middle with a minimum of 15 bases of sequence either side.

In-Situ Hybridization (ISH): Since In-Situ Hybridization (ISH) is an important tool for biological and epidemiological research, contributions of alpha-labeled oligo tags to in-situ hybridization, provides a way to selectively target DNA, RNA transcripts, proteins and other cell components.

Chromatography: Some oligonucleotides at nucleic acid products, and oligonucleotides are also useful for the purification of DNA binding proteins or specific DNA molecules by specific hybridization based affinity chromatography. The hybridized oligonucleotide can be bound to a chromatographic matrix and used for either column or spin chromatography.

Notes: Oligos 1-7 are crude unpurified. Oligos 8-11 are gel purified. Gel lanes for oligos 8-11 correspond to crude followed by gel purified.

Crude Desalted vs. SPE Purified vs. Gel Purified: Comparison table showing %GC, MW, T_m , nmols, and μg for various oligo lengths (10, 20, 30, 40, 50, 60, 70, 80, 90, 100, 150, 200, 300, 400, 500, 600, 700, 800, 900, 1000).

Biophysical Data: Table showing A_{260} units, nmols, and μg for various oligo lengths.

Gel Photo Documentation: Two gel electrophoresis images showing the synthesis results for oligos 1-7 and 8-11. The first gel shows a single band for oligo 1, and the second gel shows multiple bands for oligos 2-11.

Gene Link Logo: Gene Link logo with tagline "Gene Link™ Custom Oligo Synthesis".

Contact Information: 340 Elm Street, Westfield, MA 01095-1101. Phone: 413-754-1192. Fax: 413-754-1193. Email: info@gene-link.com. Website: www.gene-link.com.



Oligo Scale of Synthesis and Typical Yield

Scale	Crude Desalted			RPC Purified**			Gel Purified		
	20 mer oligo* Typical yield			30 mer oligo* Typical yield			50 mer oligo* Typical yield		
	A ₂₆₀ Units	nmols	mg	A ₂₆₀ Units	nmols	mg	A ₂₆₀ Units	nmols	mg
50 nmol	8-10	30+	0.2-0.3	4-5	12+	0.1-0.16	NR* [1-2]	NR* [2-4]	NR* [0.03-0.06]
200 nmol	20-25	80+	0.6-0.8	8-12	24+	0.26-0.4	4-6	8+	0.13-0.2
1 μmol	100-120	400+	3-4	40-50	30+	1.3-1.6	20-25	40+	0.6-0.8
Purity & Yield	Purity is greater than 80% depending on oligo sequence and structure. Refer to coupling efficiency table for oligo length dependent purity and yield. No further purification required for PCR and sequencing applications. Gel purification recommended for oligos above 50 mer and all applications involving cloning and mutagenesis.			Purity 85% to 95% depending on oligo sequence and structure. Yield and purity will be lower for sequences with high GC content. Not recommended for oligos longer than 35 mer. **RPC is reverse phase purification using a cartridge; a substitute for HPLC.			Purity 98% to ~100% depending on oligo sequence and structure. Yield will gradually decrease as length of oligo increases. Palindromes, hairpins and high GC content oligos and oligos containing stretches of 3 or more G's induces strong secondary structure and base stacking thus decreasing purity and yield. NR* Not Recommended		

*Yield of 30 μg/A₂₆₀ unit for oligos is calculated for an ~equimolar base composition. Long stretches of a single base or homopolymers will have variable yields. Example for homopolymeric 50 mer: A(50) = ~20/A₂₆₀ Unit; G(50) = ~28/A₂₆₀ Unit; T(50) = ~35/A₂₆₀ Unit and C(50) = ~39/A₂₆₀ Unit.

Unmodified DNA Oligo Synthesis*

Scale of Synthesis	Catalog No.	Price (\$)
50 nmol	26-6400-05	0.90
200 nmol	26-6400-02	2.00
1 μmol	26-6400-01	3.75
2 μmol	26-6400-03	6.50
10 μmol	26-6400-10	32.00
15 μmol	26-6400-15	38.00

*minimum charge for 15 mer applies. Please visit www.genelink.com for current list prices. Call for institutional discount pricing structure.

Same Day Oligo*

Design your oligos today and use them tomorrow morning! Investigators who just can not wait order our rush service (order by 12 noon EST). We ship the same day for next early morning delivery in the US and 72 hours for most international destinations.

* Turn-around time stated is for unmodified oligos. Please inquire about purified and modified oligos

Purification



Product	Catalog No.	Scale of Synthesis Price (\$)/ purification					
		50 nmol	200 nmol	1 μmol	2 μmol	10 μmol	15 μmol
Gel Purification	26-6400-XX	75.00	75.00	150.00	280.00	1500.00	1800.00
Reverse Phase Cartridge	26-6400-XX	30.00	30.00	90.00	170.00	750.00	900.00



Quality • Consistency • Confidence

Customer Name: Alyson Rodgers
Customer Number: 10532AJ1
Order Number: 136039
Date: June 13, 2004

Lane	Oligo Name	Sequence (5'-3')	Size	MW	TM	nmols	µg	A ₂₆₀ Units
1.	Primer 1	CATCCTGCAGGGCTAGCTCATAGACTTGCGCGTCAATT AGGATACCTAGG	51	15,715	74.8	48.7	765.3	26.61
2.	Primer 2	GGTGCTTAGATCAGGAGCTTGCGCAGTCCCCGTTGGG GATACC TAGTCACGTACTACTATGTCA	64	19,719	77.9	47.8	941.9	32.11
3.	Primer 3	CATCCTGCAGGGCTAGCTCATAGACTTGCGCGTCAATT AGAGCTTGG	48	14,800	74.6	45.5	672.7	23.10
4.	Primer 4	CTCAAGCAGGAAATCGGGAGCGGCACCTTCGTACGGCG CGTCC	42	12,950	77.0	48.9	633.6	21.92
5.	Primer 5	CGGAATTCGGTCCACAGGCTTGGTCA	25	7,698	63.9	48.2	370.7	12.79
6.	Primer 6	GGTCTGTCTGGGATCCCA	18	5,507	56.6	52.6	289.8	9.49
7.	Primer 7	AAGAGAAAGGTAGGAAGCAC	20	6,257	53.4	40.2	251.5	10.37
8.	Primer 8	CCAACCTCCTGTCCACCACTTTCCTTCGTTGGATGTC CATCTGGGGCGTTATGTTGGTTCTCCTGTAGGACTG GAA	78	23,869	77.7	9.3	222.7	7.26
9.	Primer 9	TGGTCAGAAATCTAGCCTTCGTGACGAAATTTAAACATA AAAGAAAGGCTTCTTGATATATTATCAAGAAACCTTCTT TCTATAAATTTACA	96	29,513	70.2	8.5	252.0	9.14
10.	Primer 10	AATTCAGTACTGTGTTTCAGCAGAAAGGAGTCTTACAT GTGATGGGGGTACAACTGAAAAGTCAAAAAGAAAGTT TGTAATACCAATTTCAATAGCAGTATAAAGGTTCTCTTT GGATCCAGTTGTTGCTGCTTACTACTCTTTCTAGTGC TTAGCT	161	49,736	76.2	6.2	309.4	10.87
11.	Primer 11	CTCAAGCAGGAAATCGAGCGGCACCTTCGTACGTAAT GCCAA	42	12,925	71.1	28.6	369.8	13.53



NOTES

Oligos 1-7 are crude unpurified. Oligos 8-11 are gel purified. Gel lanes for oligos 8-11 correspond to crude followed by gel purified.



Gene Link™

Mobility of an oligonucleotide is dependent upon the size and base composition. Oligos of the same size may not share the same mobility patterns based on the following migration rate C>A>T>G. A stretch G's and GC's induces strong secondary structure that travels as higher mobility fragments.

Custom Oligo Specifications

Gene Link custom oligonucleotides are supplied desalted and lyophilized. They are ready to use after appropriate reconstitution. Dry oligonucleotides are stable at room temperature for an extended period of time.

Storage & Reconstitution

The oligonucleotide should preferably be frozen upon receipt. TE buffer (10 mM Tris, 1 mM EDTA, pH 7.5) is recommended for dissolving the oligonucleotides. After reconstitution store the stock solution at -80°C or -20°C.

Purity & Usage

The crude, desalted oligonucleotide supplied is suitable for all amplification and sequencing protocols. Gel purification is advised for all oligos used for cloning applications and for oligos longer than 50mer.

Gel Photo Documentation

An actual gel picture of the synthesized custom oligonucleotide is supplied. A major single band represents high purity of the crude oligonucleotide.

Biophysical Data

Each oligo after desalting is quantified by recording A₂₆₀. Exact nmols and µg is determined by the extinction coefficient and molecular weight of the oligo.

Primer Design

Successful use of oligos as primers for amplification and sequencing starts with functional primer design followed by optimized PCR amplification conditions.

Fortunately, both PCR and sequencing reactions are inherently 'robust' and have been observed to tolerate wide variations in quality of primers when using unique templates. The same 'tolerance' can also lead to false priming, poor results and frustrating time loss with templates of higher complexity.

Primer specificity alone does not guarantee an optimum amplification yield. Numerous computer applications are available for primer search and design. Most of these applications do not consider the effect of hairpin structures which tend to be quite stable thermodynamically.

General guidelines for primer design are given below followed by a brief account of stable hairpin structure formation and non-Watson-Crick base pairing induced by a stretch of G's and C's interspersed with A's or C's (1-3).

General Guidelines

- 1. Specificity:** Select an 18 to 24-mer stretch with perfect specificity.
- 2. Base Composition:** Preferably maintain GC content below 60% with no stretches of more than 3G's or 4 runs of the same base.
- 3. Tm:** Select primer Tm within a few degrees of the pair.
- 4. Cross Homologies:** Perform NCBI blast to determine extent of cross homologies.
- 5. Secondary Structure:** Perform computer assisted analysis to view formation of stable dimers, loops and hairpins.

Oligo Scale of Synthesis and Typical Yield of Unmodified Oligos*

Scale	Crude Desalted		RPC Purified***		Gel Purified	
	A ₂₆₀ Units	nmols	A ₂₆₀ Units	nmols	A ₂₆₀ Units	nmols
50 nmol	8-10	30+	4-5	12+	NR* [1-2]	NR* [2-4]
200 nmol	20-25	80+	8-12	24+	4-6	8+
1 µmol	100-120	400+	40-50	30+	20-25	40+
Purity & Yield	Purity is more than 80% depending on oligo sequence and structure.		Purity 85% to 95% depending on oligo sequence and structure. Not recommended for oligos longer than 35mer.		Purity 98% to ~100% depending on oligo sequence and structure. Yield will gradually decrease as length of oligo increases.	

*The yield of modified oligos varies based on modification.

**Yield of 30 µg/A₂₆₀ unit for oligos is calculated for an -equimolar base composition. Long stretches of a single base or homopolymers will have variable yields.

Example for homopolymeric 50mer: A(50) = -20/A₂₆₀ Unit; G(50) = -28/A₂₆₀ Unit; T(50) = -35/A₂₆₀ Unit and C(50) = -39/A₂₆₀ Unit.

***RPC is reverse phase purification using a cartridge; a substitute for HPLC.

NR* Not Recommended.

Hairpin Loop Formation and Primer Design*

Sequence	5'-CAGCCGCTACAGGG-ATGAGCT-3' 22 mer; dg = -47.5; Tm(NN) = 61.6°C	5'-GTCCGCACTACAGGATGAGCT-3' 18 mer; dg = -38.4; Tm(NN) = 57.0°C	5'-GTCAOCCGCACTACGGGACAT-3' 21 mer; dg = -46.3; Tm(NN) = 61.70°C	5'-AGTAAAGCCACTACGGACTTACGAC-3' 24 mer; dg = -47.1; Tm(NN) = 58.8°C
*Dimers	5'-TGCAGTACGGAGATCAGCGGC-5' STACK AT 3 IS 4 BP LONG. dg = -4.8; Tm = -58.4°C	5'-GTCCGCACTACAGGATGAGCT-3' + 3'-TACAGGATGCGCGCTG-5' STACK AT 8 IS 6 BP LONG. dg = -5.7; Tm = -42.4°C	5'-GTCCGCACTACAGGATGAGCT-3' + 3'-TACAGGATGCGCGCTG-5' STACK AT 11 IS 6 BP LONG. dg = -4.65; Tm = -28.2°C	5'-AGTAAAGCCACTACGGACTTACGAC-3' + 3'-CAGGATTCAGGCATCAGCAATGA-5' STACK AT 2 IS 4 BP LONG. dg = -2.05; Tm = -47.3°C
Hairpin Loops	None	5'-GTCCGCGAC 3'-TACAGGCATG STEM AT 1 IS 5 BP LONG. LOOP=6. dg = -5.3; Tm = 87.3°C	5'-GTCCGCGCAC 3'-TACAGGCATG STEM AT 6 IS 3 BP LONG. LOOP=6. dg = -2.4; Tm = 68.9°C	5'-AGTAAAGCCACT 3'-CAGGATTCAGGCAT STEM AT 2 IS 4 BP LONG. LOOP=12. dg = -0.8; Tm = 13.8°C

*Secondary structure results are truncated to show the most stable structures. All thermodynamic values including Tm and secondary structures calculated and displayed solely indicate the relative stability of the secondary structures. They should only be used to compare the relative stability of the structures. dg value unit is kcal/mol.

Visit www.genelink.com/tools/gt-500.asp to design oligos or click on the 'Analyze' button while on the online oligo ordering page.


Hairpin Structures

One essential element of efficient primer design is to minimize internal secondary structure, especially hairpin loops which tend to be deceptively stable at standard annealing temperatures. Hairpins are stable with as few as 4 bases stacked in the stem and a loop size of 4 to 6 bases. The stability declines as the loop size increases. The stem and loop size are related proportionately such that longer stem sizes can tolerate longer loop sizes (4).

As a general rule, avoid hairpins with more than 3 bases in the stem. Stable hairpin loop formation drastically reduces the primer concentration available for hybridization to the target sequence.

Base Composition

Higher GC content stabilizes hybridization, but a string of G's and C's can exhibit internal Hoogsteen base pairing, non-Watson-Crick base pairing and should be avoided (3,4). Although this anom-



alous behavior is difficult to predict, these structures can disrupt stable primer binding. In general, avoid runs of more than three consecutive G's in primers. Also, examine potential primers for self-complementary and hairpin structures. Nuclear Magnetic Resonance (NMR) studies have shown that a stable hairpin can form with just four G-C basepairs in the stem and just three bases in the loop (5).

References

1. Michael Zuker (2003) *Nucleic Acids Res.*, 31, 3406-3415.
2. Santalucia, J. (1998) *Proc. Nat. Acad. Sci. USA* 95, 1460.
3. Sarochi, M-T., Courtois, Y., Guschbauer, W. 1970. *Eur. J. Biochem.* 14: 411.
4. Gene Link, Inc. internal data.
5. Summer, M.F., Byrd R.A., Gallo, K.A., Samson, C.J., Zon, G., Egan, W. 1985. *Nucleic Acids Res.* 13: 6375.