

Product Guide

Custom Oligo Synthesis, antisense oligos, RNA oligos, chimeric oligos,
Fluorescent dye labeled oligos, Molecular Beacons, TaqMan Probes
siRNA, Aptamers

Control Fluorescent Probes

For research use only. Not for use in diagnostic procedures for clinical purposes.

Applications

Real Time Quantitative PCR Analysis (QPCR) Probes
Fluorescent Genotyping
siRNA Gene Knockout Validation
Allelic Discrimination
Antisense Targeting
SNP Detection
Aptamers Detection Probes

Endogenous Control Fluorescent Probes

IMPORTANT NOTE

1. Consult product specification sheet and material supplied for specifications of product received.
2. This product guide is not specific to any particular endogenous control fluorescent probe.
3. This product guide should be used in conjunction with the particular instrument manual and specifications and is not intended to replace those specifications.

Storage Instructions:

1. Shipped lyophilized at room temperature.
 2. Store at -20°C upon receipt.
 3. Store at -20°C after reconstitution.
-

Endogenous Control Fluorescent Probes

Reconstitution, Use & Stability of Fluorescent Probes

All Gene Link custom oligo products including, molecular probes, RNA and siRNA includes a datasheet that contains the exact nmols, μg , A_{260} units(OD Units) and other physical data. This data is important for reconstituting the product. All fluorescent probes are shipped in amber tubes to prevent exposure to light and minimize photobleaching. Gene Link guarantees the stability of oligos for 1 year and fluorescent molecular probes for 6 months if reconstituted and stored appropriately as detailed below.

In our experience unmodified oligos are stable for numerous years if reconstituted and stored properly. Avoid multiple freeze thaws; do not exceed 6-10 freeze thaw cycles. If the same oligo is intended to be used repeatedly then it is prudent to make a numerous aliquot of the stock solution and stored frozen.

Reconstitution & Storage

Gene Link probes are supplied lyophilized in amber tubes to protect from light and to reduce photo bleaching. These are stable at room temperature for an extended period but should preferably be frozen upon receipt. TE buffer is recommended for dissolving the probes and oligonucleotides; EDTA inhibits the activity of the nucleases.

Preferred TE Buffer Reconstitution & Storage pH for Fluorescent Probes	
6-FAM, HEX, TET, ROX, and TAMRA	TE Buffer pH 7.5 or 8.0
Cy3, Cy3.5, Cy5, and Cy5.5	TE Buffer pH 7.0 or 7.5
AFDyes	TE Buffer pH 7.5 or 8.0
Cy dyes rapidly degrade in acidic pH	

Further dilution can be made in TE buffer. After reconstitution store the stock solution at -80°C or -20°C . Fluorescently labeled oligos should be stored protected from light.

Preparation of Stock Solution of 100 pmols/ μl [100 μM]

Gene Link provides the exact amount of nmols of each probe supplied on the tube and on the Product Specifications Sheet included with the product.

Multiply the 'nmol' amount by 10 to arrive at the volume of TE to be added.

Example: 2 nmols \times 10 = 20 μL

Dissolve the oligo in 20 μL to get 100 pmols/ μl stock solution [100 μM].

Working Stock Solution.

Dilute 10 fold to prepare a 10pmols/ μL [10 μM]. Dilute further as required.

Storage

For optimal long-term storage, it is recommended that the oligonucleotides should be stored dry at -20°C in the dark. If numerous experiments are planned using the same oligonucleotide, prepare aliquots, dry them and store the aliquots at -20°C .

Stability

Gene Link guarantees the stability of oligos for 1 year and fluorescent molecular probes for 6 months if reconstituted and stored appropriately as recommended by Gene Link. The stability can be increased several fold by instituting proper handling conditions, avoiding exposure to light and multiple freeze thaws.

Fluorescent Probes PAGE Purified Yield*			
Dye/Quencher	50 nmol scale Yield	200 nmol scale Yield	1 μ mol scale Yield
	PAGE Purified*	PAGE Purified*	PAGE Purified*
5' 6-Fam/3'-Tamra, BHQ1 or Dabcyl	10 nmol	25 nmol	60 nmol
5'-Hex, Tet, Cy3 or Cy5/BHQ2	5 nmol	15 nmol	40 nmol
5' -Cy3.5, Cy5.5,/quencher	5 nmol	12 nmol	25 nmol
AFDyes and other NHS Dyes/quencher	2 nmol	5 nmol	16 nmol

*Approximate yield for dye and quencher listed. Lower yield with additional modifications. ** Purified yield is for gel purified (PAGE, polyacrylamide gel electrophoresis). Inquire for higher scale of synthesis

Typical Real Time PCR Component Mix				
Component	Stock Solution	Final Concentration	per 25 μ L rxn	5x 25 μ L reaction
Water	Water		10 μ L	50 μ L
10X PCR buffer	10X	1X	2.5 μ L	12.5 μ L
dNTP	2 mM	200 μ M	2.5 μ L	12.5 μ L
MgCl ₂	25 mM	3 mM	3 μ L	15 μ L
Primer Mix	3 μ M	300 nM (0.3 pmol/ul)	3 μ L	15 μ L
Probe	2 μ M	200 nM (0.2 pmol/ul)	3 μ L	15 μ L
Template			1 μ L	
Taq polymerase	5u/ μ L	0.025 unit/ μ L	0.5 μ L	
Total Volume			25 μ L	


Gene Link Molecular Beacon Melt Curve Protocol

1. Prepare Molecular Beacon stock solution at 100 pmols/ μ l [100 μ M (micromolar)] in 1 X PCR Buffer. Gene Link provides the exact amount of nmoles of each oligo supplied on the tube and on the Oligo Report. Multiply the 'nmol' amount by 10 to arrive at the volume of solvent to be added.
2. Prepare Molecular Beacon working solution at 5 pmols/ μ l [5 μ M (micromolar)] in 1 X PCR Buffer.
3. Set up two 25 μ l reactions, one with probe alone, one with target + probe as follows.


Molecular Beacon Probe Alone	Molecular Beacon Probe plus Target
2.5 μ l 10X PCR buffer	2.5 μ l 10x PCR buffer
3 μ l 25 mM MgCl ₂	3 μ l 25 mM MgCl ₂
1 μ l 5 pmol/ μ l probe [0.2 pmol/ μ l final or 200 nM]	1 μ l 5 pmol/ μ l probe [0.2 pmol/ μ l final or 200 nM]
	3 μ l 5 pmol/ μ l target [0.6 pmol/ μ l final or 600 nM]
18.5 μ l H ₂ O	15.5 μ l H ₂ O

The protocol for Molecular Beacon Melt Curve ramps from 40 to 95 degrees C at 0.2 degrees/second

Screen capture of graphs from a Cepheid SmartCycler



Multiple Molecular Beacon Probe Alone Melt Curve



Molecular Beacon Probe plus Target and Probe Alone Melt Curve

Notes:

- MgCl₂ needs to be higher for MB reactions than for regular PCR as it helps to stabilize the stem structure of the probe during the high ramp rate. Final concentration of MgCl₂ should be between 2.5 and 4 mM. Here we use 3 mM final. This is the same range of concentration used in an actual amplification reaction.
- Final concentration of probe should be 200-600 nM. Here we use 200 nM. This is also the same concentration range used for the real time reaction.
- For a melt curve it is important to saturate the probe with target. Use 2-3 X the *Molar* amount of target. Here we use 3X target for a final concentration of 600 nM. For real time monitoring, 500 ng genomic DNA, diluted 10X to various concentrations can be used as a starting point.

QPCR

Once you have your melt curve you want to select an annealing temperature for your real time PCR where the probe alone is completely closed (shows no fluorescence), and the probe+target is completely open (shows maximal fluorescence). This temperature should be about 5-8 degrees below the T_m of the probe/target hybrid (red vertical line on melt curve of probe+target). It is important to test your primers at the annealing T to ensure that you will have strong, clean amplification at this temperature.

Signal-to-noise (S:N) ratios is calculated by dividing the fluorescence signal of a 25-mer in the presence of a two to five-fold excess of an exactly complementary target sequence by the fluorescence intensity of the probe alone.


Fluorophore Spectral Data & Quencher Selection Guide

Fluorophore Name	Excitation Max, nm	Emission Max, nm	Extinction coefficient*	Color**	Quencher	
AFDye-350 NHS	346	445	19,000	Blue	Dabcyl λ (max) = 453 nm Range = 380-530 nm	
AFDye-405 NHS	402	424	33,000			
PBlue-455 NHS	410	455	46,000			
MBlue-460 NHS	362	459	20,000	Blue-Green	BHQ-1 λ (max) = 534 nm Range = 480-580 nm	
AFDye-488 NHS	494	517	73,000	Yellow-Green		
FAM	495	520	75,850			
TET	521	536	99,000	Yellow		
AFDye-430 NHS	430	539	15,000			
Cal Fluor Gold 540	552	543	81,100			
JOE	520	548	75,000	Yellow-Orange		
Yakima Yellow	531	549	83,800			
AFDye-532 NHS	530	555	81,000			
HEX	535	556	98,000	Orange		BHQ-2 λ (max) = 579 nm Range = 550-650 nm
Cal Orange 560	537	558	81,000			
Cy3	550	570	150,000			
AFDye-555 NHS	555	572	155,000		Yellow-Orange	
TAMRA	555	576	65,000			
CAL Fluor Red 590	569	591	79,000		Orange	
Redmond Red	579	595	52,300			
Cy3.5	581	596	150,000		Orange-Red	
ROX NHS	575	602	82,000			
AFDye-568 NHS	578	602	88,000		Red	
Cal Red 610	590	610	108,000			
TXRed-616 NHS	589	616	69,000			
AFDye-594 NHS	590	617	92,000			
CAL Fluor Red 635	616	637	112,000	Near-IR region. Human vision is insensitive to light beyond ~650 nm; it is not possible to view near-IR fluorescent dyes.		
LC Red 640 NHS	625	640	110,000			
AFDye-647 NHS	649	671	270,000	Red		
Cy5	649	670	250,000			
Cy5.5	675	694	190,000			
AFDye-680 NHS	678	701	185,000	Near-IR region. Human vision is insensitive to light beyond ~650 nm; it is not possible to view near-IR fluorescent dyes.		BBO-650 λ (max) = 650nm Range = 550-750 nm
Cy7 NHS	750	773	199,000			
IR 750 NHS	756	776	260,000			
Cy7.5 NHS	788	808	223,000			

* Extinction coefficient at λ (max) in $\text{cm}^{-1}\text{M}^{-1}$. ** Typical emission color seen through the eyepiece of a conventional fluorescence microscope with appropriate filters.

BioTools application from Gene Link for iPhone/iPod/iPad

BioTools: An Array of Genetic Tools



The BioTools app also has advanced modules for setup of Polymerase Chain Reaction (PCR) and Quantitative Real Time PCR (Q-PCR).

The main focus of this app is to have a handy source of calculation modules and quick reference sections for designing and executing experiments involving PCR and Q-PCR.

BioTools	Reference
<ol style="list-style-type: none"> 1. Oligo Tm: A robust oligo melting temperature calculation module using three methods; it also calculates other physical attributes. 2. Oligo Resuspension 5. Oligo Dilution 6. PCR & QPCR: Convenient calculator for multiple reaction setup for PCR, TaqMan QPCR and Molecular Beacon QPCR setup. Includes stock solution information and cycling profiles 7. Molarity Calculator 8. Reagent Dilution 9. Measurement Converter: A convenient selection of calculators to convert length, area, mass, temperature and volume units. 10. Genetic Code Translator: Enter DNA sequence to see coding pattern. 	<p>A selection of topics, relevant to life scientists for quick access to basic information. This section includes the following sections and sub sections.</p> <ol style="list-style-type: none"> 1. Nucleic Acid Basics 2. Modifications Table: A list of common modifications with molecular structure and basic properties. 3. Dye & Quencher Table: A convenient list of fluorophores and quencher matching the emission max. 4. Duplex Stability & Nuclease Resistance 5. Gene Information: Simply enter the accession number and retrieve detailed gene information from NCBI database, 6. Amino Acid Table: Molecular structure and detailed physical properties of all amino acids. 7. Translation Table: Genetic code for all amino acids. 8. Codon Picker: Select codon sequence and see the corresponding amino acid and detailed information.

Document Warranty and Liability

Information in this document is subject to change without notice. This document and all information presented in this document are written as a guide. Gene Link, Inc. does not warrant this document to be free of errors and assumes no responsibility for any errors that may appear in this document.

Gene Link disclaims all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall Gene Link be liable, whether in contract, tort, warranty, or under any statute or on any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including but not limited to the use thereof.

Website

As the receipt of information on the Internet is highly dependent upon factors, including without limitations to, the user's computer, browser, operation system, etc., information may be perceived incorrectly. Therefore, Gene Link does not warrant or guarantee that the information contained on its website www.genelink.com is error free.

Product Warranty and Liability

Warranty: Gene Link makes no warranty of any kind, specifically disclaims and excludes all other warranties of any kind or nature, directly or indirectly, express or implied, including, without limitation, as to the suitability, productivity, durability, fitness for a particular purpose or use, merchantability, condition, or any other matter with respect to Gene Link products. Gene Link products are for research purposes only including custom products. There is no warranty or claim of its performance for any specific research application. All Gene Link products are guaranteed to meet or exceed the specifications stated. Each Gene Link product is shipped with documentation stating specifications and other technical information. If the product fails to meet the stated specifications the sole remedy is prompt replacement by Gene Link or within 30 days of purchase a refund of the purchased price.

Liability: Under no circumstances shall Gene Link be liable for any damages directly or indirectly related to Gene Link's products and services. Whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of Gene Link products to perform in accordance with the stated specifications.

Research Use Only. Not for use in diagnostic or clinical procedures.

Notice to Purchaser: The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact support@genelink.com.

© 2021 Gene Link Inc. All rights reserved.

The trademarks mentioned herein are the property of their respective owners.

Gene Link, Inc.

Email: support@genelink.com

www.genelink.com