



## Product Manual

### Hepatitis C Virus Genemer™ Amplification Kit

Amplification premix for HCV Viral Specific Fragment

Catalog No.: 60-2003-11K

Store at  $-20^{\circ}\text{C}$

For research use only. Not for use in diagnostic procedures for clinical purposes



- **All human and animal samples, bodily fluids including purified DNA & RNA for pathogen detection should be considered infectious.**
- **Appropriate protective gear (gloves, mask, goggles and disposable laboratory coat) and decontamination protocol should be followed.**
- **All waste should be decontaminated prior to disposal.**
- **Follow approved laboratory, state and local regulations for decontamination protocol and disposal of all samples and laboratory ware.**
- **See appendix for an example of decontamination protocol.**

#### Important Information

All Gene Link, Inc. products are for research use only. Not for use in diagnostic procedures for clinical purposes. Product to be used by experienced researchers who are appropriately trained in performing molecular biology techniques following established safety procedures. Additional qualification and certification is required for interpretation of results.

## Material supplied

# Hepatitis C Virus (HCV) Genemer™ Amplification Kit

Ready-to-use kit for reliable genotyping of a gene fragment. Includes a specific primer pair for gene or mutation specific amplification, optimized buffers and dNTPs. Ready to use; simply add sample and Taq Polymerase. The kit includes nuclease free water and 5X Orange G/XC loading buffer.

## Storage Instructions

1. Shipped on ice.
2. Store at  $-20^{\circ}\text{C}$  upon receipt.

	Qty	Catalog No.	Product	Size
<input type="checkbox"/>		<b>60-2003-11KL</b>	<b>Hepatitis C Virus (HCV) Genemer™ Amplification Kit; 100 Reactions</b>	<b>100 Reactions</b>
<input type="checkbox"/>	1	60-2003-11L	Hepatitis C Virus (HCV) Genemer™ PCR Premix 1 mL	100 Reactions
<input type="checkbox"/>	1	60-2003-06C	Hepatitis C Virus (HCV) Genemer™ Control DNA 20 $\mu\text{L}$ (1ng/ $\mu\text{L}$ )	20 $\mu\text{L}$
<input type="checkbox"/>	1	40-3001-01	Nuclease Free Water	1 ml
<input type="checkbox"/>	1	40-3004-05	Loading buffer 5X Orange G/XC non-denaturing	0.5 ml

	Qty	Catalog No.	Product	Size
<input type="checkbox"/>		<b>60-2003-11KS</b>	<b>Hepatitis C Virus (HCV) Genemer™ Amplification Kit; 50 Reactions</b>	<b>50 Reactions</b>
<input type="checkbox"/>	1	60-2003-11S	Hepatitis C Virus (HCV) Genemer™ PCR Premix 0.5mL	50 Reactions
<input type="checkbox"/>	1	60-2003-06C	Hepatitis C Virus (HCV) Genemer™ Control DNA 20 $\mu\text{L}$ (1ng/ $\mu\text{L}$ )	20 $\mu\text{L}$
<input type="checkbox"/>	1	40-3001-01	Nuclease Free Water	1 ml
<input type="checkbox"/>	1	40-3004-05	Loading buffer 5X Orange G/XC non-denaturing	0.5 ml

\*The polymerase chain reaction (PCR) process is covered by patents owned by Hoffmann-La Roche. A license to perform is automatically granted by the use of authorized reagents.

## Introduction

---

Hepatitis C virus (HCV) is mainly transmitted via blood or blood products. However, sexual, oral and perinatal infections are also possible. Following a general malaise including appetite loss, vomiting and abdominal problems, about 10-20% of patients develop fever, exanthema (skin rash) as well as rheumatoid joint and muscle problems. Two to fourteen days later jaundice develops which may be accompanied by itching. Fulminant hepatitis occurs in about 1% of all infected patients and is frequently fatal. Some of hepatitis C patients develop chronic liver inflammation, which can progress to cirrhosis of liver or primary liver cell carcinoma.

Detection of HCV, involves serologic screening utilizing an enzyme immunoassay (EIA). Serologic assays for antibodies to HCV, i.e., anti-HCV antibodies, are 97% specific. However, these assays cannot be used to distinguish an acute infection from a chronic infection. In three successive versions of EIA, sensitivity has increased progressively. With older tests, some HCV infections may have been missed six to nine months after infection. The most recent third-generation EIA involves core protein and nonstructural proteins 3, 4, and 5; these can be used to detect antibodies within four to ten weeks after the onset of infection. False-negative results for the presence of HCV antibody can occur in persons with compromised immune systems, such as those with HIV-1 infection, patients with renal failure, and those with HCV-associated essential mixed cryoglobulinemia. False-positive EIA results can occur in persons without risk factors and in those without signs of liver disease, such as blood donors or health care workers. Recombinant immunoblot assay (RIBA-2) is used to confirm HCV infection. A positive immunoblot assay result is defined as the detection of antibodies against two or more antigens, and an indeterminate assay result is defined as the detection of antibodies against a single antigen. RIBA-2 is useful to confirm positive EIA results in low-risk populations.

## Detection Methods

---

Gene Link gene detection system uses the Polymerase Chain Reaction (PCR) to amplify the amount of a specific gene fragment DNA or cDNA derived from DNA present in a sample. Amplification of RNA virus templates is based on reverse transcriptase PCR (RT-PCR). The method is highly sensitive and is capable of detecting a few copies of template. Detection systems offered by Gene Link are for research use only and should not be used for clinical diagnosis.

### Qualitative Detection System

This system detects the presence or absence of the specific target DNA or DNA template. It is highly sensitive and gives a positive or negative result. In someone who has never been exposed to the specific pathogen a negative result would occur. On occasion, a positive result is obtained in individuals who have overcome an infection and are clinically asymptomatic. This possibly indicates that minute quantities of pathogen are sometimes present.

### Quantitative Detective System

Quantitative PCR (QPCR) also termed as Real Time PCR determines the quantity of pathogen template in the sample. This system requires the use of fluorescent probe and a real time fluorescent detection system. A specialized PCR system is employed that measures the quantity of amplification at each cycle of the process. The relative amplification per cycle is directly related to the initial amount of DNA or DNA template. This system provides an indication of the number of copies of the pathogen present. This result is used to determine what is known as "viral load" or the number of copies present.

## Procedure



- All human and animal samples, bodily fluids including purified DNA & RNA for pathogen detection should be considered infectious.
- Appropriate protective gear (gloves, mask, goggles and disposable laboratory coat) and decontamination protocol should be followed.
- All waste should be decontaminated prior to disposal.
- Follow approved laboratory, state and local regulations for decontamination protocol and disposal of all samples and laboratory ware.
- See appendix for an example of decontamination protocol.

## Genemer™ Amplification Premix

The Genemer™ Amplification premix is supplied as a 2X concentrated ready to use reagent. Keep on ice and immediately freeze after use. Simply add template and Taq polymerase.

## Amplification and Detection

This protocol assumes that amplification is being performed from cDNA/DNA.

Set up the following amplification files on a thermal cycler. Please refer to the instrument manufacturer's manual for setting up of the program.

Amplification Profile		
Step	Temperature & Time	Cycles
Initial Denaturation	5 min at 95°C	1
Denaturation	15 seconds at 94°C	40
Annealing	30 seconds at 55°C	
Extension	30 seconds at 72°C	
Fill up	7 minutes at 72°C	1
Hold	Hold for infinity at 4°C	Hold

### 1. PCR\* Setup

It is a good strategy initially to amplify multiple amounts for optimization. Use known positive control DNA as one sample. Genemer™ positive control template is provided in this kit. Use 2 µl of this control DNA as one sample. Use one amplification reaction as blank with no added template DNA to determine background amplification due to contamination. The blank reaction should have no amplification product.

PCR Setup		
Component	1 X 20 µl Rxn.	10 X 20 µl Rxns.
Sterile Water	Up to 8 µl	Up to 80 µl
2X Genemer™ PCR Premix	10 µl	100 µl
Taq DNA Polymerase Enzyme or Enzyme mix (EM)	0.5 µl	3 µl
Template DNA	1-5 µl	Add 1-5 µl DNA to each tube
Total Volume	20 µl	
<b>Keep on ice during set up. After adding template start PCR File</b>		

### Recipe

TE Buffer pH 7.5 Composition
1 X TE Buffer pH 7.5
10 mM Tris-HCl pH 7.5
1 mM EDTA



Always use filter barrier pipette tips to prevent cross contamination



Stringent laboratory conditions should be followed to avoid cross contamination. Generally a few copies of template are sufficient for a successful amplification.



Program your thermal cycler instrument with an amplification profile prior to starting the amplification protocol. Consult your appropriate instrument manufacturer's manual.



The PCR premix preparation protocol is written considering that more than one amplification reaction will be performed at the same time. If only one reaction is planned then there is no need to prepare the Taq Enzyme Mix (EM).

### Recipe

Gene Link PCR Buffer
1 X PCR Buffer
10 mM Tris-HCl pH 8.3
50 mM KCl
1.5 mM MgCl <sub>2</sub>
0.001% Gelatin

**2. Taq Polymerase mix Preparation (EM). Label tube "EM"**

Taq Enzyme Mix Preparation (EM)	
Component	10 X 20 µl Rxns.
Sterile Water	33 µl
10 X PCR Buffer	4 µl
Taq Polymerase	3 µl
Add 3 µl to each reaction or to the premix as required	

**Recipe**

1 X TAE Buffer
40 mM Tris-Acetate pH 7.8
1 mM EDTA

**Detection of amplification fragment by agarose gel electrophoresis.**

Prepare a 2 % agarose gel containing 1 µg/ml ethidium bromide in TAE buffer.

1. Add 2 µl of 5X Orange G loading buffer to the amplified samples. Any other non-denaturing loading buffer can be substituted.
  2. Load 10 µl of the amplified samples to the gel.
  3. Electrophorese at 8 volts/cm for approximately 1 hour.
- Visualize under UV. Obtain documentation by taking a gel picture.



Ethidium bromide is a carcinogen. Follow Health and Safety Procedures established by your institution.

Follow proper Hazardous Material Disposal procedures established by your institution.

**Recipe**


6X Orange G loading buffer
10 mM Tris-HCl pH 7.6
60 mM EDTA
60% Glycerol
0.15% Orange G
0.03% Xylene Cyanol

## Results and Interpretation

Refer to the table below to determine the expected size of amplified fragment.

Pathogen Genemer™ Expected Fragment Size		
Pathogen	Control DNA Catalog No.	Fragment Size*
WNV	60-2001-06	154 bp
HIV	60-2002-06	125 bp
HCV	60-2003-06	214 bp
MTB	60-2004-06	151 bp
HBV	60-2007-06	138 bp
MTB 16s RNA	60-2012-06	175 bp

\*Fragment size obtained after PCR amplification



## Troubleshooting

---

### PCR/RT-PCR

#### Precautions:

1. Infectious agents should be handled as per cautious laboratory practices protocol
2. Purified nucleic acids should be amplified immediately, otherwise stored at  $-20\text{ }^{\circ}\text{C}$  for DNA and  $-70\text{ }^{\circ}\text{C}$  for RNA
3. Always use filter barrier pipette tips to prevent cross contamination
4. There should be a physical separation between the Pre-PCR and Post-PCR areas and the flow of movement of reagents and personal should always be unidirectional, i.e. always from Pre-PCR to Post-PCR and not otherwise. This will prevent the possibility of the laboratory being filled with amplicons as aerosol in air or equipment used, which can contaminate potential negative samples and give false positive results.
5. The items of the system should be stored as recommended.

#### Troubleshooting:

- No band in the positive control
  - a. Check the PCR conditions
  - b. Check the post-reconstitution storage of the reagents
  - c. Check the post purification storage of the nucleic acids
  - d. Check the sample collection protocol
  - e. Repeat the PCR reaction after checking the above
- Broad smear in place of sharp bands
  - a. Check the PCR conditions
  - b. Check the voltage in the agarose gel running apparatus
  - c. Repeat the PCR reaction after checking the above
- Two bands in the negative control
  - a. Cross contamination of samples/reagents
  - b. Clean the pre-PCR area with bleach followed by ethanol and water
  - c. Repeat the PCR reaction with fresh reagents

## References

---

1. Yeh SH, Tsai CY, Kao JH, Liu CJ, Kuo TJ, Lin MW, Huang WL, Lu SF, Jih J, Chen DS, Chen PJ. Quantification and genotyping of hepatitis B virus in a single reaction by real-time PCR and melting curve analysis. *J Hepatol.* 2004 Oct; 41(4):659-66.
2. Hanazaki K. Antiviral therapy for chronic hepatitis B: a review. *Curr Drug Targets Inflamm Allergy.* 2004 Mar; 3(1): 63-70. Review.
3. Zeuzem S. Overview of commercial HBV assays systems. *Methods Mol Med.* 2004; 95:3-13. Review.
4. Hu KQ. Occult hepatitis B virus infection and its clinical implications. *J Viral Hepat.* 2002 Jul; 9(4):243-57.

## Appendix

### Pathogen DNA Extraction

1. Add 250 µl of DNA lysis buffer and 10 µl of Glass Beads. Vortex for 1 minute and incubate at RT for 5 minutes.
2. Transfer all the contents into an elution tube. Incubate at 65°C for 15 minutes.
3. Incubate at 95°C for 5 minutes.
4. Incubate at room temperature (RT) for 15 minutes.
5. Centrifuge at 5,000 rpm for 1 minute at RT.
6. Discard supernatant and add 500 µl of DNA wash buffer and vortex for 1 minute. Note: It is important to completely re-suspend the glass bead pellet at this step and in all subsequent steps.
7. Centrifuge at 5,000 rpm for 1 minute at RT. Discard supernatant and add 500 µl of DNA wash buffer and vortex for 1 minute.
8. Centrifuge at 5,000 rpm for 1 minute at RT.
9. Discard supernatant and add 35 µl of RNase free water and vortex for 1 minute.
10. Centrifuge at 5,000 rpm for 1 minute at RT to elute purified DNA.
11. Use eluted DNA for the PCR reaction.
12. The eluted DNA can be stored at this stage in -20°C.



Proper pathogen sample containment and decontamination protocols should be followed. DNA extraction should be performed in a biological safety cabinet with unidirectional workflow for all procedures.



The pathogen RNA purification protocol is specific for the Omni-Pure™ Viral RNA Purification System; Gene Link catalog No: 40-3650-XX.

Other reliable purification protocol and or product can be substituted.



## Appendix

### Decontamination of Bodily Fluids and Tissue Samples

All human and animal samples used for purification of DNA & RNA should be considered infectious and proper decontamination protocol should be followed for eventual disposal. The following protocol is an easy and tested decontamination protocol.

#### Bodily Fluids

1. Prepare 1 L of 1x bleach solution in a large narrow mouth bottle. Keep the bottle capped. See recipe.
2. Transfer all liquid waste to this bottle. You can add up to 300 ml waste to this 1 L bleach solution.
3. At the end of the DNA purification protocol and after at least 1 hour decontamination, this bleach solution can be safely discarded in a regular sink/sewer. Precipitates appear after longer storage.
4. Let cold water run for 3-5 minutes to completely rinse, dilute and wash the sink.



- All bodily fluids and tissue samples are to be considered infectious and hazardous.
- Wear gloves and protective clothing to prevent any exposure.
- All waste materials should be properly decontaminated and disposed following institutional guidelines.
- The decontamination protocol given here is for information only and is not a substitute for any other protocol established by your institution or OSHA.



- Household bleach is a readily available and effective disinfectant.
- Common household bleach contains 5% sodium hypochlorite. This is a convenient 10X solution.
- Extended heating at 80°C to 100°C for 20 minutes or longer denatures and inactivates most pathogens.

#### Solid Waste & Laboratory Plastic Disposables

1. All solid wastes should be disposed of in orange biohazard bags for eventual autoclaving and disposal.
2. All sharps should be disposed in sharps container and disposed of after autoclaving.
3. Paper towels, pipet tips and disposable plastic ware should be treated as solid waste.

#### Recipe

1x Bleach Solution	
Dilution of household bleach	
10x Bleach	Water
100 ml	900 ml

## Ordering Information

**Genemer™:** Primer pair for specific amplification of a gene fragment. Special optimized conditions may be required for certain amplifications.

Catalog No.	Product	Size*	Price \$
60-2001-10	West Nile Virus (WNV) Genemer™	2 nmols	60.00
60-2002-10	Human Immunodeficiency Virus-1 (HIV-1) Genemer™	2 nmols	60.00
60-2003-10	Hepatitis C Virus (HCV) Genemer™	2 nmols	60.00
60-2004-10	Mycobacterium tuberculosis (MTB) Genemer™	2 nmols	60.00
60-2007-10	Hepatitis B Virus (HBV) Genemer™	2 nmols	60.00

\*The quantity supplied is sufficient for 200 regular 20 µl PCR reactions

\*Please visit [www.genelink.com](http://www.genelink.com) for other Genemer™ not listed here

**Genemer™ Amplification Kit:** Ready-to-use kit for reliable genotyping of a gene fragment. Includes a specific primer pair for gene or mutation specific amplification, optimized buffers and dNTPs. Ready to use; simply add sample and Taq Polymerase. The kit includes sterile water and 5X Orange G/XC loading buffer.

Catalog No.	Product	Size	Price \$
60-2001-11KL	West Nile Virus (WNV) Genemer™ Amplification Kit; 100 detections (1 ml)	100 Rxns	140.00
60-2001-11KS	West Nile Virus (WNV) Genemer™ Amplification Kit; 50 detections (0.5 ml)	50 Rxns	80.00
60-2002-11KL	Human Immunodeficiency Virus-1 (HIV-1) Genemer™ Amplification Kit; 100 detections (1 ml)	100 Rxns	140.00
60-2002-11KS	Human Immunodeficiency Virus-1 (HIV-1) Genemer™ Amplification Kit; 50 detections (0.5 ml)	50 Rxns	80.00
60-2003-11KL	Hepatitis C Virus (HCV) Genemer™ Amplification Kit; 100 detections (1 ml)	100 Rxns	140.00
60-2003-11KS	Hepatitis C Virus (HCV) Genemer™ Amplification Kit; 50 detections (0.5 ml)	50 Rxns	80.00
60-2004-11KL	Mycobacterium tuberculosis (MTB) Genemer™ Amplification Kit; 100 detections (1 ml)	100 Rxns	140.00
60-2004-11KS	Mycobacterium tuberculosis (MTB) Genemer™ Amplification Kit; 50 detections (0.5 ml)	50 Rxns	80.00
60-2007-11KL	Hepatitis B Virus (HBV) Genemer™ Amplification Kit; 100 detections (1 ml)	100 Rxns	140.00
60-2007-11KS	Hepatitis B Virus (HBV) Genemer™ Amplification Kit; 50 detections (0.5 ml)	50 Rxns	80.00

\*20 µl PCR reaction volume

\*Please visit [www.genelink.com](http://www.genelink.com) for other Genemer™ Control DNA not listed here

**Genemer™ Control DNA:** Cloned fragment of a particular gene for use with gene or mutation specific Genemer™ products. They are ideal genotyping templates for optimizing and performing control amplification with unknown DNA.

Catalog No.	Product	Size	Price \$
60-2001-06	West Nile Virus (WNV) Genemer™ Control DNA	500 ng	175.00
60-2002-06	Human Immunodeficiency Virus-1 (HIV-1) Genemer™ Control DNA	500 ng	175.00
60-2003-06	Hepatitis C Virus (HCV) Genemer™ Control DNA	500 ng	175.00
60-2004-06	Mycobacterium tuberculosis (MTB) Genemer™ Control DNA	500 ng	175.00
60-2007-06	Hepatitis B Virus (HBV) Genemer™ Control DNA	500 ng	175.00

\*Please visit [www.genelink.com](http://www.genelink.com) for other Genemer™ Kits not listed here