



# Premium Oligonucleotide Synthesis

**QUALITY • CONSISTENCY • CONFIDENCE**

Gene Link oligos are for demanding applications and consistent results. We believe that investigators who value time and have no room for an experiment to fail due to oligo quality should consider Gene Link. Our numerous quality control steps for each oligo assure confidence.



**GOLD STANDARD**

## **Actual Gel Photo**

An actual gel photo of each oligo is affixed on the oligo report. An absolute testimony of quality.

Gene Link has raised the standard since inception over a decade ago.

*We have the pictures to prove it!*

## Superior to “Mass-Produced Factory Oligos”

Gene Link is not an oligo factory. Each oligo is synthesized, processed and quality assured to Gene Link’s absolute standards. This includes coupling efficiency monitoring of each base during synthesis and electrophoretic analysis of each oligo on a polyacrylamide gel to visually assess quality.



### Coupling Efficiency

We maintain a coupling efficiency threshold of greater than 99.5% for all oligos by using premium reagents of exacting specifications, membrane synthesis, state-of-the-art instruments and optimized software-driven protocols. This may not be evident when comparing short oligos, as PCR and sequencing reactions are very robust and can tolerate up to 50% failure/truncated sequence oligos. However, you are clearly taking a chance by using long oligos synthesized at anything below 99.5% coupling efficiency.

### Trityl Monitoring

All Gene Link DNA synthesizers are equipped with trityl monitors for monitoring coupling efficiency of each added base. The instruments are programmed to halt when it falls below the threshold.

See example of routine trityl bars.



Actual trityl coupling efficiency of a 210 mer.

## The Gene Link Advantage

- Stringent Quality Control Measures
- Specializing in Long Oligos up to 250 mer
- Trityl Monitoring of All Oligos
- Polyacrylamide Gel Photograph of Each Oligo
- All Modifications Available
- All Oligo Types Available
- Easy Online Ordering System
- Online Design and Analysis Tools
- Knowledgeable Technical Support
- Personalized, Friendly Customer Service



### Long Oligos

Ask our competitors how often they synthesize 200 to 250 mer oligonucleotides. Gene Link specializes in long oligos.

*You are invited to compare.*



### Oligo Scale of Synthesis and Typical Yield

Scale	Crude Desalted			RPC Purified**			Gel Purified		
	20 mer oligo* Typical yield			30 mer oligo* Typical yield			50 mer oligo* Typical yield		
	A <sub>260</sub> Units	nmols	mg	A <sub>260</sub> Units	nmols	mg	A <sub>260</sub> Units	nmols	mg
50 nmol	8-10	30+	0.2-0.3	4-5	12+	0.1-0.16	NR* [1-2]	NR* [2-4]	NR* [0.03-0.06]
200 nmol	20-25	80+	0.6-0.8	8-12	24+	0.26-0.4	4-6	8+	0.13-0.2
1 μmol	100-120	400+	3-4	40-50	30+	1.3-1.6	20-25	40+	0.6-0.8
<b>Purity &amp; Yield</b>	<b>Purity is greater than 80%</b> depending on oligo sequence and structure. Refer to coupling efficiency table for oligo length dependent purity and yield.  No further purification required for PCR and sequencing applications.  Gel purification recommended for oligos above 50 mer and all applications involving cloning and mutagenesis.			<b>Purity 85% to 95%</b> depending on oligo sequence and structure.  Yield and purity will be lower for sequences with high GC content.  Not recommended for oligos longer than 35 mer.  **RPC is reverse phase purification using a cartridge; a substitute for HPLC.			<b>Purity 98% to ~100%</b> depending on oligo sequence and structure.  Yield will gradually decrease as length of oligo increases. Palindromes, hairpins and high GC content oligos and oligos containing stretches of 3 or more G's induces strong secondary structure and base stacking thus decreasing purity and yield.  NR* Not Recommended		

\*Yield of 30 μg/A<sub>260</sub> unit for oligos is calculated for an ~equimolar base composition. Long stretches of a single base or homopolymers will have variable yields. Example for homopolymeric 50 mer: A(50) = ~20/A<sub>260</sub> Unit; G(50) = ~28/A<sub>260</sub> Unit; T(50) = ~35/A<sub>260</sub> Unit and C(50) = ~39/A<sub>260</sub> Unit.

#### Unmodified DNA Oligo Synthesis\*

Scale of Synthesis	Catalog No.	Price (\$)
50 nmol	26-6400-05	0.90
200 nmol	26-6400-02	2.00
1 μmol	26-6400-01	3.75
2 μmol	26-6400-03	6.50
10 μmol	26-6400-10	32.00
15 μmol	26-6400-15	38.00

\*minimum charge for 15 mer applies. Please visit [www.genelink.com](http://www.genelink.com) for current list prices. Call for institutional discount pricing structure.

#### Same Day Oligo\*

**Design your oligos today and use them tomorrow morning!** Investigators who just can not wait order our rush service (order by 12 noon EST). We ship the same day for next early morning delivery in the US and 72 hours for most international destinations.

\* Turn-around time stated is for unmodified oligos. Please inquire about purified and modified oligos

#### Purification

Product	Catalog No.	Scale of Synthesis Price (\$)/ purification					
		50 nmol	200 nmol	1 μmol	2 μmol	10 μmol	15 μmol
Gel Purification	26-6400-XX	75.00	75.00	150.00	280.00	1500.00	1800.00
Reverse Phase Cartridge	26-6400-XX	30.00	30.00	90.00	170.00	750.00	900.00

## Synthesis of Long Oligos

Synthesis of long oligos up to 250 mer requires greater than 99.5% coupling efficiency. This can only be attained by using reagents of exacting specifications, optimized protocols and state-of-the-art instruments. Gene Link has perfected and maintains these standards. *You are invited to compare.*

**PCR and sequencing reactions are very robust and can tolerate up to 50% failure/truncated sequence oligos.**

**However, you are clearly taking a chance by using long oligos synthesized at anything below 99.5% coupling efficiency. See the coupling efficiency table and graph.**


Gene Link specializes in long oligos. Our description of a long oligo is 180 mer to 250 mer. *You are invited to compare.*

### Purification

Crude oligo is the total yield after chemical synthesis; this contains the full-length product as well as all truncated n-1 sequences. For example, at 99% coupling efficiency the crude yield of a 70 mer is ~50% full-length and ~50% truncated sequences. Gel purification is strongly recommended for all oligos above 50 mer.

### Coupling Efficiency and Full Length Oligo Yield

Oligo Size	99.50%	99.00%	98.00%
20	90.916	82.617	68.123
40	82.243	67.573	45.48
60	74.398	55.268	30.363
80	67.301	45.204	20.27
100	60.881	36.973	13.533
120	55.074	30.24	9.034
140	49.821	24.734	6.031
160	45.068	20.23	4.027
180	40.769	16.546	2.688
200	36.88	13.533	1.795
220	33.36	11.07	1.19
240	30.18	9.05	0.8
250	28.7	8.19	0.65



PCR and sequencing reactions are very robust and can tolerate up to 50% failure/truncated sequence oligos. However, you are clearly taking a chance by using long oligos synthesized at anything below 99.5% coupling efficiency.

### Coupling Efficiency

Chemical DNA synthesis comprises of multiple reactions to complete a cycle of the appropriate base coupling. Thus the use of reagents of exacting specifications, state-of-the-art instruments and optimized software driven protocols are necessary to maintain the highest possible

coupling efficiency. This becomes enormously important when synthesizing a long oligo. Coupling efficiency of 99% or 98% seems very good but on closer examination the yield is almost half for a 40 mer! See the coupling efficiency table.

### Long Oligo Scale of Synthesis and Typical Yield

Gel Purified 150 mer oligo typical yield			
Scale	A <sub>260</sub> Units	nmols	mg
1 μmol	4–6	4+	0.13–0.2
2 μmol	8–12	8+	0.26–0.8

**Purity & Yield 98% to ~100%** depending on oligo sequence and structure.

Yield will gradually decrease as length of oligo increases. Palindromes, hairpins and high GC content oligos and oligos containing stretches of 3 or more G's induces strong secondary structure and base stacking thus decreasing purity and yield.

### Oligo Size and Purification Recommendations

Scale	Synthesis Scale	Recommended Purification
1-49 mer	50 nmol	No purification required. Purification dependent upon desired application.
50-99 mer	200 nmol	Gel purification
100-199 mer	1 μmol	Gel purification
200-250 mer	2 μmol	Gel purification



**Customer Name:** Helen Estrada  
**Customer Number:** 1063643  
**Order Number:** 116194  
**Date:** June 11, 2006

Lane	Oligo Name	Sequence (5'-3')	Size	MW	%G	%C	%A	%T	App. Info.
1	PNBA-F1	ATCTCCGACGCTGGCTCCCTCGAAGGCTCTG... TAC TGGTTCAGACCTGGAGAGGCTTGGATCTC	71	14,168	18.4	9.4	26.6	45.6	871
2	PNBA-F2	TTAAAGCAGCAGCAGCAGCAGCAGCAGCAGCAG... TGGTTCAGACCTGGAGAGGCTTGGATCTC	80	16,201	19.0	7.4	20.6	53.0	928
3	PNBA-F3	GGCTCCGCTCGAAGGCTTGGATCTCAGAGGCT... TAAAGCAGCAGCAGCAGCAGCAGCAGCAGCAG	78	15,866	19.1	9.8	29.6	41.5	836
4	PNBA-F4	AGGCTCCGCTCGAAGGCTTGGATCTCAGAGGCT... TAAAGCAGCAGCAGCAGCAGCAGCAGCAGCAG	88	18,180	19.0	9.0	28.0	44.0	1040
5	PNBA-F5	TACAGCTCCGCTCGAAGGCTTGGATCTCAGAGG... TAAAGCAGCAGCAGCAGCAGCAGCAGCAGCAG	100	21,494	19.0	4.1	18.0	68.9	1100
6	PNBA-F6	CGCTCCGCTCGAAGGCTTGGATCTCAGAGGCT... TAAAGCAGCAGCAGCAGCAGCAGCAGCAGCAG	101	21,526	19.1	2.1	12.0	76.8	1100
7	PNBA-F7	AGGCTCCGCTCGAAGGCTTGGATCTCAGAGGCT... TAAAGCAGCAGCAGCAGCAGCAGCAGCAGCAG	170	32,887	19.0	2.8	19.1	78.1	1800
8	PNBA-F8	CGCTCCGCTCGAAGGCTTGGATCTCAGAGGCT... TAAAGCAGCAGCAGCAGCAGCAGCAGCAGCAG	204	43,628	19.0	1.8	11.0	77.2	436

**NOTE:**  
Gel purified oligos: Gel lane represents crude followed by gel purified.

## Purity and Yield

Gel purified oligo purity is generally between 98% to ~99.9% depending on oligo sequence and structure. Yield will gradually decrease as oligo length increases. Palindromes, hairpins, high GC content oligos and oligos containing stretches of 3 or more G's induce strong secondary structure and base stacking. These are not completely denatured and travel as broad bands on a polyacrylamide gel thus decreasing purity and yield.

## Unmodified DNA Oligo Synthesis\*

Scale of Synthesis	Price (\$)/base
200 nmol scale	2.00
1 µmol scale	3.75
2 µmol scale	6.50

\*minimum charge for 15 mer applies. Please visit [www.genelink.com](http://www.genelink.com) for current list prices. Call for institutional discount pricing structure.

## Purification

All Gene Link oligos shorter than 40 mer usually do not require any further purification if the application is for PCR or sequencing. Gene Link recommends purification of oligos longer than 50 mer and all oligos destined to be cloned.

Product	Scale of Synthesis Price (\$)/purification					
	50 nmol	200 nmol	1 µmol	2 µmol	10 µmol	15 µmol
Gel Purification	75.00	75.00	150.00	280.00	1500.00	1800.00



**GOLD STANDARD**