



Product Specification

Catalog No. 26-3000-08
Lambda gt11 forward primer 24 mers

Sequencing Primers

For research use only. Not for use in diagnostic procedures for clinical purposes. Commercial licenses may be obtained directly from Gene Link.

Product	Lambda gt11 forward primer 24 mers
Catalog No.	26-3000-08
Quantity	25 µg
Shipping Condition	Ambient
Storage	-20°C
Sequence	5'-GGTGGCGACGACTCCTGGAGCCCG-3'
MW	7395
Lot No.	78542

Description

This lambda gt11 forward primer is upstream of the unique EcoR cloning site present in the phage imm434 gene region. This primer can be used to sequence inserts, and can be used with the Lambda gt10 reverse primer (catalog # 26-3000-11) to amplify the insert using polymerase chain reaction.

The product is supplied as a lyophilized powder, after reconstitution store at -20 C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.



Gene Link™

140 Old Saw Mill River Road Hawthorne, NY 10532
Tel: 914.769.1192 www.genelink.com Fax: 914.769.1193