



Product Specification

Catalog No. 26-3000-12
SP6 Promotor primer; 18 mers

Sequencing Primers

For research use only. Not for use in diagnostic procedures for clinical purposes. Commercial licenses may be obtained directly from Gene Link.

Product	SP6 Promotor primer; 18 mers
Catalog No.	26-3000-12
Quantity	25 µg
Shipping Condition	Ambient
Storage	-20°C
Sequence	5'-ATTTAGGTGACACTATAG-3'
MW	5537
Lot No.	78542

Description

This SP6 RNA polymerase primer is complementary to the conserved class class III sequence. Sequencing data will be obtained for all templates containing the Class III T7 RNA polymerase promoter sequences.

The product is supplied as a lyophilized powder, after reconstitution store at -20°C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.



Gene Link™

140 Old Saw Mill River Road Hawthorne, NY 10532
Tel: 914.769.1192 www.genelink.com Fax: 914.769.1193