



# Product Specification Summary

## PCR Amplification & Analysis

Catalog Number	40-3010-50
Product Name	Agarose LE Molecular Biology Grade; 500 g
Size	500 g
Description	Agarose LE Molecular Biology Grade; 500 gms
Component/Note	Agarose GLG is certified molecular biology grade DNase and RNase-free agarose powder. It is suitable for routine DNA and RNA agarose gel electrophoresis.
Component/Note	Agarose GLG is certified molecular biology grade DNase and RNase-free agarose powder. It is suitable for routine DNA and RNA agarose gel electrophoresis. Electrophoretic separation of 200 bp to 13 kb fragment is ideally achieved by preparing a 0.8% to 1% gel in a TAE buffer system.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=1071>



Product Manual

<http://www.genelink.com/Literature/ps/PS40-3010-XX Ver 1.2.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL
Loading buffer 5X BPB/XC non-denaturing; 1 mL	40-3002-10	1 mL
Loading Buffer 5X Orange G/XC non-denaturing; 1 mL	40-3004-10	1 mL
Loading Buffer 2X BPB/XC Denaturing for Sequencing; 1 mL	40-5027-10	1 mL
Omni-Clean™ Gel DNA Beads Purification System; 100 purification	40-4110-10	100 Purifications
Omni-Clean™ Gel DNA Beads Purification System; 500 purifications	40-4110-50	500 Purifications
Omni-Clean™ Gel DNA Spin Column Purification System; 100 purifications	40-4120-10	100 Purifications
Omni-Clean™ Gel DNA Spin Purification; 500 purifications	40-4120-50	500 Purifications
Omni-Clean™ DNA Concentration System; 100 purifications	40-4130-10	100 Purifications
Omni-Clean™ DNA Concentration System; 500 purifications	40-4130-50	500 Purifications
Omni-Clean™ DNA Spin Column Concentration System; 100 purification	40-4140-10	100 Purifications
Omni-Clean™ DNA Spin Column Concentration System; 500 purification	40-4140-50	500 Purifications
Taq DNA Polymerase 400 units 5 u/ul; 80 uL; 400 units	40-5200-40	400 each