



# Product Specification Summary

## Huntington Disease

Catalog Number	40-2025-15FMS
Product Name	Huntington Disease GScan™ 6-FAM Kit (20rxns); 1 Kit
Size	1 Kit
Description	Huntington Disease GScan™ 6-FAM Kit (20rxns); 1 Kit
Component/Note	6-Fam Dye
Component/Note	Fluorescent Detection

Gene Link's GScan™ gene detection kits are safe, convenient and sensitive, and afford automated compilation of data. The kits contain optimized PCR amplification reagents and a wide array of fluorescent-labeled primers for genotyping after PCR using fluorescent genetic analyzer instrument(s). Included in these kits are ready-to-run control samples of various repeats of the triple repeat disorder kit. These control samples are for calibration with the molecular weight markers for accurate size determination of the amplified fragments. The GScan kits are simple and robust for routine triple-repeat detection of greater than 100 repeats of all triple repeat disorders listed, except Fragile X. The CGG repeat in Fragile X can be detected up to ~50 repeats.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=11645>



Product Manual

[http://www.genelink.com/Literature/ps/M40-2025-15\\_Ver7.2.pdf](http://www.genelink.com/Literature/ps/M40-2025-15_Ver7.2.pdf)



Product MSDS

[http://www.genelink.com/Literature/ps/SDS\\_40-20XX-15.pdf](http://www.genelink.com/Literature/ps/SDS_40-20XX-15.pdf)



## Related Products

Product	Catalog No	Size
GLHD 7 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-05HX	25 uL
GLHD 18 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-01HX	25 uL
GLHD 31 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-07HX	25 uL
GLHD 34 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-02HX	25 uL
GLHD 37 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-08HX	25 uL
GLHD 44 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-03HX	25 uL
GLHD 89 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-04HX	25 uL
GLHD 49 CAG repeat Scan Hex labeled ready to run Control DNA; 25 uL	40-2025-09HX	25 uL
GLHD 116 CAG repeat GScan Genemer Control DNA; 25 uL	40-2025-06HX	25 uL
GLHD 134 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-61HX	25 uL
GLHD 182 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-62HX	25 uL
Huntington Disease GScan™ TET Kit	40-2025-15TT	1 Kit
Huntington Disease GScan™ HEX Kit	40-2025-15HX	1 Kit
Huntington Disease GScan™ Cy3 Kit	40-2025-15C3	1 Kit
Huntington Disease GScan™ Cy5 Kit	40-2025-15C5	1 Kit