



Product Specification Summary

Huntington Disease

| | |
|----------------|--|
| Catalog Number | 40-2025-06 |
| Product Name | Huntington Disease 116 CAG Genemer Control DNA; 500 ng |
| Size | 500 ng |
| Description | GLHD Genemer |
| Component/Note | GLHD ~116 CAG repeats |
| Component/Note | Huntington Disease |

A tube containing 500 ng of lyophilized control DNA segment. The control DNA's contains varying triple repeats. These are specific DNA templates for Gene Link's Genemer and PCRProber products. The above control DNA is an ideal genotyping template for optimizing and performing control amplification with unknown DNA. The size of the triple repeats has been determined by sequencing and gel electrophoresis. The stability of size repeats upon cloning and amplification has NOT been determined. Thus, the size should be considered approximate and there is no claim for each fragment to contain the exact number of triple repeats. The quantity supplied is sufficient for 1000 regular 50ul PCR reaction.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=754>



Product Manual

http://www.genelink.com/Literature/ps/PS40-2025-05_Ver3.1.pdf



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|---|--------------|-------|
| GLHD 7 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-05HX | 25 uL |
| GLHD 18 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-01HX | 25 uL |
| GLHD 31 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-07HX | 25 uL |
| GLHD 34 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-02HX | 25 uL |
| GLHD 37 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-08HX | 25 uL |
| GLHD 44 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-03HX | 25 uL |
| GLHD 89 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-04HX | 25 uL |
| GLHD 49 CAG repeat Scan Hex labeled ready to run Control DNA; 25 uL | 40-2025-09HX | 25 uL |
| GLHD 116 CAG repeat GScan Genemer Control DNA; 25 uL | 40-2025-06HX | 25 uL |
| GLHD 134 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-61HX | 25 uL |
| GLHD 182 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-62HX | 25 uL |
| Huntington Disease GScan™ TET Kit | 40-2025-15TT | 1 Kit |
| Huntington Disease GScan™ HEX Kit | 40-2025-15HX | 1 Kit |
| Huntington Disease GScan™ Cy3 Kit | 40-2025-15C3 | 1 Kit |
| Huntington Disease GScan™ Cy5 Kit | 40-2025-15C5 | 1 Kit |