



Product Specification Summary

RAPD Decamer Sets

| | |
|----------------|---------------------------|
| Catalog Number | 40-0001-04 |
| Product Name | GL RAPD Decamer 04-GH Set |
| Size | 1 set |
| Description | GL RAPD Decamer 04-GH Set |
| Component/Note | Supplied Lyophilized |
| Component/Note | 40 Decamers 10 nmols each |

Arbitrary sequence decamer sets for rapid discerning of DNA polymorphism. There are 30 sets available each containing 40 primers with a 60-70% GC content. These primers are used in amplification reactions for determining Random Amplified Polymorphic DNA (RAPD) fragments. Each tube contains 10 nmols of lyophilized primer. This is sufficient for thousands of amplification reaction. Please see Protocol section for details. Please see the pdf file for detailed information

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=724>



Product Manual

<http://www.genelink.com/Literature/ps/PS40-0001-04.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|-----------------------------|------------|-------|
| GL RAPD Decamer Set 1- 30 | 40-0001-00 | 1 set |
| GL RAPD Decamer 10-ST Set | 40-0001-10 | 1 set |
| GL RAPD Decamer 11-UV Set | 40-0001-11 | 1 set |
| GL RAPD Decamer 13-YZ Set | 40-0001-13 | 1 set |
| GL RAPD Decamer 16-AEAF Set | 40-0001-16 | 1 set |
| GL RAPD Decamer 17-AGAH Set | 40-0001-17 | 1 set |
| GL RAPD Decamer 18-AIAJ Set | 40-0001-18 | 1 set |
| GL RAPD Decamer 19-AKAL Set | 40-0001-19 | 1 set |
| GL RAPD Decamer 20-AMAN Set | 40-0001-20 | 1 set |
| GL RAPD Decamer 21-AOAP Set | 40-0001-21 | 1 set |
| GL RAPD Decamer 22-AQAR Set | 40-0001-22 | 1 set |
| GL RAPD Decamer 23-ASAT Set | 40-0001-23 | 1 set |
| GL RAPD Decamer 24-AUAV Set | 40-0001-24 | 1 set |
| GL RAPD Decamer 25-AWAX Set | 40-0001-25 | 1 set |
| GL RAPD Decamer 26-AYAZ Set | 40-0001-26 | 1 Set |