



## Product Specification Summary

### RAPD Decamer Sets

|                |                           |
|----------------|---------------------------|
| Catalog Number | 40-0001-05                |
| Product Name   | GL RAPD Decamer 05-IJ Set |
| Size           | 1 set                     |
| Description    | GL RAPD Decamer 05-IJ Set |
| Component/Note | Supplied Lyophilized      |
| Component/Note | 40 Decamers 10 nmols each |
| Storage        | Store at -20°C            |

Arbitrary sequence decamer sets for rapid discerning of DNA polymorphism. There are 30 sets available each containing 40 primers with a 60-70% GC content. These primers are used in amplification reactions for determining Random Amplified Polymorphic DNA (RAPD) fragments.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=725>



Product Manual

<http://www.genelink.com/Literature/ps/PS40-0001-05.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

| Product                     | Catalog No | Size  |
|-----------------------------|------------|-------|
| GL RAPD Decamer Set 1- 30   | 40-0001-00 | 1 set |
| GL RAPD Decamer 10-ST Set   | 40-0001-10 | 1 set |
| GL RAPD Decamer 11-UV Set   | 40-0001-11 | 1 set |
| GL RAPD Decamer 13-YZ Set   | 40-0001-13 | 1 set |
| GL RAPD Decamer 16-AEAF Set | 40-0001-16 | 1 set |
| GL RAPD Decamer 17-AGAH Set | 40-0001-17 | 1 set |
| GL RAPD Decamer 18-AIAJ Set | 40-0001-18 | 1 set |
| GL RAPD Decamer 19-AKAL Set | 40-0001-19 | 1 set |
| GL RAPD Decamer 20-AMAN Set | 40-0001-20 | 1 set |
| GL RAPD Decamer 21-AOAP Set | 40-0001-21 | 1 set |
| GL RAPD Decamer 22-AQAR Set | 40-0001-22 | 1 set |
| GL RAPD Decamer 23-ASAT Set | 40-0001-23 | 1 set |
| GL RAPD Decamer 24-AUAV Set | 40-0001-24 | 1 set |
| GL RAPD Decamer 25-AWAX Set | 40-0001-25 | 1 set |
| GL RAPD Decamer 26-AYAZ Set | 40-0001-26 | 1 Set |