



Product Specification Summary

Huntington Disease

| | |
|----------------|---|
| Catalog Number | 40-2025-15FMS |
| Product Name | Huntington Disease GScan™ 6-FAM Kit (20rxns); 1 Kit |
| Size | 1 Kit |
| Description | Huntington Disease GScan™ 6-FAM Kit (20rxns); 1 Kit |
| Component/Note | 6-Fam Dye |
| Component/Note | Fluorescent Detection |

Gene Link's GScan™ gene detection kits are safe, convenient and sensitive, and afford automated compilation of data. The kits contain optimized PCR amplification reagents and a wide array of fluorescent-labeled primers for genotyping after PCR using fluorescent genetic analyzer instrument(s). Included in these kits are ready-to-run control samples of various repeats of the triple repeat disorder kit. These control samples are for calibration with the molecular weight markers for accurate size determination of the amplified fragments. The GScan kits are simple and robust for routine triple-repeat detection of greater than 100 repeats of all triple repeat disorders listed, except Fragile X. The CGG repeat in Fragile X can be detected up to ~50 repeats.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=11645>



Product Manual

http://www.genelink.com/Literature/ps/M40-2025-15_Ver7.2.pdf



Product MSDS

http://www.genelink.com/Literature/ps/SDS_40-20XX-15.pdf



Related Products

| Product | Catalog No | Size |
|---|--------------|-------|
| GLHD 7 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-05HX | 25 uL |
| GLHD 18 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-01HX | 25 uL |
| GLHD 31 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-07HX | 25 uL |
| GLHD 34 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-02HX | 25 uL |
| GLHD 37 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-08HX | 25 uL |
| GLHD 44 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-03HX | 25 uL |
| GLHD 89 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-04HX | 25 uL |
| GLHD 49 CAG repeat Scan Hex labeled ready to run Control DNA; 25 uL | 40-2025-09HX | 25 uL |
| GLHD 116 CAG repeat GScan Genemer Control DNA; 25 uL | 40-2025-06HX | 25 uL |
| GLHD 134 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-61HX | 25 uL |
| GLHD 182 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL | 40-2025-62HX | 25 uL |
| Huntington Disease GScan™ TET Kit | 40-2025-15TT | 1 Kit |
| Huntington Disease GScan™ HEX Kit | 40-2025-15HX | 1 Kit |
| Huntington Disease GScan™ Cy3 Kit | 40-2025-15C3 | 1 Kit |
| Huntington Disease GScan™ Cy5 Kit | 40-2025-15C5 | 1 Kit |