



Product Specification Summary

Huntington Disease

Catalog Number	40-2025-06
Product Name	Huntington Disease 116 CAG Genemer Control DNA; 500 ng
Size	500 ng
Description	GLHD Genemer
Component/Note	GLHD ~116 CAG repeats
Component/Note	Huntington Disease

A tube containing 500 ng of lyophilized control DNA segment. The control DNA's contains varying triple repeats. These are specific DNA templates for Gene Link's Genemer and PCRProber products. The above control DNA is an ideal genotyping template for optimizing and performing control amplification with unknown DNA. The size of the triple repeats has been determined by sequencing and gel electrophoresis. The stability of size repeats upon cloning and amplification has NOT been determined. Thus, the size should be considered approximate and there is no claim for each fragment to contain the exact number of triple repeats. The quantity supplied is sufficient for 1000 regular 50ul PCR reaction.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=754>



Product Manual

http://www.genelink.com/Literature/ps/PS40-2025-05_Ver3.1.pdf



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
GLHD 7 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-05HX	25 uL
GLHD 18 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-01HX	25 uL
GLHD 31 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-07HX	25 uL
GLHD 34 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-02HX	25 uL
GLHD 37 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-08HX	25 uL
GLHD 44 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-03HX	25 uL
GLHD 89 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-04HX	25 uL
GLHD 49 CAG repeat Scan Hex labeled ready to run Control DNA; 25 uL	40-2025-09HX	25 uL
GLHD 116 CAG repeat GScan Genemer Control DNA; 25 uL	40-2025-06HX	25 uL
GLHD 134 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-61HX	25 uL
GLHD 182 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-62HX	25 uL
Huntington Disease GScan™ TET Kit	40-2025-15TT	1 Kit
Huntington Disease GScan™ HEX Kit	40-2025-15HX	1 Kit
Huntington Disease GScan™ Cy3 Kit	40-2025-15C3	1 Kit
Huntington Disease GScan™ Cy5 Kit	40-2025-15C5	1 Kit