



# Product Specification Summary

## Genetic Tools & Reagents

Catalog Number	40-5024-20
Product Name	Lumisol III,Hybridization Solution for Oligo Probe; 200 mL
Size	200 mL
Description	Lumisol III, Hybridization Solution for Oligo Probes
Component/Note	Product Description Lumisol™ II & III are ready-to-use hybridization solutions specifically formulated for non-radioactive chemiluminescent hybridization.
Component/Note	Lumisol™ III hybridization solution is specifically formulated for oligonucleotide probes that require only 30 minute pre-hybridization and 30 minute hybridization. Longer hybridization is not recommended. The hybridization temperature should be calculated for proper results and should not be more than 55°C for alkaline phosphatase labeled oligo probes.
<b>Product Description</b>	

Lumisol™ II & III are ready-to-use hybridization solutions specifically formulated for non-radioactive chemiluminescent hybridization.

**Lumisol™ II** Lumisol™ II is for use with digoxigenin or biotin labeled DNA or RNA probes in Southern & Northern blot protocols. It can also be used with radioactive labeled probes.

**Lumisol™ III** Lumisol™ III is specifically formulated for use with oligonucleotide probes labeled with biotin, digoxigenin, alkaline phosphatase or other detection ligand or enzymes.

**Applications** Lumisol™ II hybridization solution can be used for all stringent types of nucleic acid blot hybridization conditions, particularly using non-radioactive digoxigenin labeled probes. Lumisol™ II hybridization solution can be used for prehybridization and hybridization. The hybridization temperature should be calculated for proper results. Duration of hybridization can be reduced to 6 hours and overnight hybridization can be used for high sensitivity requirements and for convenience. Lumisol™ III hybridization solution is specifically formulated for oligonucleotide probes that require only 30 minute pre-hybridization and 30 minute hybridization. Longer hybridization is not recommended. The hybridization temperature should be calculated for proper results and should not be more than 55°C for alkaline phosphatase labeled oligo probes.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=95>



Product Manual

[http://www.genelink.com/Literature/ps/PS40-5023\\_Lumisol\\_Ver\\_1.5.pdf](http://www.genelink.com/Literature/ps/PS40-5023_Lumisol_Ver_1.5.pdf)



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni- Marker™ Universal unlabeled; 500 uL	40-3005-05	500 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni- Marker™ Low unlabeled; 100 uL	40-3006-01	100 uL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL
Omni- Marker™ Low unlabeled; 1 mL	40-3006-10	1 mL
GScan Marker Tamra Labeled 50 bp - 600 bp; 100 uL	40-3061-01	100 uL
GScan Marker Tamra Labeled 50 bp - 600 bp; 500 uL	40-3061-05	500 uL
Omni-Marker™ GScan-1 Tet labeled 50 bp-1kb	40-3071-01	100 ul
Omni-Marker™ GScan-1 Tet labeled 50 bp-1kb	40-3071-05	500 ul
Omni-Marker™ GScan-2 Tet labeled 50 bp- 600 bp	40-3072-01	100 ul
Omni-Marker™ GScan-2 Tet labeled 50 bp- 600 bp	40-3072-05	500 ul
GScan Marker Hex Labeled 50 bp - 600 bp; 100 uL	40-3081-01	100 uL