

# Directions to Gene Link, Inc.

140 Old Saw Mill River Road, Hawthorne, NY 10532

## **From Long Island and Queens**

From Whitestone Bridge or Throgs Neck Bridge, continue straight onto the Hutchinson River Parkway to exit 13-Cross County Parkway West. At end of Cross County Parkway, bear right onto Saw Mill Parkway North. Proceed to exit 25-9A North/Hawthorne. Off the exit ramp make a left at traffic light onto 9A North. Proceed about 3/10 mile to Old Saw Mill River Road-a soft right off 9A. Follow about 3/10 mile to driveway on right.

## **From George Washington Bridge**

Take I-87 North-Major Deegan Expressway through Yonkers to Saw Mill Parkway North. Follow directions above.

## **From Yonkers and Upper Westchester**

Saw Mill Parkway to exit 25. Follow directions above.

## **From Connecticut**

Follow Cross Westchester Expressway/287 West to exit 2-9A North. Follow for 3 miles. Take soft right after traffic light at Skyline Drive onto Old Saw Mill River Road. Follow directions above.

## **From West Side Manhattan**

From West Side take West Side Highway which becomes Henry Hudson Parkway which becomes Saw Mill Parkway to exit 25. Follow directions above.

## **From Upstate New York/Tappan Zee Bridge**

Off Bridge take exit 8-Cross Westchester Expressway/287 East to exit 1-Saw Mill Parkway North. Follow directions above.

## **From East Side Manhattan**

From East Side take FDR Drive to Triboro Bridge to Bruckner Expressway to Bronx River Parkway North. Continue to Taconic/Sprain Brook Parkways North. Take route 141/Hawthorne exit. Turn left at end of ramp and left again at stop sign onto route 141. Follow signs to route 9A South. From 9A South, make a soft left after Citgo gas station onto Old Saw Mill River Road. Gene Link is third driveway on left.

## **From Upstate New York/East of Hudson River**

Take Taconic State Parkway South to route 100 exit/Briarcliff Manor. Turn left at traffic light onto 100 South. Follow to 9A South for about 2 miles. Follow directions above.

