



Product Specification

Catalog No. 26-3100-18
Bam HI/Not I Adaptor

Adaptors

For research use only. Not for use in diagnostic procedures for clinical purposes.
Commercial licenses may be obtained directly from Gene Link.

Product	Bam HI /Not I Adaptor
Catalog No.	26-3000-18
Quantity	20 µg
Shipping Condition	Ambient
Storage	-20°C
Sequence	5'- GATCCGCGGCCGC-3' 3'- CGCCGGCGp -5'
MW	6,348
Lot No.	132852

Description

Gene Link supplied adaptors are short synthetic oligonucleotide preannealed duplexes with 5' blunt end. These have an internal restriction endonuclease site, which is created by ligation to fragments with complementary overhangs. The duplexes have an overhang and a blunt end. Adaptors are supplied as lyophilized powder.

The product is supplied as a lyophilized powder, after reconstitution store at -20°C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.

