



## Product Specification

Catalog No. 26-3100-22  
Eco RI/Bam HI Adaptor

# Adaptors

**For research use only.** Not for use in diagnostic procedures for clinical purposes. Commercial licenses may be obtained directly from Gene Link.

<b>Product</b>	<b>Eco RI/Bam HI Adaptor</b>
<b>Catalog No.</b>	<b>26-3100-22</b>
<b>Quantity</b>	<b>20 µg</b>
<b>Shipping Condition</b>	<b>Ambient</b>
<b>Storage</b>	<b>-20°C</b>
<b>Sequence</b>	<b>5'- AATTCGGATCC-3'</b> <b>3'-GCCTAGGp-5'</b>
<b>MW</b>	<b>5597.6</b>
<b>Lot No.</b>	<b>219317</b>

### Description

Gene Link supplied adaptors are short synthetic oligonucleotide pre-annealed duplexes with 5' blunt end. These have an internal restriction endonuclease site, which is created by ligation to fragments with complementary overhangs. The duplexes have an overhang and a blunt end. Adaptors are supplied as lyophilized powder.

The product is supplied as a lyophilized powder, after reconstitution store at -20°C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.