



Product Specification

Catalog No. 26-3100-23
Sal I/Bam HI Adaptor

Adaptors

For research use only. Not for use in diagnostic procedures for clinical purposes. Commercial licenses may be obtained directly from Gene Link.

Product	Sal I/Bam HI Adaptor
Catalog No.	26-3100-23
Quantity	20 µg
Shipping Condition	Ambient
Storage	-20°C
Sequence	5'- TCGACGGATCC-3' 3'-GCCTAGGp-5'
MW	5598.6
Lot No.	219317

Description

Gene Link supplied adaptors are short synthetic oligonucleotide pre-annealed duplexes with 5' blunt end. These have an internal restriction endonuclease site, which is created by ligation to fragments with complementary overhangs. The duplexes have an overhang and a blunt end. Adaptors are supplied as lyophilized powder.

The product is supplied as a lyophilized powder, after reconstitution store at -20°C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.