



Product Specification

Catalog No. 26-3000-02
M13/pUC Reverse (-24); 16mers

Sequencing Primers

For research use only. Not for use in diagnostic procedures for clinical purposes. Commercial licenses may be obtained directly from Gene Link.

Product	M13/pUC Reverse (-24); 16mers
Catalog No.	26-3000-02
Quantity	25 µg
Shipping Condition	Ambient
Storage	-20°C
Sequence	5'AACAGCTATGACCATG-3'
MW	4874
Lot No.	78542

Description

All M13/pUC primers are complementary to sequence in the M13 and pUC vectors (1). These primers will yield sequencing data for inserts cloned in these vectors. In addition to these vectors sequencing data will be obtained for all templates/vectors containing the *lacZ* sequences.

The product is supplied as a lyophilized powder, after reconstitution store at -20 C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.

1. Messing J. (1983) Meth. Enzymol. 101:20.



Gene Link™

140 Old Saw Mill River Road Hawthorne, NY 10532
Tel: 914.769.1192 www.genelink.com Fax: 914.769.1193