



Product Specification

Catalog No. 26-3000-06
T3 Promotor Primer; 20 mers

Sequencing Primers

For research use only. Not for use in diagnostic procedures for clinical purposes. Commercial licenses may be obtained directly from Gene Link.

Product	T3 Promotor Primer; 20 mers
Catalog No.	26-3000-06
Quantity	25 µg
Shipping Condition	Ambient
Storage	-20°C
Sequence	5'-TAATACGACTCACTATAGGGAGA-3'
MW	7080
Lot No.	78542

Description

This T3 RNA polymerase primer is complementary to the conserved class III sequence. Sequencing data will be obtained for all templates containing the Class III T 7 RNA polymerase promoter sequences.

The product is supplied as a lyophilized powder, after reconstitution store at -20°C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.



Gene Link™

140 Old Saw Mill River Road Hawthorne, NY 10532
Tel: 914.769.1192 www.genelink.com Fax: 914.769.1193