



## Product Specification

Catalog No. 26-3000-22  
T3 Promotor/Bluescript primer 20mers

# Sequencing Primers

**For research use only.** Not for use in diagnostic procedures for clinical purposes. Commercial licenses may be obtained directly from Gene Link.

<b>Product</b>	<b>T3 Promotor/Bluescript primer 20mers</b>
<b>Catalog No.</b>	<b>26-3000-22</b>
<b>Quantity</b>	<b>25 µg</b>
<b>Shipping Condition</b>	<b>Ambient</b>
<b>Storage</b>	<b>-20°C</b>
<b>Sequence</b>	<b>5'-AATTAACCCTCACTAAAGGG-3'</b>
<b>MW</b>	<b>6094</b>
<b>Lot No.</b>	<b>78542</b>

### Description

This T3 RNA polymerase primer is complementary to the conserved class III sequence. Sequencing data will be obtained for all templates containing the Class III T3 RNA polymerase promotor sequences. This primer in addition contains sequence complementary to Stratagenes's pBluescript plasmid flanking the T3 sequence.

The product is supplied as a lyophilized powder, after reconstitution store at -20°C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.



**Gene Link**™

140 Old Saw Mill River Road Hawthorne, NY 10532  
Tel: 914.769.1192 www.genelink.com Fax: 914.769.1193