

Product Specification



Primers

Product	Xho I /Not I
Catalog No	26-3100-06
Quantity	20 µg
Ship	Ambient
Storage	-20°C
Sequence	5'-TCGAGGCGGCCGC-3' -----3'-CCGCCGGCGP 5'

MW

Description

Gene Link supplied adaptors are short synthetic oligonucleotide preannealed duplexes with 5' blunt end phosphorylated. These have an internal restriction endonuclease site which is created by ligation to fragments with complementary overhangs. The duplexes have an overhang and a blunt end. Adaptors are supplied as lyophilized powder.

The product is supplied as a lyophilized powder, after reconstitution store at -20°C. Oligo purity is greater than 98% as determined by denaturing polyacrylamide gel electrophoresis.



Gene Link™

140 Old Saw Mill River Road Hawthorne, NY 10532
Tel: 914.769.1192 www.genelink.com Fax: 914.769.1193