



## Product Specification Summary

### Fragile X Syndrome

Catalog Number	40-2004-11
Product Name	Fragile X Genemer Kit; 1 Kit
Size	1 Kit
Component/Note	repeat region
Component/Note	Fragile X CGG triple

Genemer™ are optimized primer pairs for amplification of the gene fragment of interest, particularly those spanning a mutation. A wide range is available from single base pair mutation sites, the complex triple repeat disorders and various pathogens. Genemer™ kits are complete, easy-to-use systems for reliable genotyping of large repeats in certain triple repeat disorder. The Genemer™ product contains 10 nmols of the primer pair that is sufficient for 400 regular 50 µl amplifications. The Genemer™ Kit contains reagents as well as control DNA sufficient for 100 amplifications.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=1008>



Product Manual

<http://www.genelink.com/Literature/ps/M40-2004-11Ver5.4.pdf>



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Fragile X GeneScan TET Kit; 1 Kit	40-2004-15TT	1 Kit
Fragile X GScan V2 HEX Kit; 1 Kit	40-2004-15HX	1 Kit
Fragile X GeneScan Cy3 Kit; 1 Kit	40-2004-15C3	1 Kit
Fragile X GScan™ V2 6-FAM Kit; 1 Kit	40-2004-15FM	1 Kit
Fragile X GScan™ V2 6-FAM Kit 20 rxns; 1 Kit	40-2004-15FMS	1 Kit
Fragile X GScan™ Cy5 Kit (20rxns); 1 Kit	40-2004-15C5S	1 Kit
Fragile X GeneScan Cy5 Kit; 1 Kit	40-2004-15C5	1 Kit
Fragile X Genemer Kit; 1 Kit	40-2004-11	1 Kit
FX TaqMan 6 CGG 5'FAM; 25 ug	44-0130-06	25 ug
FX TaqMan 6 CGG 5' Fluorescein; 25 ug	44-0140-06	25 ug
FX CGG 5-MB-CY3; 25 ug	45-0110-05	25 ug
FX CGG 5-MB-CY5; 25 ug	45-0120-05	25 ug
FX CGG 5-MB-FAM; 25 ug	45-0130-05	25 ug
FX CGG 5-MB-FI; 25 ug	45-0140-05	25 ug
FX CGG 5-MB-HEX; 25 ug	45-0150-05	25 ug