



# Product Specification Summary

## siRNA:RNA Interference

Catalog Number	27-6411-20FL
Product Name	siRNA non-silencing FL labeled; 20 nmols
Size	20 nmols
Description	siRNA non-silencing FI labeled. Negative control
Component/Note	siRNA non-silencing FL labeled duplex

A negative control should also be included when performing siRNA experiments to eliminate the possibility of nonspecific silencing effects. For this purpose, Gene Link provides a siRNA sequence that has been shown to have no effect on gene silencing. The negative control siRNA is also available labeled with fluorescein, which allows transfection efficiency to be tracked by fluorescent microscopy.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1011>



Product Manual [http://www.genelink.com/Literature/ps/RNAi\\_Controls\\_COA\\_PS\\_V6.1.pdf](http://www.genelink.com/Literature/ps/RNAi_Controls_COA_PS_V6.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit	27-6402-01	1 Kit
3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each	27-6402-02	1 Kit

Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit  
8 Westchester Plaza, Suite 130, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | [www.genelink.com](http://www.genelink.com) |  
SS27-6411-20FL.pdf Print Date Version : April 7, 2026

27-6402-06  
Page 1 of 1

