

Product Specification Summary

siRNA:RNA Interference

Catalog Number 27-6413-10FL

Product Name Human Beta Actin control siRNA duplex FL; 10 nmols

Size 10 nmols

Description Human Beta Actin siRNA; NM 001101 Fl labeled

A negative control should also be included when performing siRNA experiments to eliminate the possibility of nonspecific silencing effects. For this purpose, Gene Link provides a siRNA sequence that has been shown to have no effect on gene silencing. The negative control siRNA is also available labeled with fluorescein, which allows transfection efficiency to be tracked by fluorescent microscopy.

Scan the QR Code or visit the following links

Product Infomation http://www.genelink.com/geneprodsite/product.asp?p=1017



Product Manual http://www.genelink.com/Literature/ps/RNAi_Controls_COA_PS_V6.1.pdf



Product MSDS http://www.genelink.com/Literature/ps/MSDSNH.pdf



Related Products

Product	Catalog No	Size
Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit	27-6402-01	1 Kit
3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each	27-6402-02	1 Kit
Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit	27-6402-06	1 Kit

