



# Product Specification Summary

## siRNA:RNA Interference

Catalog Number 27-6413-10FL  
Product Name Human Beta Actin control siRNA duplex FL; 10 nmols  
Size 10 nmols  
Description Human Beta Actin siRNA; NM\_001101 FI labeled

A negative control should also be included when performing siRNA experiments to eliminate the possibility of nonspecific silencing effects. For this purpose, Gene Link provides a siRNA sequence that has been shown to have no effect on gene silencing. The negative control siRNA is also available labeled with fluorescein, which allows transfection efficiency to be tracked by fluorescent microscopy.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1017>



Product Manual [http://www.genelink.com/Literature/ps/RNAi\\_Controls\\_COA\\_PS\\_V6.1.pdf](http://www.genelink.com/Literature/ps/RNAi_Controls_COA_PS_V6.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

| Product  | Catalog No | Size  |
|--|------------|-------|
| Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit | 27-6402-01 | 1 Kit |
| 3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each         | 27-6402-02 | 1 Kit |
| Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit                                 | 27-6402-06 | 1 Kit |