Product Specification Summary



Omni-PureTM RNA Purification System

Catalog Number 40-4091-10

Product Name Omni-Pure Universal RNA Purification System; 100 purifications

Size 100 Purifications

Description Omni-Pure™ Universal RNA Purification System

Gene Link provides a rapid purification system for extraction of viral DNA or RNA from human bodily fluids including blood. Viral DNA or RNA is captured on a special membrane and then eluted in a low volume for direct use in qualitative and quantitative amplification protocols for detection of a pathogen. The Viral RNA purification system is ideal for small volumes of human bodily fluid samples, i.e., serum, plasma and CSF. Using the easy spin column format, purification of HIV, HCV and other RNA viruses is easily accomplished in less than 30 minutes and ready for RT-PCR amplification.

Scan the QR Code or visit the following links

Product Infomation http://www.genelink.com/geneprodsite/product.asp?p=1083



Product Manual http://www.genelink.com/Literature/ps/M40-4081-XX V4.1 RNA Fluid&Spin.pdf



Product MSDS http://www.genelink.com/Literature/ps/MSDS 40-Omni-Pure.pdf



Related Products

Product	Catalog No	Size
Omni-Pure Viral RNA Spin Column Purification System; 100 purifications	40-3650-01	100 Purifications
Omni-Pure Viral RNA Spin Column Purification System 500 purifications	40-3650-05	500 Purifications

Omni-Pure Blood & Bodily fluids RNA Purification System; 50 purifications 40-4081-05

8 Westchester Plaza, Suite 130, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | www.genelink.com |
SS40-4091-10.pdf Print Date Version: October 16, 2025 Page 1 of 1

