



Product Specification Summary

Omni-Pure™ RNA Purification System

Catalog Number	40-4091-50
Product Name	Omni-Pure Universal RNA Purification System; 500 purifications
Size	500 Purifications
Description	Omni-Pure™ Universal RNA Purification System

Gene Link provides a rapid purification system for extraction of viral DNA or RNA from human bodily fluids including blood. Viral DNA or RNA is captured on a special membrane and then eluted in a low volume for direct use in qualitative and quantitative amplification protocols for detection of a pathogen. The Viral RNA purification system is ideal for small volumes of human bodily fluid samples, i.e., serum, plasma and CSF. Using the easy spin column format, purification of HIV, HCV and other RNA viruses is easily accomplished in less than 30 minutes and ready for RT-PCR amplification.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=1084>



Product Manual <http://www.genelink.com/Literature/ps/M40-4081-XX V4.1 RNA Fluid&Spin.pdf>



Product MSDS http://www.genelink.com/Literature/ps/MSDS_40-Omni-Pure.pdf



Related Products

Product	Catalog No	Size
Omni-Pure Viral RNA Spin Column Purification System; 100 purifications	40-3650-01	100 Purifications
Omni-Pure Viral RNA Spin Column Purification System 500 purifications	40-3650-05	500 Purifications

Omni-Pure Blood & Bodily fluids RNA Purification System; 50 purifications
8 Westchester Plaza, Suite 130, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | www.genelink.com |
SS40-4091-50.pdf Print Date Version : June 13, 2026 Page 1 of 1

