



# Product Specification Summary

## siRNA:RNA Interference

Catalog Number 27-6438-05  
Product Name Custom siRNA duplex PAGE Purified; 50 nmols  
Size 50 nmols  
Description Custom siRNA duplex PAGE Purified, 50 nmols

Custom siRNA are synthesized at different scales of synthesis, the yield varies between these synthesis scales, individual synthesis, modifications and purification methods. Your custom siRNA duplex is accompanied by a Certificate of Analysis stating the exact yield in nmols and  $\mu\text{g}$ . The exact specifications are also printed on the tube label. Please refer to the Certificate of Analysis or the tube label for the exact specifications. Place your custom siRNA order online at the following link <http://www.genelink.com/newsite/allonlineorders.asp> or you can email your sequences at [oligos@genelink.com](mailto:oligos@genelink.com)

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17114>



Product Manual [http://www.genelink.com/Literature/ps/siRNA\\_Custom\\_Duplex\\_PS\\_V1.1.pdf](http://www.genelink.com/Literature/ps/siRNA_Custom_Duplex_PS_V1.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit	27-6402-01	1 Kit
3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each	27-6402-02	1 Kit

Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit 27-6402-06  
8 Westchester Plaza, Suite 130, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | [www.genelink.com](http://www.genelink.com) |  
SS27-6438-05.pdf Print Date Version : June 9, 2026 Page 1 of 1

