



# Product Specification Summary

## siRNA:RNA Interference

Catalog Number	27-6438-01
Product Name	Custom siRNA duplex PAGE Purified; 100 nmols
Size	100 nmols
Description	Custom siRNA duplex PAGE Purified, 100 nmols

Custom siRNA are synthesized at different scales of synthesis, the yield varies between these synthesis scales, individual synthesis, modifications and purification methods. Your custom siRNA duplex is accompanied by a Certificate of Analysis stating the exact yield in nmols and  $\mu\text{g}$ . The exact specifications are also printed on the tube label. Please refer to the Certificate of Analysis or the tube label for the exact specifications. Place your custom siRNA order online at the following link <http://www.genelink.com/newsite/allonlineorders.asp> or you can email your sequences at [oligos@genelink.com](mailto:oligos@genelink.com)

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17115>



Product Manual [http://www.genelink.com/Literature/ps/siRNA\\_Custom\\_Duplex\\_PS\\_V1.1.pdf](http://www.genelink.com/Literature/ps/siRNA_Custom_Duplex_PS_V1.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit	27-6402-01	1 Kit
3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each	27-6402-02	1 Kit

Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit  
8 Westchester Plaza, Suite 130, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | [www.genelink.com](http://www.genelink.com) |  
SS27-6438-01.pdf Print Date Version : April 11, 2026 Page 1 of 1

