



Product Specification Summary

Huntington Disease

Catalog Number	40-2025-15S
Product Name	HD GScan™ Fam amplified control, 44 CAG Repeats; 10 uL
Size	10 uL
Description	HD 6-Fam Control Sample
Component/Note	~44 CAG; ~318bp
Component/Note	Ready to Run

HD GScan™ Fam amplified control. Ready to run control sample, 44 CAG Repeats for fragment analysis

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17187>



Product Manual http://www.genelink.com/Literature/ps/M40-2025-15_Ver5.3.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
GLHD 7 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-05HX	25 uL
GLHD 18 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-01HX	25 uL
GLHD 31 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-07HX	25 uL