



## Product Specification Summary

### Genetic Tools & Reagents

Catalog Number	40-5014-05
Product Name	RNA Reconstitution Solution (1 mM Sodium Citrate pH 6.4); 50 mL
Size	50 mL
Description	RNA Reconstitution Solution (1 mM Sodium Citrate pH 6.4); 50 mL

#### **RNA and RNA Oligonucleotide Reconstitution and Storage Solution**

Sodium citrate, 1 mM pH 6.4 is prepared in RNase free water and sealed in bottles after autoclaving. It is guaranteed RNase free and is the recommended solution for reconstitution and storage of RNA and RNA oligonucleotides. The low pH of 6.4 and low ionic strength of 1 mM Sodium citrate reduces base hydrolysis considerably and is an efficient chelating agent. Sodium citrate, 1 mM pH 6.4 as a reconstitution and storage solution is compatible with all RNA based applications. RNA Reconstitution Solution (1 mM Sodium Citrate pH 6.4) is also available in convenient 10 tubes 1.6 mL aliquot packaging (Catalog No.: 40-5014-16) to meet the rigorous requirements of single use applications.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=17306>



Product Manual

[http://www.genelink.com/Literature/ps/PS40-5014-RNA\\_ReconstSoln\\_Ver2.1.pdf](http://www.genelink.com/Literature/ps/PS40-5014-RNA_ReconstSoln_Ver2.1.pdf)



Product MSDS

[http://www.genelink.com/Literature/ps/SDS\\_40-DNA&RNA\\_solutions.pdf](http://www.genelink.com/Literature/ps/SDS_40-DNA&RNA_solutions.pdf)



## Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL
Omni- Marker™ Universal unlabeled; 500 uL	40-3005-05	500 uL
Omni-Marker™ Universal unlabeled; 1 mL	40-3005-10	1 mL
Omni- Marker™ Low unlabeled; 100 uL	40-3006-01	100 uL
Omni-Marker™ Low unlabeled; 500 uL	40-3006-05	500 uL
Omni- Marker™ Low unlabeled; 1 mL	40-3006-10	1 mL
GScan Marker Tamra Labeled 50 bp - 600 bp; 100 uL	40-3061-01	100 uL
GScan Marker Tamra Labeled 50 bp - 600 bp; 500 uL	40-3061-05	500 uL
Omni-Marker™ GScan-1 Tet labeled 50 bp-1kb	40-3071-01	100 ul
Omni-Marker™ GScan-1 Tet labeled 50 bp-1kb	40-3071-05	500 ul
Omni-Marker™ GScan-2 Tet labeled 50 bp- 600 bp	40-3072-01	100 ul
Omni-Marker™ GScan-2 Tet labeled 50 bp- 600 bp	40-3072-05	500 ul
GScan Marker Hex Labeled 50 bp - 600 bp; 100 uL	40-3081-01	100 uL