



## Product Specification Summary

### GeneProber™ Digoxigenin Labeled Probes

Catalog Number 40-2058-41  
Product Name Mouse 18S rDNA genomic GeneProber Digoxigenin Labeled  
Size 110 uL  
Description Mouse 18S rDNA genomic GeneProber Digoxigenin Labeled  
Component/Note Mouse 18S rDNA genomic GeneProber Digoxigenin Labeled  
Mouse 18S rDNA genomic GeneProber Digoxigenin Labeled for non-radioactive Southern blot analysis.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17361>



Product Manual [http://www.genelink.com/Literature/ps/40-2058-41\\_V1.pdf](http://www.genelink.com/Literature/ps/40-2058-41_V1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



#### Related Products

Product	Catalog No	Size
Kennedy Disease (SBMA) GeneProber Digoxigenin labeled Probe; 100 uL	40-2032-41	110 uL
SCA7 GeneProber Digoxigenin labeled Probe; 110 uL	40-2041-41	110 uL
DRPLA GeneProber Digoxigenin labeled; 110 uL	40-2042-41	110 uL