



Product Specification Summary

Genetic Tools & Reagents

Catalog Number	40-4100-05
Product Name	Omni-Mag™ DNA Purification System Kit; 500 purifications.
Size	500 Purifications
Description	Omni-Mag DNA Purification System Kit; 500 purifications.

Omni-Mag™ DNA Purification System Kit An optimized kit for DNA extractions from trace samples and specifically touch DNA samples for forensic and biometric applications.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17368>



Product Manual http://www.genelink.com/Literature/ps/M40-4100_Omni-Mag_V1.4.pdf



Product MSDS http://www.genelink.com/Literature/ps/SDS_40-4100_Omni-Mag.pdf



Related Products

Product	Catalog No	Size
Omni-Array™ mRNA Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL