



Product Specification Summary

Genetic Tools & Reagents

Catalog Number	40-4118-08
Product Name	Omni-Mag™ PMP (Paramagnetic particles); 800 µL
Size	800 µL
Description	Omni-Mag™ Lysis Buffer (OCC-3); 50 mL
Component/Note	Omni-Mag™ PMP (Paramagnetic particles); 800 µL
Component of Omni-Mag™ DNA purification System (OM). An optimized kit for trace DNA extraction for forensic and biometric applications.	

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17370>



Product Manual http://www.genelink.com/Literature/ps/M40-4100_V1.2_Omni-Mag.pdf



Product MSDS http://www.genelink.com/Literature/ps/SDS_40-4118_PMP_V1.2.pdf



Related Products

Product	Catalog No	Size
Omni-Array™ mRNA Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 µL	40-3005-01	100 µL