



# Product Specification Summary

## siRNA:RNA Interference

Catalog Number	27-6430-05G
Product Name	Pre-designed 3 siRNA PAGE purified; 5 nmols
Size	5 nmols
Description	Pre-designed 3 siRNA PAGE purified; 5 nmols
Component/Note	Negative control siRNA duplex 2 nmols
Component/Note	Pre-designed 3 siRNA RPC purified; 5 nmols each

### **Predesigned RNAi Explorer. \*Minimum 70% silencing.**

You simply give us the accession number or sequence for the gene to down regulate. Gene Link designs siRNA's using published criteria and proprietary algorithms. Gene Link guarantees a minimum of 70% silencing of your gene with at least one of the siRNA supplied. Three target and a negative control siRNA will be supplied.

Guaranteed RNAi Explorer Kit is only available for Human and Mouse Reference Sequence (RefSeq) Accession Numbers.  
- 3 specific target siRNAs custom synthesized by Gene Link, plus 1 negative control. Each siRNA is supplied annealed and ready to use after reconstitution. 2 nmols of negative control siRNA duplex is supplied.

- Guarantee: Gene Link guarantees a minimum of 70% silencing of your gene with at least one of the siRNA supplied. Gene Link will synthesize another set of 3 siRNA if evidence is provided that none of the provided siRNA was functional in achieving a minimum of 70% silencing.

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17431>

Product Manual  [http://www.genelink.com/Literature/ps/RNAi\\_Guar\\_PS\\_Ver5.2.pdf](http://www.genelink.com/Literature/ps/RNAi_Guar_PS_Ver5.2.pdf)

Product MSDS  <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit	27-6402-01	1 Kit
3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each	27-6402-02	1 Kit
Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit	27-6402-06	1 Kit
siRNA non-silencing FL labeled; 20 nmols	27-6411-20FL	20 nmols
Human Vimentin custom siRNA duplex; 20 nmols	27-6412-20	20 nmols
Mouse Beta Tubulin control siRNA duplex; 20 nmols	27-6414-20	20 nmols
Rat Chromogranin-A control siRNA duplex; 20 nmols	27-6415-20	20 nmols
Human Vimentin control siRNA duplex FL label; 10 nmols	27-6412-10FL	10 nmols
Human Beta Actin control siRNA duplex FL; 10 nmols	27-6413-10FL	10 nmols
Rat Chromogranin-A control siRNA duplex; 10 nmols	27-6415-10FL	10 nmols
Custom siRNA duplex standard purity; 2 nmols	27-6436-20	2 nmols
Custom siRNA duplex standard purity; 10 nmols	27-6436-10	10 nmols
Custom siRNA duplex standard purity; 50 nmols	27-6436-05	50 nmols
Custom siRNA duplex standard purity; 100 nmols	27-6436-01	100 nmols
Custom siRNA duplex Reverse Phase Purified; 2 nmols	27-6437-20	2 nmols