



## Product Specification Summary

### siRNA:RNA Interference

Catalog Number	27-6736-15S
Product Name	SmartBase custom siRNA single strand standard purity 150 nmols (~1 mg)
Size	150 nmols
Description	SmartBase custom siRNA single strand standard purity 150 nmols (~1 mg)
Component/Note	SmartBase custom siRNA single strand standard purity 150 nmols (~1 mg)

Custom siRNA are synthesized at different scales of synthesis, the yield varies between these synthesis scales, individual synthesis, modifications and purification methods. Your custom siRNA duplex is accompanied by a Certificate of Analysis stating the exact yield in nmols and  $\mu\text{g}$ . The exact specifications are also printed on the tube label. Please refer to the Certificate of Analysis or the tube label for the exact specifications. Place your custom siRNA order online at the following link <https://www.genelink.com/newsite/allonlineorders.asp> or you can email your sequences at [oligos@genelink.com](mailto:oligos@genelink.com)

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=17527>

Product Manual  [http://www.genelink.com/Literature/ps/siRNA\\_Custom\\_Duplex\\_PS\\_V1.1.pdf](http://www.genelink.com/Literature/ps/siRNA_Custom_Duplex_PS_V1.1.pdf)

Product MSDS  <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



## Related Products

Product	Catalog No	Size
Guaranteed RNAi Explorer Kit w/ Fluorescein Molecular Beacon 3 siRNA duplex; 1 Kit	27-6402-01	1 Kit
3 Target siRNA's & 1 negative control duplex. Standard purity, 2 nmol each	27-6402-02	1 Kit
Guaranteed RNAi Explorer Kit 3 siRNA duplex; 1 Kit	27-6402-06	1 Kit
siRNA non-silencing FL labeled; 20 nmols	27-6411-20FL	20 nmols
Human Vimentin custom siRNA duplex; 20 nmols	27-6412-20	20 nmols
Mouse Beta Tubulin control siRNA duplex; 20 nmols	27-6414-20	20 nmols
Rat Chromogranin-A control siRNA duplex; 20 nmols	27-6415-20	20 nmols
Human Vimentin control siRNA duplex FL label; 10 nmols	27-6412-10FL	10 nmols
Human Beta Actin control siRNA duplex FL; 10 nmols	27-6413-10FL	10 nmols
Rat Chromogranin-A control siRNA duplex; 10 nmols	27-6415-10FL	10 nmols
Custom siRNA duplex standard purity; 2 nmols	27-6436-20	2 nmols
Custom siRNA duplex standard purity; 10 nmols	27-6436-10	10 nmols
Custom siRNA duplex standard purity; 50 nmols	27-6436-05	50 nmols
Custom siRNA duplex standard purity; 100 nmols	27-6436-01	100 nmols
Custom siRNA duplex Reverse Phase Purified; 2 nmols	27-6437-20	2 nmols