



Product Specification Summary

Huntington Disease

Catalog Number 40-2025-32
Product Name Huntington PCRProber Kit Huntington PCR Detection Kit; 1 Kit
Size 1 Kit
Component/Note
Component/Note

PCRProber Products

Gene Link's PCRProber™ Kit is based on PCR amplification followed by Southern blot chemiluminescent detection using an alkaline phosphatase labeled oligonucleotide probe. This kit is a safe and sensitive alternate to radioactive-based detection methods. The amplified products are resolved on a sequencing polyacrylamide gel, and then blotted and processed for chemiluminescent detection. The PCRProber™ Kit is simple and robust for routine triple repeat detection of greater than 100 repeats of all triple repeat disorders Gene Link products.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=200>



Product Manual

http://www.genelink.com/Literature/ps/M40-2025-32_Ver5.3.pdf



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
GLHD 7 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-05HX	25 uL
GLHD 18 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-01HX	25 uL
GLHD 31 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-07HX	25 uL
GLHD 34 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-02HX	25 uL
GLHD 37 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-08HX	25 uL
GLHD 44 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-03HX	25 uL
GLHD 89 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-04HX	25 uL
GLHD 49 CAG repeat Scan Hex labeled ready to run Control DNA; 25 uL	40-2025-09HX	25 uL
GLHD 116 CAG repeat GScan Genemer Control DNA; 25 uL	40-2025-06HX	25 uL
GLHD 134 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-61HX	25 uL
GLHD 182 CAG repeat GScan Hex labeled ready to run Control DNA; 25 uL	40-2025-62HX	25 uL
Huntington Disease GScan™ TET Kit	40-2025-15TT	1 Kit
Huntington Disease GScan™ HEX Kit	40-2025-15HX	1 Kit
Huntington Disease GScan™ Cy3 Kit	40-2025-15C3	1 Kit
Huntington Disease GScan™ Cy5 Kit	40-2025-15C5	1 Kit