



Product Specification Summary

Myotonic Dystrophy

Catalog Number	40-2026-31
Product Name	Myotonic Dystrophy PCRProber (Oligo AP labeled); 5 blots
Size	5 blots
Description	Gene Detection Systems
Component/Note	labeled oligo probe
Component/Note	Myotonic Dystrophy

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=201>



Product Manual http://www.genelink.com/Literature/ps/M40-2026-31_Ver3.1.pdf



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng	40-2026-39	500 ng
Myotonic Dystrophy GScan™ TET Kit	40-2026-15TT	1 Kit
Myotonic Dystrophy GScan™ HEX Kit	40-2026-15HX	1 Kit