



# Product Specification Summary

## Myotonic Dystrophy

Catalog Number 40-2026-32  
Product Name Myotonic Dystrophy PCRProber kit; 1 Kit  
Size 1 Kit  
Component/Note  
Component/Note

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=202>



Product Manual [http://www.genelink.com/Literature/ps/M40-2026-32\\_Ver3.1.pdf](http://www.genelink.com/Literature/ps/M40-2026-32_Ver3.1.pdf)



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng	40-2026-39	500 ng
Myotonic Dystrophy GScan™ TET Kit	40-2026-15TT	1 Kit
Myotonic Dystrophy GScan™ HEX Kit	40-2026-15HX	1 Kit