



# Product Specification Summary

## Sry, X & Y

Catalog Number	44-6430-01
Product Name	SRY TaqMan 5'FAM; 25 ug
Size	25 ug
Description	Fam/Tamra TaqMan
Component/Note	Sry Genemer Specific

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=345>



Product Manual <http://www.genelink.com/Literature/ps/PS44-6430-01.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
X alphoid TaqMan 5'FAM; 25 ug	44-6130-01	25 ug
X alphoid TaqMan 5' Fluorescein; 25 ug	44-6140-01	25 ug
SRY TaqMan 5' Fluorescein; 25 ug	44-6440-01	25 ug