



Product Specification Summary

Genetic Tools & Reagents

Catalog Number 40-5032-00
Product Name 6% Polyacrylamide Gel mix; 125 mL
Size 125 mL

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=41>



Product Manual <http://www.genelink.com/Literature/ps/ps40-5032-00.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Omni-Array™ mRN Antisense strand amplification kit	08-0021-02	1 Kit
Omni-Array™ Sense Strand mRNA Amplification Kit	08-0011-02	1 Kit
Omni-Marker™ Universal unlabeled; 100 uL	40-3005-01	100 uL