



Product Specification Summary

Sry, X & Y

| | |
|----------------|----------------------|
| Catalog Number | 45-6430-01 |
| Product Name | SRY MB-FAM; 25 ug |
| Size | 25 ug |
| Description | Sry Fam/Dabcyl MB |
| Component/Note | Sry Genemer Specific |
| Component/Note | Sry Molecular Beacon |

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=441>



Product Manual <http://www.genelink.com/Literature/ps/PS44-6430-01.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

| Product | Catalog No | Size |
|--|------------|-------|
| X alphoid TaqMan 5'FAM; 25 ug | 44-6130-01 | 25 ug |
| X alphoid TaqMan 5' Fluorescein; 25 ug | 44-6140-01 | 25 ug |
| SRY TaqMan 5' Fluorescein; 25 ug | 44-6440-01 | 25 ug |