



# Product Specification Summary

## RT-PCRmer™

Catalog Number	40-1002-00
Product Name	beta actin control PCR mix (human & Rat); 200 uL
Size	200 uL
Description	Human, Rat and Mouse
Component/Note	Ready to use
Component/Note	Add Template &Taq Pol.

html head meta http-equiv="Content-Type" content="text/html; charset=windows-1252" meta name="GENERATOR" content="Microsoft FrontPage 4.0" meta name="ProgId" content="FrontPage.Editor.Document"> title>New Page 1/title> /head> /title> '>

Scan the QR Code or visit the following links

Product Infomation <http://www.genelink.com/geneprodsite/product.asp?p=52>



Product Manual <http://www.genelink.com/Literature/ps/ps40-1002-00.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



### Related Products

Product	Catalog No	Size
beta actin D Human RT-PCRmer™ (F2/R2) Exon 4-3' U; 10 nmols	40-1023-10	10 nmols
beta actin G Human RT-PCRmer™ (F4/R1) Exon 2-3; 10 nmols	40-1025-10	10 nmols

beta actin B Mouse RT-PCRmer (F2/R2) Exon 6-7: 10 nmols  
8 Westchester Plaza, Suite 130, Hawthorne, NY 1032-1515 | Tel: 914-769-1192 | Fax: 914-769-1193 | www.genelink.com |  
SS40-1002-00.pdf Print Date Version : April 15, 2026 Page 1 of 1

