



Product Specification Summary

Myotonic Dystrophy

Catalog Number	40-2026-41
Product Name	Myotonic Dystrophy CTG Repeat Genotyping GeneProber™ Digoxigenin labeled ; 110 uL
Size	110 uL
Description	Myotonic Dystrophy CTG Repeat Genotyping GeneProber™ Digoxigenin labeled ; 110 uL
Component/Note	CTG Triple repeat region
Component/Note	Myotonic Dystrophy CTG Repeat Genotyping GeneProber™ Digoxigenin labeled ; 110 uL

Gene Link is the leader in triple repeat disorder genotyping using non-radioactive based methods. We have more than a decade of expertise and have developed facile non-radioactive detection methods for safe, sensitive and reliable genotyping of human genetic disorders specifically fragile X, Huntington Disease, Myotonic Dystrophy etc. Take a look at the simple agarose and polyacrylamide gel based systems, the chemiluminescent Southern blot detection methods and the fluorescent systems for genotyping of triple repeat disorders.

Scan the QR Code or visit the following links

Product Information

<http://www.genelink.com/geneprodsite/product.asp?p=632>



Product Manual

http://www.genelink.com/Literature/ps/M40-2026-41_Ver6.2.pdf



Product MSDS

<http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Myotonic Dystrophy GeneProber GLDM2 (Pst I) unlabeled Probe; 500 ng	40-2026-39	500 ng
Myotonic Dystrophy GScan™ TET Kit	40-2026-15TT	1 Kit
Myotonic Dystrophy GScan™ HEX Kit	40-2026-15HX	1 Kit
Myotonic Dystrophy GScan™ Cy3 Kit	40-2026-15C3	1 Kit
Myotonic Dystrophy GScan™ Cy5 Kit	40-2026-15C5	1 Kit
Myotonic Dystrophy GeneProber GLDM3 unlabeled Probe; 500 ng	40-2026-38	500 ng
Myotonic Dystrophy GScan™ 6-FAM Kit; 1 Kit	40-2026-15FM	1 Kit
Myotonic Dystrophy GScan™ 6-FAM Kit (20rxns); 1 Kit	40-2026-15FMS	1 Kit
Myotonic Dystrophy Genemer Kit; 1 Kit	40-2026-11	1 Kit
DM CTG-6-MB-CY3; 25 ug	45-0210-06	25 ug
DM CTG-6-MB-CY5; 25 ug	45-0220-06	25 ug
DM CTG-6-MB-FAM; 25 ug	45-0230-06	25 ug
DM CTG-6-MB-FI; 25 ug	45-0240-06	25 ug
DM CTG-6-MB-HEX; 25 ug	45-0250-06	25 ug
DM CTG-6-MB-TET; 25 ug	45-0260-06	25 ug