



Product Specification Summary

GeneProber™ Digoxigenin Labeled Probes

Catalog Number 40-2032-41
Product Name Kennedy Disease (SBMA) GeneProber Digoxigenin labeled Probe; 100 uL
Size 110 uL

Scan the QR Code or visit the following links

Product Information <http://www.genelink.com/geneprodsite/product.asp?p=655>



Product Manual <http://www.genelink.com/Literature/ps/ps40-2032-41.pdf>



Product MSDS <http://www.genelink.com/Literature/ps/MSDSNH.pdf>



Related Products

Product	Catalog No	Size
Kennedy Disease (SBMA) GeneProber Digoxigenin labeled Probe; 100 uL	40-2032-41	110 uL
SCA7 GeneProber Digoxigenin labeled Probe; 110 uL	40-2041-41	110 uL
DRPLA GeneProber Digoxigenin labeled; 110 uL	40-2042-41	110 uL